





Instructions for Use

Power Knee™
PKA01

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	16
FR	Notice d'utilisation	28
ES	Instrucciones para el uso	41
IT	Istruzioni per l'uso	53
NO	Bruksanvisning	65
DA	Brugsanvisning	77
SV	Bruksanvisning	88
EL	Οδηγίες Χρήσης	100
FI	Käyttöohjeet	113
NL	Gebruiksaanwijzing	125
PT	Instruções de Utilização	137
PL	Instrukcja użytkowania	150
TR	Kullanım Talimatları	162
RU	Инструкция по использованию	174
JA	取扱説明書	187
ZH	中文说明书	199
KO	사용 설명서	210
CS	Návod k použití	222
SK	Indikácie použitia	234
RO	Indicații de utilizare	246
HR	Upute za uporabu	259
HU	Használati javallat	271
BG	Указания за употреба	283
SL	Indikácie pre použitie	296

ÖSSUR
POWER KNEE™ PKA01

CIN: XXXXXX
Input 24V
1.25A(max)
Battery Operated
33.6V 15A(max)
FCC ID:2ADHKATWINC1500
IC:20266-WINC1500PB
FCC ID:2AASB04; IC:12208A-04

2797

IP54

Össur Iceland ehf, Grjótháls 5, 110 Reykjavík, Iceland

ÖSSUR

POWER KNEE™ BATTERY

Input 24V 1.25A MAX
Output 33.6V

Rechargeable Li-Ion Battery - 8ICR1966
Nominal Voltage 28.8V
3.3Ah 95.04 Wh

REF PKA01102 SN HF55XXXX

Össur Iceland ehf, Grjótháls 5, 110 Reykjavík, Iceland

IP67

CAUTION:
DO NOT
DISASSEMBLE,
CRUSH, HEAT
ABOVE 75°C OR
INCINERATE.

Input: 100-240V~/50-60Hz/600-300mA
Output: 24V ~ 1/250mA
Type: FW8030M/24
Made in XXXXXX
FRUWO Gerätebau GmbH
Von-Liebig-SträÙe 11
48346 Ostbevern
Oscur Japan G.K.

CAUTION:
RISK OF ELECTRICAL SHOCK
DRY LOCATION USE ONLY

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DESCRIPTION

The device is a microprocessor-controlled, and motor-powered prosthetic knee. It enables patients to regain mobility and perform activities of daily living.

The device consists of the following main components (**Fig. 1**):

1. Proximal pyramid adapter;
2. Motor;
3. Power button and status indicator;
4. Security tab;
5. Battery Release Clip;
6. Battery;
7. Charging port;
8. Battery indicator and battery button;
9. Distal pyramid Connector;

Device labels (**Fig. 5**) can be found at the rear of the device, on the power supply and on the internal face of the battery pack.

This document is for patients (intended device operator) and healthcare professionals.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

Unilateral transfemoral or knee disarticulation amputation.

Unilateral hip-disarticulation or hemipelvectomy amputation.

Bilateral amputations that combine unilateral amputation listed above on one side with transfemoral level amputation or any amputation below that level on the contralateral side.

No known contraindications.

The device is for low to moderate impact use, e.g., walking.

The weight limit for the device is between 50 kg and 116 kg (110 lbs and 256 lbs).

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: In case of an unexpected power loss the device will stop moving and may not return to full extension.

Warning: Avoid placing hands or fingers near moving joints.

Warning: If the device gives a warning signal – sound or vibration, immediately stop walking and proceed with caution.

Warning: If device emits a warning signal, stop using device and refer to **Table 1** for more information. Note that certain errors can cause the device to shut down.

Caution: The device is motor powered and can move on its own.

Note: High intensity uses or operation in high ambient temperature will cause the device actuator to heat up and become hot to the touch.

Note: Improper handling and/or adjustment of the device can cause malfunction. Patient should avoid:

- Direct impact to the device;
- Excessive shocks or vibrations;
- High-impact activities, sports, excessive loading, and heavy-duty use.

Note: Do not attempt to remove covers or modify device in any way.

The device is for single patient use.

REQUIRED COMPONENTS

The device is delivered with a power supply and a Power Knee™ Battery.

Össur Logic App

For the initial fitting, the healthcare professional must use the Össur Logic app. It is available on Apple App Store. See Össur Logic on the App Store for mobile device compatibility.

The patient can use the Össur Logic app to connect to the device and access some functions.

DEVICE SELECTION

See Össur catalog for a list of recommended components.

Use Össur titanium adapters to achieve weatherproof performance.

Note: Adapter installation and set-screw tightening shall be performed according to instructions provided in their respective accompanying documents.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass 0-5 mm posterior to knee axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes 0-5 mm posterior to knee axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Static Alignment

Turn the device on.

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

1. Keep device powered ON.
2. Ensure the patient is familiar with how the device functions.
3. Tell the patient to walk between parallel bars with usual gait. This lets the patient get used to the behavior of the device.
4. Adjust dynamic alignment as required.
 - Make sure that the motion in swing and stance phases stays in the line of progression.
 - Make sure to achieve step-length symmetry.
5. The Össur Logic App must be used for setting device parameters during dynamic alignment.

Note: Use the "Initial Setup" function in the Össur Logic app to adjust the device for basic use. Use the "Advanced Settings" function to fine tune the device.

USAGE

Battery

Initial charge

The battery pack is in sleep mode during shipping and is charged to approximately 30%. To activate it, charge the battery pack. Fully charge the battery pack before first use.

Charging

1. Plug the power supply into the wall socket.
2. Plug the power supply into the charging port. If the battery pack is in the device when it is connected to the charging port, the device will turn off.
3. After a short while, the battery indicator will show the battery charge level with solid green lights and indicate charging with a flashing green light. If the power supply is disconnected during charging, the battery indicator may continue flashing for a short while.
4. After charging is complete, the battery indicator lights will turn off.
5. After charging, disconnect the power supply from the charging port.
6. Disconnect the power supply from the wall socket to save energy.

Warning: A damaged power supply can lead to electric shock.

Warning: Always disconnect the power supply from the device before donning.

Caution: The Battery Pack charge connector is Nickel-plated which may cause allergic reactions.

Note: Charge the battery pack indoors.

Note: It takes approximately 5 hours to charge the device battery pack when fully depleted.

Note: A fully charged battery pack can last for 4 to 20 hours, depending on activity level.

Note: The battery pack cannot be overcharged.

Note: The battery pack can be charged while inserted in the device or removed for charging.

Note: The battery pack can only be charged in temperatures between 10°C and 45°C. If charge is attempted outside these limits, it will be delayed until the battery pack temperature is within the specified range.

Battery Charge Level Indication

Push the battery indicator button. The green lights on the battery indicator show approximate battery level.

Removing and inserting the battery pack

Caution: Never remove the battery pack unless the device is powered off as data could be lost.

How to Remove the Battery Pack:

1. Turn the device off to avoid losing data collected during device operation.
2. Pull out the security tab at the top of the battery pack (**Fig. 4a**).
3. Push down the battery release clip at the top of the battery pack until you hear a click (**Fig. 4b**). The battery pack will release from the battery compartment.
4. Remove the battery pack.

Caution: Only remove the battery pack in a dirt-free environment.

How to Insert the Battery Pack:

1. Lift the security tab up and out of the way of the battery compartment.
2. Align the battery pack with the battery compartment (**Fig. 4c**).
3. Push the battery pack in at the top until you hear the battery release clip click (**Fig. 4d**).
4. Check if the battery pack is securely mounted and flush with battery compartment.
5. Push the security tab into position flush with the battery pack (**Fig. 4e**).

Warning: The device can turn off unexpectedly if battery pack is not inserted correctly. Make sure the battery pack and security tab are fully in place before using device.

Device operation

Turning the Device On

1. Push and hold the power button for longer than 3 seconds.
2. The device will play three ascending tones and vibrate.
3. When the status indicator flashes green, the device is ready to use.

Turning the Device Off

1. Push and hold the power button for longer than 3 seconds.
2. The device will play three descending tones and turn off.

Connecting to the mobile app

1. Turn on the device.

2. Open the app and follow the instructions.

Note: The required PIN code can be found on the device label at the back of the device.

In case of warnings, stop using device and troubleshoot as follows:

1. Check battery status and recharge if necessary.
2. Check if the device is hot and allow to cool down.

If this does not help, contact your healthcare professional.

The Össur Logic App can also give information on device warnings.

Supported activities

The device automatically adapts to activities such as walking, standing up and more. It has several different states that the patient must learn to recognize and activate.

Standing State

The device is in standing state by default and it will revert to standing state if it does not recognize a movement pattern.

The device will provide support when the patient applies weight to the prosthesis. The knee will release/swing freely when the prosthesis is unloaded/no weight is applied to it. The patient can take small steps and turns in a natural and physiological way.

Standing Lock

Standing lock has been modified to reduce incidence of overheating after activation of the lock.

Standing lock allows standing with the knee locked in a partially flexed position.

To enter standing lock:

1. Remove weight from the prosthesis.
2. Flex the knee up to the desired position.
3. Put some weight on the device and hold in a steady position for 3 seconds.
4. The knee is now in a standing locked state and full weight may be applied to it.

To exit standing lock take weight off the device, this will return the knee into standing state.

Note: When engaging the stance lock with high degree of flexion with full load, the internal temperature of the actuator will increase.

To prevent overheating, the knee will transition from standing lock to stance yielding function following delivery of a device temperature warning. It is recommended that the standing lock is used for short-term standing support up to 30 seconds.

Walking on Level Ground

The device adapts to the patients walking speed and style. It will give adaptive support in stance phase and powered extension and flexion in swing phase.

If using a cane or a crutch, the patient must take care to put weight on the knee during stance phase or the knee will not provide support.

Sitting Down

Sitting down has been modified to make it easier and more natural to move from standing to sitting.

The device automatically detects sitting down motion. Support is provided while the patient sits down. The patient can use this to distribute weight equally to both legs and control how fast to sit down.

To sit down:

1. Stand in front of the chair with even weight on both feet.
2. Lean back toward the chair. The knee will provide support during knee flexion until seated.
3. The patient can now take the weight off the leg to release the knee and move the knee freely while seated.

After a few seconds of sitting without moving the knee, the knee will enter standby. Standby turns off the motor to save power.

To exit standby, flex or extend the knee slightly.

Standing Up

The device will recognize when the patient starts to stand up. The device will provide a powered extension to let the patient apply equal weight to both legs.

To stand up:

1. Ensure that the prosthetic foot is placed directly under the knee.
2. Put equal weight on both feet.
3. Start to stand up. The patient will feel the knee assist standing up naturally.

Note: Placing hands on knees helps applying weight on both legs.

Walking Up Stairs

The device will provide a powered knee extension when it detects stair ascent.

Initial training:

1. Stop in front of the first step of the stairs.
2. Flex the hip to lift the prosthesis off the ground. The knee will swing freely.
3. If needed, swing the leg outwards to get the leg onto the first step to prevent the toe from catching the step.
4. Put the prosthetic foot flat onto the first step.
5. Put weight on the prosthesis and extend the knee fully. The knee will detect the stairs and will provide a powered extension.
6. When the knee has fully extended, place the other foot on the same step.
7. Take the weight off the prosthesis and if stair ascent has been triggered, the knee will flex.
8. After a brief pause the knee will automatically extend.
9. Practice going up and down the first step, repeating the above steps until the patient is familiar with the knee action.

To walk up stairs, the patient should:

1. Stop in front of the first step of the stairs.
2. Flex the hip to lift the prosthesis off the ground. The knee will swing freely.
3. If needed, swing the leg outwards to get the leg onto the first step to prevent the toe from catching the step.
4. Put the prosthetic foot flat onto the first step.
5. Put weight on the prosthesis and extend the knee fully. The knee will detect the stairs and will provide powered extension.
6. When the prosthetic knee has fully extended, place the other foot on the step above.
7. Flex the hip on the prosthetic side to bring the prosthesis to the next step. The prosthesis will provide knee flexion to provide clearance and to prepare for the next step.
8. Continue climbing the stairs, with a step over step pattern.
9. At the top of the stairs:
 - a. If prosthetic foot is leading, place the other foot next to it, and keep equal weight on both feet for 3 seconds before continuing.
 - b. If the other foot is leading, flex the hip to place the prosthetic knee slightly in front of the other knee until the prosthetic knee extends.

Walking Up Ramps

Walking up ramps requires no change from level ground walking.

Caution: Always use the handrail when walking up ramps or stairs.

Walking Down Stairs and Ramps

Caution: Always use the handrail when you walk down ramps and stairs.

The knee will give support while it flexes when walking down ramps and stairs.

To walk down stairs or ramps:

1. Take the first step down with the prosthesis.
 1. Put weight on the prosthesis. Lean slightly backwards to flex the knee. The knee will flex and give support.
 2. Lower the other foot onto the next step or ramp.
 3. Continue to walk down the stairs or ramp. The knee will adapt the support to the walking speed.

Note: Leaning backwards keeps from falling forward if balance is lost.

Kneeling

The device automatically detects kneeling motion. The device will provide support until the knee reaches the ground.

To kneel:

1. Take a step forward with the sound leg. If needed, use a hand support to maintain balance.
2. Take most of the weight off the prosthesis momentarily.
3. Put weight on the prosthesis.
4. Within 3 seconds, slightly flex the hip on the prosthetic side and hold pressure until the knee starts to flex.
5. Kneel. The device gives support while kneeling.

Driving

When driving a vehicle, the device must be turned off.

Device warning signals

The device may give warnings in the form of audible beeps and vibration. In case of warnings, refer to **Table 1** for a list of possible causes and take action.

In case of warnings, stop using device and troubleshoot as follows:

1. Check battery status and recharge if necessary.
2. Check if the device is hot and allow to cool down.

If this does not help, contact your healthcare professional.

The Össur Logic App can also give information on device warnings.

Table 1 - Warning Signals

Warning type	Device state	Device Feedback	Action
Low battery warning	Low battery (15% charge)	Beeps twice. Vibration for 0.5 seconds at next load cycle. Device status indicator flashing orange. Full gait function active.	Charge battery or replace battery.
	Critical low battery (10% charge)	Beeps 10 times every 7.5 seconds. Vibration 3 times every 10 seconds. Device status indicator flashing red. Full gait function active for 60 seconds until device triggers shutdown with 3 descending tones.	Charge battery or replace battery.
Device temperature warning	High motor/device temperature	Beeps twice. Vibration for 0.5 seconds at next load cycle. Device status indicator flashing orange. Full gait function active.	Stop using device or reduce activity, allow to cool down.
	Critical high motor/device temperature	Beeps 10 times every 7.5 seconds. Vibration 3 times every 10 seconds. Device status indicator flashing red. Full gait function active for 60 seconds until device triggers shutdown with 3 descending tones.	Stop using device or reduce activity, allow to cool down. Send to Össur for service.
Device temperature warning while using flexed device for support	High motor temperature	Beeps once. Vibration 2 times. Gait function switches from locked to yielding.	Stop using device for full support at an angle.
Non-gait related sensor error	Sensor error.	Beeps 10 times every 7.5 seconds. Vibration 3 times every 10 seconds. Device status indicator flashing red. Full gait function active for 60 seconds until device triggers shutdown with 3 descending tones.	Stop using device. Send to Össur for service.
Gait related sensor error	Sensor error.	Beeps 10 times every 7.5 seconds. Vibration 3 times every 10 seconds. Device status indicator flashing red. Reduced gait function (device extends under load) active for 60 seconds until device triggers shutdown with 3 descending tones.	Stop using device. Send to Össur for service.
Battery critical error	Sensor error	Device stops operating. Device triggers shutdown with 3 descending tones.	Stop using device. Send to Össur for service.

Cleaning and care

It is recommended to maintain the device in good outward condition by cleaning the surface of the device regularly.

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Caution: Do not use solvents stronger than isopropyl alcohol as it might degrade the material.

Caution: Do not use compressed air or pressurized water jets to clean device.

Note: Device is not delivered in a sterilized state or intended to be sterilized.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Caution: If the device is exposed to salt water or chlorinated water it must be sent to Össur for servicing.

See which environmental conditions apply to the device in **Table 2**.

Table 2 - Environmental Conditions

	Use	Charging	Shipping	Extended storage
Temperature	-10 °C to 45 °C (14 °F to 113 °F)	10 °C to 40 °C (50 °F to 105 °F)	-25 °C to 70 °C (-13 °F to 158 °F)	0 °C to 45 °C (32 °F to 113 °F)
Relative humidity	0% to 100%, condensing	0% to 100%, relative humidity (non-condensing)	10% to 100%, condensing	10% to 90%, non-condensing
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

If the device is used outside the allowed temperature range, device performance may be unexpected and battery warnings may occur (See **Table 1**). If the device is stored below 0°C (32°F) for an extended time, device will not perform as expected. When the device warms up, performance will be restored.

Warning: Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.

MAINTENANCE

Examine the device regularly for signs of damage. If you find signs of damage, consult your healthcare professional.

Warning: Do not attempt device repairs or remove covers. Risk of injury or damage to device. Send device to Össur for maintenance or repair.

Expected life duration under normal usage is 6 million steps. Depending on patient activity, this corresponds to a duration of use of 4 to 6 years.

SPECIFICATIONS

Table 3 - Device Specifications

Device weight	2.65 kg (5.8 lbs)
Device build height (Fig. 2)	27.7 cm (10 7/8")
Battery pack specification	Lithium-Ion / 3300 mAh / 95.04Wh
Battery pack weight	521 g (1.15 lbs)
Energy to charge the battery pack	251.7 VAh at 230 V AC, 50 Hz 216.5 VAh at 110 V AC, 60 Hz
Primary packaging weight	3.2 kg (7.1 lbs)
Primary packaging materials	Polypropylene (PP 5) plastic and Polyurethane Foam

Note: The primary packaging is made from recyclable materials.

Power supply specifications

Table 4 - Power Supply Specifications

Manufacturer	FRIWO
Model no.	FW8030M/24
Input	100–240 VAC, 50–60 Hz, 0.6–0.3 A
Output	24 VDC, 1.25A

Caution: Only use the provided power supply. Do not use a different power supply.

ELECTROMAGNETIC COMPATIBILITY

Warning: Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The device is suitable for use in any environment except where immersion in water or any other fluid is possible, or in flammable / explosive atmospheres, or where exposure to highly electrical and/or magnetic fields can occur (e. g., electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

The device may be susceptible to electromagnetic interference from portable and mobile RF communications equipment such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

Table 5 - Electromagnetic Emissions Compliance Information

Emissions test	Standard	Compliance level
Conducted and radiated RF emissions	CISPR 11	Group 1 – Class B
Harmonics emissions	IEC 61000-3-2	Class A
Voltage fluctuations / Flicker emissions	IEC61000-3-3	Complies

Table 6 - Electromagnetic Immunity Compliance Information

Immunity Test	Standard	Compliance Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 15 kV air
Electrical fast transient/burst	IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines
Surge	IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage dips	IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field	IEC 61000-4-8	30 A/m
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
	EN ISO22523	12 V/m 26 MHz – 1 GHz
Immunity to proximity magnetic fields	IEC 61000-4-39	134.2kHz @ 65A/m 3.56 MHz @ 7.5A/m 30kHz @8A/m

Table 7 - RF wireless communications Compliance Information

Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

Wireless modules regulated information

Table 8 - Wireless Modules

Model	Regulatory certificates	Type and frequency characteristics	Effective radiated power
WiFi Module Model: ATWINC1500- MR210PB	<p>FCC Contains FCC ID: 2ADHKATWINC1500</p> <p>Canada Contains transmitter module IC: 20266-WINC1500PB</p> <p>Japan Contains transmitter with certificate number  005-101762</p> <p>Korea Contains transmitter with certificate number  R-CRM-mcp-WINC1510MR210P</p> <p>Taiwan Contains transmitter with certificate number  CCAN18LP0321T2</p>	2412 - 2462 MHz (5 MHz Interval, 11 Channels)	N/A

Model	Regulatory certificates	Type and frequency characteristics	Effective radiated power
Bluetooth 4.2 Module Model: BMD-300	FCC Contains FCC ID: 2AA9B04	F1D 2402-2480 MHz	0.35 mW/MHz
	Canada Contains IC: 12208A-04		
	Japan Contains transmitter with certificate number  R210-106799		
	Korea Contains transmitter with certificate number  R-CRM-Rgd-BMD-300		
	Brazil Contains ANATEL approved module # 00820-21-05903		
	Mexico Este equipo contiene el módulo con IFT #: NYCE/CT/0146/17/TS		

USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation – This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- this device may not cause interference, and
- this device must accept any interference, including interference that may cause undesired operation of this device.

Caution: Exposure to Radio Frequency Radiation – The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website www.hc-sc.gc.ca/rpb

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance


The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

The device complies with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment (ROHS3).

The device complies with Directive 2014/53/EU of the European Parliament and of the Council of 16 April 2014 on the making available on the market of radio equipment.


This device has been tested according to ISO 10328 standard to three million load cycles.

The load level specification is ISO10328 - P5 - 116 kg

SYMBOLS

Symbol	Description
	Medical Device
	Caution Symbol
	Type B Applied Part
	Consult instructions for use
	Device contains electronic components and/or batteries that should not be disposed of in regular waste
	Serial number
	Year of manufacture
	Model number
	See operating instructions

Symbol	Description
	Fragile, handle with care
	Keep Dry
	For indoor use only
	ON / OFF
	Non-ionizing electromagnetic radiation
	Alternating current
	Direct current
IP54	Protected from harmful effects of dust and against the effects of splashing water
IP67	Protected against harmful effects of dust and temporary immersion in water Note: This only applies to the battery pack.
	QR Code containing device specific wireless communication PIN Code
	Positive center connection
	Surface may be hot
 Li-ion	Recyclable Lithium Ion Battery
	This way up
	Recyclable symbol
	Class II Equipment
	Energy Efficiency Level VI
 CAUTION: DO NOT DISASSEMBLE, CRUSH, HEAT ABOVE 75°C OR INCINERATE.	Do not disassemble, crush, heat above 75 °C or incinerate battery pack

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