
Instructions for Use




---

PROPRIO FOOT

PSX01

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	23
FR	Notice d'utilisation	43
ES	Instrucciones para el uso	62
IT	Istruzioni per l'uso	82
NO	Bruksanvisning	101
DA	Brugsanvisning	119
SV	Bruksanvisning	138
EL	Οδηγίες Χρήσης	157
FI	Käyttöohjeet	176
NL	Gebruiksaanwijzing	195
PT	Instruções de Utilização	215
PL	Instrukcja użytkowania	235
CS	Návod k použití	255
TR	Kullanım Talimatları	274
RU	Инструкция по использованию	293
JA	取扱説明書	313
ZH	中文说明书	332
KO	사용 설명서	350

## DESCRIPTION OF SYMBOLS

Caution Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

ON / OFF

Type BF Applied Part (IEC60601-1)

Non-ionizing electromagnetic radiation

Manufacturer and year of manufacturing

Direct current

Consult instructions for use  
**Note:** symbol can appear in white on blue background.

This marking on the product, packaging, accessories or literature indicates that the product contains electronic components and/or batteries that should not be disposed of in regular waste at the end of its usable life. To prevent possible harm to the environment or human health from uncontrolled waste disposal users are requested to separate these items from other types of waste and recycle them responsibly to support sustainable reuse of material resources. Users should contact their local government office for information on how these items can be recycled or disposed of in an environmentally sound manner.  
To protect natural resources and to promote material reuse, please separate batteries and electrical components from other types of waste and recycle them through your local, free electronic parts return system.

Protected against solid foreign objects of 2.5 mm diameter and greater and protected against splashing water

Serial number

CE Mark with Notified Body Identification Number

Device specific wireless communication PIN Code

## DESCRIPTION OF SYMBOLS ON PRIMARY PACKAGING

This way up.  
Indicates correct upright position of the transport package.

Consult instructions for use  
Indicates the need for the user to consult the instructions for use.

Fragile, handle with care  
Indicates a medical device that can be broken or damaged if not handled carefully.

Recyclable Packaging  
Indicates that the packaging can be recycled.

Caution  
Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.

Keep Dry.  
Indicates that the packaged medical device needs to be protected from moisture.

This marking on the product, packaging, accessories or literature indicates that the product contains electronic components and/or batteries that should not be disposed of in regular waste at the end of its usable life.

## PROPRIO FOOT®

The PROPRIO FOOT® is referred to as “*the device*” in the following document. This document provides important information on the indications for use and handling of the device. It is intended for a certified prosthetist and the user of the device. The device may only be fitted by a qualified practitioner authorized by Össur after completing the corresponding training.

## PRODUCT DESCRIPTION

The device is a microprocessor prosthetic ankle and foot system, consisting of the following components:

1. Charge connector
2. Actuator
3. Pyramid adapter
4. Power button
5. LED Indicator Ring
6. Unit serial number
7. Carbon blades

The device increases toe clearance in swing phase by raising the forefoot momentarily in mid-swing. It adjusts the ankle angle to changing terrain, providing stable foot placement on inclined surfaces. Adjustments only occur in swing phase when the ankle is not under load.

---

## INDICATIONS FOR USE

---

The device is to be used exclusively for the exoprosthetic fitting of amputations of the lower limb.

- Unilateral trans-tibial amputation
- Bilateral trans-tibial amputation
- Unilateral trans-femoral amputation

**Note:** Use of the device by patients with the following conditions is subject to a case-by-case assessment, and advanced fitting performed by a qualified Össur practitioner:

- Bilateral transfemoral amputation
- Limited residual limb control, due to e.g. unilateral hip disarticulation and lower-extremity dysmelia patients

**Note:** User should have sufficient control over residual limb to support device activities.

---

## CONTRAINDICATIONS

---

None known.

---

## INTENDED USE

---

The device is a non-invasive, reusable prosthetic device designed for single patient use.

The device supports low to moderate impact levels.

Maximum user weight is 125 kg (275 lbs.).

---

## SPECIFICATIONS (Figure 2)

---

Device build heights and weights:

Size 22–24: 174 mm (6 <sup>7</sup> / <sub>8</sub> ”)	Weight: 1.3 kg (2.9 lbs)	Size 24, Category 2.
Size 25–27: 180 mm (7 <sup>1</sup> / <sub>16</sub> ”)	Weight: 1.5 kg (3.3 lbs)	Size 27, Category 5.
Size 28–30: 188 mm (7 <sup>3</sup> / <sub>8</sub> ”)	Weight: 1.6 kg (3.5 lbs)	Size 30, Category 8.

Device ankle motion range is 33°.

---

## CATEGORY SELECTION

---

The table below defines category selection as a function of user weight and impact level.

Weight kg	45–52	53–59	60–68	69–77	78–88	89–100	101–116	117–125
Weight lbs.	99–115	116–130	131–150	151–170	171–194	195–221	222–256	257–275
Low Impact Level	1	1	2	3	4	5	6	7
Moderate Impact Level	1	2	3	4	5	6	7	8

---

## PACKAGING

---

The primary packaging of the device is made of recyclable materials: Polypropylene (**PP 5**) plastic and polyurethane foam.

Primary packaging weights 3.1 kg / 6.8 lbs.

---

## REQUIRED DEVICES

---

The Össur Logic mobile app requires an Apple iOS device supported by the manufacturer, e.g. an iPhone or iPad.

---

## SAFETY INSTRUCTIONS

---

### **Warnings**

Össur products and components are designed and tested according to applicable official standards, or an in-house defined standard should no official standard apply. Compatibility and compliance with these standards is achieved only if Össur products are used with recommended components. Please refer to the Össur catalog for recommended device components.

Immediately cease using the device if unusual movement or product wear is detected and consult a practitioner.

In cases of damage, malfunction, changes in alignment, or functional loss, the user should immediately cease using the device, and return it to their practitioner.

Improper handling and/or adjustment of the device may lead to malfunction.

The device is battery-operated. In cases of power failure or depletion, the device may malfunction and subject the user to risk of falling.

Listen to audio warning signals. Attention to audio signals emitted by the device is critical for safe device operation. Take the measures provided in this document, specifically in the section “*Additional Device Feedback*”.

The device has been designed and tested based on single patient usage. **DO NOT** use with multiple patients.

Adhere to local ordinance or regulations for proper disposal of the device, accessories, and packaging.

### **Precautions**

The device does not support use in high impact activities, sports, excessive loading, and heavy-duty use. Do not wear boots or shoes restricting movement of the ankle area (i. e. due to height or stiffness of material). Ankle motion restriction will impair device function.

Avoid impact to the interface panel on front and back of device.

Extended exposure to temperatures below 0°C (32°F) may prevent the ankle from moving as expected.

Do not sterilize the device.

Do not expose the device to environments containing flammable anesthetic mixtures.

Avoid excessive mechanical shocks and vibrations.

Beware fingers can be caught and pinched the back of the device when in operation.

Do not insert wedges between the carbon blades as this may cause premature failure of the carbon blades.

### ***Weatherproof***

A weatherproof device can be used in a wet and/or humid environment but cannot be submerged. Freshwater splashing against the enclosure from any direction will have no harmful effect.

Dry thoroughly after contact with freshwater. Freshwater: includes tap water. Excludes salt and chlorinated water.

Weatherproof components must be used with the device for a weatherproof system.

### **Precautions:**

Do not submerge device. Do not expose device to salt water, chlorinated water, or liquids other than fresh water.

Do not expose device to pressurized water jets.

Do not use device as a shower prosthesis.

Immediately return the device to Össur if exposed to any of the conditions above.

**Note:** Full weatherproof performance is only achieved by using Össur titanium adapters. No other adapters or screws can be used. Adapter installation and set-screw tightening must be performed according to instructions provided in respective accompanying documents.

---

## SUPPORTED ACTIVITIES

---

The device automatically adjusts its behavior to activities including standing, walking, ramp ascent/descent, and stair ascent/descent.

**Note:** The device only moves in swing phase when not under load. Excessive loading in swing phase, caused by short swing phase duration (lower than 0.4 sec), heavy shoes, interrupted swing cycle or foot contact with the ground, can result in an incomplete movement.

The device is equipped with a safety mechanism to cease motion in case it is unable to move to the desired position.

### **Level Walking**

- When walking on level ground above 1.5 km/h (0.9 mph), the device will raise the forefoot in swing phase to provide more ground clearance to reduce risk of stumbling.
- This motion occurs only after the first two prosthetic steps and in each prosthetic step thereafter.

**Note:** Circumduction or any abnormal side movement may interrupt initiation of dorsiflexion during swing phase.

### **Incline – Ascent**

**Warning:** Inclines of more than 15° (ascending and descending) may trigger the stair response function, resulting in unexpected device alignment.

- The device will gradually raise forefoot to accommodate the angle of the incline.
- During swing phase, the device will raise forefoot to provide additional ground clearance, lowering risk of stumbling.

- The user is not required to lead with a specific limb during ascent.

### **Incline - Descent**

- The device will gradually lower forefoot during descent to accommodate the angle of the incline.
- The user is not required to lead with a specific limb during descent.

### **Stair Ascent**

- The device will raise the forefoot during stair ascent.
- Leading with the sound side upstairs will provide stair adaptation after first prosthetic step.
- Leading with the prosthetic side upstairs will provide stair adaptation after second prosthetic step.
- On completion of stair ascent and stepping onto level ground, device will adapt to level ground after one stride.

**Note:** User should consult practitioner if abnormal socket pressure is felt during stair ambulation. Stair response can be altered using the Össur Logic mobile app (see “*Accompanying Software*” section) to minimize socket pressures.

### **Stair Descent**

- User must be aware that device will raise forefoot during stair descent.
- User should place entire foot on step. Doing so allows more stable descent compared to non-microprocessor-controlled ankle systems.
- When descending, stair adaptation occurs after second prosthetic step.
- When reaching bottom of stairs, device will adapt to level ground after one stride.

**Note:** User should consult practitioner if abnormal socket pressure is felt during stair ambulation. Stair response can be altered using the Össur Logic mobile app (see “*Accompanying Software*” section) to minimize socket pressures.

### **Sitting**

- Use the Össur Logic mobile app to configure the Relax Mode function (see “*Accompanying Software*” section).
- Relax Mode fully lowers forefoot, improving body symmetry to sound limb when sitting.
- Tap foot on ground in Relax Mode to revert the device to a neutral state.
- Relax Mode is OFF by default.

### **Chair-Exit**

- Chair Exit Mode function can be configured using the Össur Logic mobile app (see “*Accompanying Software*” section).
- With Chair Exit Mode enabled, exiting Relax Mode will cause the device to move automatically into Chair Exit Mode. The device will raise the forefoot, allowing the user to bring the leg further back, for more controlled transition from sitting to standing. The device will return to a neutral state after the first step.
- With Chair Exit Mode enabled, flexing the knee and placing the foot beneath the chair will allow ankle angle adjustment to take place.
- Tap foot on ground in Relax or Chair Exit Mode to revert it to a neutral state.

**Note:** With Chair Exit Mode enabled, but Relax Mode disabled, user must first extend the leg to enter Chair Exit Mode.

### **Standing**

- The device operates like a normal energy-storing and energy-returning foot while standing.
- Flexing the hip and placing the device in front of the sound foot will activate Relax Mode, lowering the forefoot when standing.
- Extending hip from the flexed position will revert ankle to a neutral state.

### **Kneeling**

- Kneeling will activate Relax Mode. Both knee and shin will rest in a less stressful position with reduced knee flexion.

## Driving

- Before driving, turn the device OFF or set to Standby mode in the Össur Logic mobile app (see “*Accompanying Software*” section).

**Warning:** Take caution when using the device in a motor vehicle. Unintended ankle movement may interfere with brake and/ or accelerator pedal function.

---

## POWER

---

The device contains a lithium-ion battery of 2900 mAh/20.88Wh.

### Warnings:

- Temperatures below -10°C (14°F) or over 45°C (113°F) will affect battery performance. This may result in battery warnings or prevent the device from moving as expected.
- Using a damaged power supply can lead to electric shock.
- In case of unexpected power loss, the device will stop moving and functioning like a mechanical foot, resulting in incorrect ankle alignment.

**Note:** Only use accompanying power supply with the device; failure to do so will invalidate the device warranty. The device comes with the following power supply:

<b>Manufacturer:</b>	FRIWO
<b>Model No.:</b>	FW8030M/24
<b>Input:</b>	100–240 VAC, 50–60 Hz, 0.6–0.3 A
<b>Output:</b>	24 VDC, 1.25 A

## CHARGING

### Warnings:

- Do not use the device when charging the battery.
- Disconnect the device from the power supply before donning.

The device can be charged turned either ON or OFF.

- Plug power supply DC connector into the charge connector at the back of the device and plug charger into the wall.
- The LED Indicator Ring on front of the device will display the device is charging, and show when it is fully charged.
- LED Indicator Ring flashing green at last clockwise position shows charging is in progress.
- LED Indicator Ring displays battery charge level through increasing number of solid green lights (see “*Battery Status*” section).
- Charging the device to 90% takes approximately 3.5 hours. A full charge (100%) takes approximately 4 hours.

A fully charged device can be used for 18 to 36 hours, depending on use.

The device cannot be overcharged.

Charging the battery from fully depleted to fully charged state consumes 57.9 VAh at 230 VAC / 50 Hz & 49.3 VAh at 110 VAC / 60 Hz.

**Note:**

- After disconnecting the power supply DC connector, the device will turn ON and be ready for operation.
- Charge device daily for optimal performance.
- disconnect power supply from power outlet to save energy when not in use.

**BATTERY STATUS**

**Warnings:**

- Device returns to last pre-set Ankle Alignment position and stops moving when battery charge is critically low. A warning signal will be emitted.
- Walk with caution if an alarm signal is heard
- Press power button once to display battery charge. 
- LED Indicator Ring will light up clockwise, showing the battery charge level.
- The number of lights on the LED Indicator Ring shows the battery charge level, as in the table below:

LED Indicator Ring Appearance	Battery Charge Level
	0% to 15% battery charge level Blinks yellow and emits 3 beeps
	16% to 29% battery charge level
	30% to 39% battery charge level
	40% to 49% battery charge level
	50% to 59% battery charge level
	60% to 69% battery charge level
	70% to 79% battery charge level
	80% to 89% battery charge level
	90% to 94% battery charge level
	Fully charged 95% to 100% battery charge level

## TURNING ON & OFF

### To turn the device ON:

- Hold the power button for longer than 2 second.
- Device emits 3 ascending tones.
- Start-up sequence will initiate, shown on the LED Indicator Ring by a green pattern rotating clockwise from the top.

- A single green light will then flash at top of LED Indicator Ring, showing the device is functional and ready to use.

### To turn the device OFF:

- Hold down power button for longer than 2 seconds.
- Device emits 3 descending tones and turns OFF.

---

## ANKLE ALIGNMENT

---

The device can be adjusted to variable heel heights, optimizing ankle alignment. This ensures correct alignment of the prosthesis, essential for proper function.

### NOTE:

- Heel height adjustment range is 0–50 mm (2") for a size 27, category 5 device.
- Ensure selected footwear allows free ankle motion.
- Remain aware of audio feedback.

### ANKLE ALIGNMENT WITHOUT WEARING PROSTHESIS

Ankle Alignment is performed without the user wearing the device.

1. Fit the device with suitable footwear and ensure free ankle motion.
2. Place the device on a level surface; make sure device is fully-charged device turned ON.
3. Initiate Ankle Alignment by pressing power button quickly twice successively.
4. Initiation confirmed: the device will beep once, and a single white LED will light up on the LED Indicator Ring.

5. Ankle Alignment in progress: LED Indicator Ring will light up white clockwise.

6. The actuator will move prosthesis into a vertical position.

*Ankle Alignment completed:* if successful, the device will emit an ascending audio tone, and LED Indicator Ring will flash white.

*Ankle Alignment failed:* if unsuccessful, LED Indicator Ring will flash yellow, showing process needs to be repeated.

**Note:** If heel height exceeds 50 mm (2") the device will emit an audio signal and not adjust to height. The user must then change shoes to a heel height under 50 mm (2") and repeat the adjustment.

### **ANKLE ALIGNMENT WEARING PROSTHESIS**

After changing shoes, the user can perform Ankle Alignment while wearing the prosthesis.

**Note:** The user should always perform User Ankle Alignment when changing shoes or going barefoot.

1. The user should be seated.
2. Ensure device is positioned flat, with both the heel and forefoot area on the ground.
3. Ensure the device is fully-charged and turned ON.
4. Initiate Ankle Alignment by quickly pressing the power button twice successively.
5. Initiation confirmed: the device will beep once, and a single white LED will light up.

6. Wait for the 2nd beep and lift the device from the ground to allow for ankle adaptation.
7. Ankle Alignment will initiate as described above from step (5).

**Users are asked to skip the "Initial Set-Up" section. These instructions are for the practitioner only. The user should not attempt these adjustments.**

---

## INITIAL SET-UP

---

The device set-up consists four steps, described in the following sections. Perform the steps in the following order:

1. Bench Alignment including Ankle Alignment
2. Dynamic Alignment
3. Auto-Adjustment
4. Advanced Set-Up

After turning the device ON, two orange lights on the LED Indicator Ring will show initial set-up is required.

One orange light on the right side of the LED Indicator Ring shows Auto-Adjustment is required.

One orange light on the left side of the LED Indicator Ring shows Ankle Alignment is required.

After successfully performing these activities, the orange lights will turn off and disappear from the LED Indicator Ring.

### **Warnings:**

- The device will not function correctly if Auto-Adjustment and Ankle Alignment are not carried out properly.
- The user must use the support of a handrail when descending stairs or ramps during initial fitting.

### **BENCH ALIGNMENT (Figure 3)**

**Precaution:** Always use the device with the recommended foot cover and Spectra Sock to ensure the device works properly. Don the Spectra Sock in a way that keeps it from interfering with ankle movement.

1. Fit the device with the selected foot cover, Spectra Sock, and shoe.
2. Turn the device ON.
3. Adapt preferred shoes by performing Ankle Alignment (see above).
4. Establish the correct height of the prosthetic system using the appropriate components.
5. Introduce suitable socket angles; flexion/extension, and abduction/adduction using the adapter set-screws.
6. Divide the foot cover length into 3 equal portions (as set out below).
7. The alignment reference line should fall at the junction of the posterior and the middle third portions (as displayed below).

### **DYNAMIC ALIGNMENT**

The heel stores energy after initial contact, slowly releasing it at mid-stance.

Heel-to-toe action can be influenced by anterior/posterior positioning of the device.

**Note:** Align the device to minimize medial-lateral motion during swing phase. Circumduction or any abnormal swing phase side motion may affect the initiation of dorsiflexion during swing phase.

**Note:** Do not use the Manual Ankle Alignment in the Össur Logic mobile app (see “*Accompanying Software*” section) to adjust flexion/extension. Use the adapter set-screws to make adjustments.

### **HEEL STIFFNESS**

#### **Heel too soft**

##### *Symptoms*

Device will come to a flat position too early (user feels he/she is sinking into a hole)

Extra energy is required to climb up over the forefoot

Toe will feel too stiff

Knee may hyperextend

##### *Actions*

Shift socket anteriorly

Consider changing device blades to higher category

#### **Heel too hard**

##### *Symptoms*

Rapid heel-to-toe movement

User has poor control of prosthesis

Toe will feel too stiff (minimal energy return)

Knee may become unstable

##### *Actions*

Shift socket posteriorly

Consider changing device blades to lower category

### **AUTO-ADJUSTMENT**

Perform Auto-Adjustment during initial fitting to calibrate the device to the user’s individual gait.

During Auto-Adjustment, the device recognizes the user’s specific gait parameters, allowing for accurate and consistent terrain detection. To complete Auto-Adjustment, the user must take 15 adequate level ground walking strides, preferably in a consecutive and uninterrupted manner.

## NOTE:

- Make sure the user has access to a level ground walkway long enough to complete at least 15 adequate strides without interruption, at a consistent self-selected speed. If Auto-Adjustment is not successful after 15 strides, continue until successful.
- Make sure the walkway surface is completely flat. It is recommended to perform Auto-Adjustment on an indoor surface. The user should walk in a straight line, while avoiding walking in a semi-circle or turning.
- If the walkway is not long enough for 15 strides, instruct the user to stop, turn around, and continue walking in a straight line. This will increase the number of strides needed.
- the device will not raise the forefoot during the Auto-Adjustment procedure.

1. Ask the user to stand still.
2. Enable Auto-Adjustment using the Össur Logic mobile app (see “*Accompanying Software*” section).
3. A single white light on the LED Indicator Ring will light up to confirm Auto-Adjustment has been initiated.

4. Ask the user to walk on level ground at a self-selected walking speed. The LED Indicator Ring will light up white in a clockwise direction, showing successful step registration.

5. The LED Indicator Ring will light white and will turn off for an instant, before showing the outcome of the Auto-Adjustment process.

*Auto-Adjustment completed:* if successful, an ascending 3-tone audio sequence will be emitted, and the LED Indicator Ring will flash white.

*Auto-Adjustment failed:* if unsuccessful, the LED Indicator Ring will flash yellow, indicating the process needs to be repeated. Please refer to the “*Troubleshooting*” subsection below.

## Troubleshooting

If Auto-Adjustment is not successful:

- Make sure alignment is correct.
- Make sure there are no obvious gait deviations such as significant circumduction or rotation, as these may prevent successful Auto-Adjustment.
- Make sure the user is walking at a consistent speed.
- Make sure the user is not walking in a semi-circle or making turns.

**Note:** The device will not operate normally until successful Auto-Adjustment is confirmed.

**Note:** Auto-Adjustment must be repeated following changes in alignment.

## ADVANCED SET-UP

Stair and ramp adaptation can be fine-tuned using the Össur Logic mobile app (see “*Accompanying Software*” section). Refer to the in-app instructions for information on how to adjust these advanced settings and other user preferences.

---

## ADDITIONAL DEVICE FEEDBACK

---

Indicator Display	Audio Feedback	Device Status
	None	Normal Single light flashing green. Device is ON and operational.
	3 beeps	Low Battery Single light flashing yellow. Battery is below 15% charge, functionality has been reduced.
	None	Error Single light flashing red. Device has registered an error and needs to be restarted. If error persists, device must be returned to Össur for service.
	10 beeps	Critical Error LED Indicator Ring flashing red. Device is faulty and non-operational. Device must be returned to Össur for service.
	2 ascending tones	Connecting Blue light rotating clockwise. Establishing connection to the Össur Logic mobile app.
No visual feedback	2 descending tones	Disconnecting Disconnecting from the Össur Logic mobile app.
	None	Connected Single light flashing blue. Device is ON, operational, and connected to the Össur Logic mobile app.

---

## ACCOMPANYING SOFTWARE

---

The Össur Logic mobile app can be downloaded on the Apple App Store and be accessed by both the user and practitioner to connect to the device and access functions relevant to fitting and use.

---

## MAINTENANCE

---

The user should periodically inspect the device for hardware integrity and consult their practitioner if any degradation is observed.

### CLEANING

Wipe the device with a soft cloth. Moisten cloth with mild soap or 30% isopropyl alcohol solution.

Debris caught between the carbon foot blades can result in noise or premature wear. Remove debris or dust without using compressed air.

**Warning:** Do not attempt to repair the device. Send the device to Össur for repairs or service.

**Warning:** Never dismantle or remove device covers.


## COMPLIANCE

The device has been tested by an independent accredited laboratory and found to comply with:  
EN60601-1:2006/A1:2013  
EN60601-1-2:2015  
IEC62133 2nd ed.

### ISO10328 Standard

The device has been tested according to ISO 10328 standard to 2 million load cycles. The expected life duration under normal usage is 2 million steps. Depending on user activity, this corresponds to a duration of use of 1 to 3 years.

### Test Levels

ISO 10328 - Label			
Category	Weight (kg)	Weight (lbs)	Label Text
1	52	115	ISO 10328-P3-52-kg
2	59	130	ISO 10328-P3-59-kg
3	68	150	ISO 10328-P3-68-kg
4	77	170	ISO 10328-P4-77-kg
5	88	194	ISO 10328-P4-88-kg
6	100	220	ISO 10328-P5-100-kg
7	116	256	ISO 10328-P5-116-kg
8	125	276	ISO 10328-P5-125-kg

## ENVIRONMENTAL CONDITIONS

	Use	Charging	Shipping	Extended storage
Temperature	-10 °C to 45 °C (14 °F to 113 °F)	10 °C to 40 °C	-25 °C to 70 °C (-13 °F to 158 °F)	0 °C to 45 °C (32 °F to 113 °F)
Humidity	0% - 100% (including condensing)	0% to 100% relative humidity (non- condensing)	10% - 100% relative humidity (including condensing)	10% - 90% relative humidity
Atmospheric pressure	700hPa to 1060hPa	700hPa to 1060hPa	700hPa to 1060hPa	700hPa to 1060hPa

---

## ELECTROMAGNETIC COMPATIBILITY

---

### Warnings

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

### Precaution

The device may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

The PROPRIO FOOT® is suitable for use in any environment except where immersion in water or any other fluid is possible, or where exposure to highly electrical and/or magnetic fields can occur (e. g. electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Electromagnetic Emissions Compliance Information		
Emissions Test	Standard	Compliance Level
Conducted and Radiated RF Emissions	CISPR 11	Group 1 – Class B
Harmonics Emissions	IEC 61000-3-2	Class A
Voltage Fluctuations / Flicker Emissions	IEC61000-3-3	Complies

Electromagnetic Immunity Compliance Information		
Immunity Test	Standard	Compliance Level
Electrostatic Discharge	IEC61000-4-2	± 8 kV contact ± 15 kV air
Electrical Fast Transient/Burst	IEC61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines
Surge	IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage Dips	IEC61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage Interruptions	IEC61000-4-11	0 % UT; 250/300 cycle
Power Frequency (50/60 Hz) Magnetic Field	IEC61000-4-8	30 A/m
Conducted Disturbances Induced by RF Fields	IEC 61000-4-6	3 V 0.15 MHz – 80 MHz 6 V in ISM and amateur bands between 0.15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF EM Fields	IEC 61000-4-3	10 V/m 80 MHz – 2.7 GHz 80 % AM at 1 kHz

RF Wireless Communications Compliance Information				
Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380 - 390	TETRA 400	Pulse modulation 18 Hz	27
450	430 - 470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704 - 787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800 - 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700 - 1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100 - 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

## BLUETOOTH MODULES REGULATED INFORMATION

This device contains the following radio frequency transmitters:

Model	Regulatory Certificates	Type and Frequency Characteristics	Effective Radiated Power
WiFi Module Model: ATWINC1500-MR210PB	FCC Contains FCC ID: 2ADHKATWINC1500  Canada Contains transmitter module IC: 20266-WINC1500PB	2412–2462 MHz (5 MHz Interval, 11 Channels)	
Bluetooth 4.2 Module Model: BMD-300	FCC Contains FCC ID: 2AA9B04  Canada Contains IC: 12208A-04  Japan Contains transmitter with certificate number   R210-106799	FID 2402–2480 MHz	0.35 mW/MHz

### USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

### Caution: Exposure to Radio Frequency Radiation.

This device must not be co-located or operating in conjunction with any other antenna or transmitter.

### Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

1. this device may not cause interference, and
2. this device must accept any interference, including interference that may cause undesired operation of this device.

**Caution: Exposure to Radio Frequency Radiation.**

The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website [www.hc-sc.gc.ca/rpb](http://www.hc-sc.gc.ca/rpb).

**DISPOSAL**

All components of the product and packaging should be disposed of in accordance with respective national environmental regulations. Users should contact their local governmental office for information on how these items can be disposed of in an environmentally sound manner.

**LIABILITY**

The manufacturer recommends using the device only under the specified conditions and for the intended purposes. The device must be maintained according to the instructions for use. The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

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**EN – Caution:** Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards is achieved only when Össur products are used with other recommended Össur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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