Instructions for Use

PYLONS AND TUBES

EN		Instructions for Use	3
DE		Gebrauchsanweisung	6
FR		Notice d'utilisation	9
ES		Instrucciones para el uso	12
IT		Istruzioni per l'uso	15
NO		Bruksanvisning	18
DA		Brugsanvisning	21
SV		Bruksanvisning	24
EL		Οδηγίες Χρήσης	27
FI		Käyttöohjeet	30
NL		Gebruiksaanwijzing	33
PT		Instruções de Utilização	36
PL		Instrukcja użytkowania	39
TR		Kullanım Talimatları	42
RU		Инструкция по использованию	45
JA		取扱説明書	48
ZH		中文说明书	51
KO		사용 설명서	54
AR		تعليمات الاستخدام	57

ENGLISH

Medical Device

INTENDED USE

The device is intended as a connection component of a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

For more detailed device information, see **Table 1**.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

For weight limit and impact level, see **Table 2**.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Only use Össur replacement parts for the device for ensured durability.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Cutting the tube:

1. Clamp the tube with a cutting jig or put into a tube cutter.

Caution: Clamping too hard may damage the device.

2. Cut the tube straight and perpendicular.

Warning: Do not cut adapter shorter than the minimum cutting length, see **Table 1**.

3. Remove all burrs and fully insert the tube into adapter.

Only applicable for female pylons:

Warning: For secure fastening use medium strength threadlocker and torque twice to 15 Nm.

Warning: When extreme alignment is necessary, use extended set screws to achieve secure fastening.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

See which environmental conditions apply to the device in **Table 1**.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour. It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

SPECIFICATIONS

Table 1 - Device Specifications

Part #	Description	Product Length	Diameter	Material	Environmental Conditions
AJ-712010	Tube Junior	250 mm	22 mm	Aluminum	Weatherproof
AJ-712020					Waterproof
A-712010	Tube	400 mm	30 mm		Weatherproof
A-712020					
A-742010	Female	250 mm			
A-742020	Pylon AL	400 mm			
A-748010	Female	250 mm		Aluminum, Titanium	Waterproof
A-748020	Pylon AL TI	400 mm			
A-748030	Female	250 mm			
A-748040	Pylon TI WP	400 mm			

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

Table 2 - Usage Specifications

Part #	Impact Level	Weight Limit	ISO 10328 - Label Text
AJ-712010	The device is intended for children and adolescent patients.	45 kg / 100 lbs	ISO10328 – P3 – 45 kg
AJ-712020			
A-712010	The device is for moderate impact use, e.g., walking.	100 kg / 220 lbs	ISO10328 – P5 – 100 kg
A-712020			
A-742010			
A-742020			
A-748010			
A-748020			
A-748030	The device is for high impact use, e.g., walking and occasional running.	166 kg / 366 lbs	ISO 10328 – P7 – 166 kg
A-748040			

ISO 10328 - “P” - “m”kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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