





Instructions for Use

---

REBOUND<sup>®</sup> AIR WALKER LT  
REBOUND<sup>®</sup> AIR WALKER

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	15
NO	Bruksanvisning	17
DA	Brugsanvisning	20
SV	Bruksanvisning	22
EL	Οδηγίες Χρήσης	25
FI	Käyttöohjeet	28
NL	Gebruiksaanwijzing	30
PT	Instruções de Utilização	33
PL	Instrukcja użytkowania	36
HR	Upute za uporabu	39
JA	取扱説明書	41
ZH	中文说明书	44
KO	사용 설명서	47
RO	Indicații de utilizare	49
BG	Указания за употреба	52

# ENGLISH

---

## SYMBOLS

Medical Device

## INTENDED USE

The device is intended for immobilization of the foot and ankle

### *Indications for use*

#### **Rebound Air Walker:**

Foot and ankle conditions that may benefit from motion restriction, such as:

- Grade 2 and 3 ankle sprains
- Stable fractures
- Soft tissue injuries
- Post-operative, trauma or rehab use

If used with the Achilles Wedge Kit, conservative and post-operative treatment of Achilles tendon ruptures.

#### **Rebound Air Walker Low Top:**

Foot and ankle conditions that may benefit from motion restriction, such as:

- Stable foot fractures
- Soft tissue injuries
- Post-operative, trauma or rehab use

No known contraindications.

### **Warnings and Cautions:**

Do not remove the device unless instructed by a healthcare professional.

- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- This device should not be applied on patients who are not able to communicate discomfort.
- Patient should use caution when walking on slippery or wet surfaces to avoid injury.
- Ensure the patient performs circulation checks. If a loss of circulation is felt, or if the walker feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.
- Care should be taken not to overtighten the device.
- Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.
- Patient maximum weight limit is 300 lbs / 136 kg (200 lbs / 90 kg for the XS version).

Air version only:

- Ensure not to overinflate the air chamber.
- Over inflation of air chamber may cause skin compromise in patients with diabetic neuropathy. Deflate air chamber to reduce compression with patients reporting discomfort or showing skin discoloration.

- Ensure patient performs circulation checks. If a loss of circulation is felt, or if the device feels too tight, deflate the air chamber to reduce compression to a comfortable level. If discomfort continues, the patient should contact their healthcare professional immediately.
- In high altitude conditions, the air chamber will expand past its optimal level. Deflate air chamber to reduce compression to a comfortable level.

## GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## USAGE

### *Device Application*

1. Open the liner by first unfastening straps and removing the front shell.  
**Note:** Do not remove the air bulb from the liner.
2. While seated, place the leg into the device. The lower leg should be kept at a 90° angle with the foot and ensure that no gap is in between the liner heel and the back of the device.
3. Wrap liner securely around the foot and the leg (**Fig. 1**).
4. Place front shell over ankle making sure that the air bulb is aligned with the opening (**Fig. 2**).
5. In the following order, fasten the ankle strap (only applies to the high-top version), the foot strap, and then the calf strap. Tighten straps until snug and comfortable. Strap adjustments may be necessary with changes in lower limb volume
6. Inflate air chambers by pressing “+” on the air bulb, and deflate by pressing “-“. Inflate or deflate until snug (**Fig. 3**).  
**Caution:** Do not overinflate!

### *Device Adjustments*

Fitting adjustments can be made to the device if desired:

- The toe guard can be trimmed off (**Fig. 4, 5**) for a better fit or if the patient has an apical lesion or wound.
- The device circumference around calf can be widened by cutting Flex Edge® (**Fig. 6, 7**) to accommodate patients with a large calf.  
**Note:** Ensure to remove the calf strap before cutting the Flex Edge®.
- If the front shell is too long it can be trimmed along a trim line if present. Remove the front shell from the device, turn it over and cut along the trim line with a strong pair of scissors. Be careful not to leave any sharp edges.

### *Device Removal*

Before removing the device, deflate air chambers by pressing “-“.

Undo the straps in any order and attach the straps back to themselves.

This will make it easier when fitting the device next time and will prolong the life of the straps. Open the Liner to lift the foot and lower leg out of the device.

## **Accessories and Replacement Parts**

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## **Cleaning and care**

### **Washing Instructions**

Device plastic can be washed with a damp cloth and mild soap. Remove liner from the device before washing.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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