





Instructions for Use







REBOUND® HIP



		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	11
ES	Instrucciones para el uso	14
IT	Istruzioni per l'uso	17
NO	Bruksanvisning	20
DA	Brugsanvisning	23
SV	Bruksanvisning	25
EL	Οδηγίες Χρήσης	28
FI	Käyttöohjeet	31
NL	Gebruiksaanwijzing	34
PT	Instruções de Utilização	37
PL	Instrukcja użytkowania	40
CS	Návod k použití	42
TR	Kullanım Talimatları	45
RU	Инструкция по использованию	48
JA	取扱説明書	51
ZH	中文说明书	54
KO	사용 설명서	56
RO	Indicații de utilizare	59



11


12

13

14

SYMBOLS

Medical Device

INTENDED USE

The device is intended for motion restriction of the hip

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Hip conditions that could benefit from motion restriction such as:

- Postoperative after (total) hip arthroplasty
- Postoperative after revision of (total) hip arthroplasty
- Postoperative hip arthroscopy
- Post-op following FAI arthroscopic surgery
- Post-op following Gluteus medius repair surgery
- Post-op following hamstring repair surgery
- Inoperable hip abnormalities
- Preventively when there is risk of hip dislocation
- Conventional treatment after dislocation/reduction

No known contraindications.

Warnings and Cautions:

- Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.
- Hip reconstructions should be stable prior to treatment with the medical device.
- Patients are advised not to operate a vehicle during treatment.
- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- Care should be taken not to overtighten the device.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is a supportive brace only and is not intended or guaranteed to prevent injury.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

1. Locate the patient's Greater Trochanter, align the center of the hinge at the hip, and have your patient hold the hinge in place (**Fig. 2**). This allows them to feel the appropriate location for the Hip Plate (A).

2. If needed, use bending irons to adapt the strut for necessary abduction/adduction angles or adjust to the patient's leg shape (for units without abduction hinge).
Note: For units with Abduction Hinge (B): Using the included Allen Wrench, adjust the hinge from -10 to 20 degrees as desired (**Fig. 3**). (Device ships at 0 degrees).
3. Adjust the length of the strut such that the distal Thigh Strap (C) will not interfere with the popliteal region of the knee. Use the Push Buttons (D) to adjust the length and click them into place (**Fig. 4**). Position the lower thigh cuff in such a way that the bottom edge is located approximately 5 cm (2 in) above the top of the patella.
Note: On the shortest setting both Push Buttons may need to be engaged to allow freedom of movement.
4. Secure the Thigh Strap into the keyhole receiver on the distal Thigh Plate (E) (**Fig. 5**). If needed, adjust and/or trim the Thigh Strap to the appropriate length (**Fig. 6**).
5. Ensure the Belt Compression System is stretched to its full width before fitting. Position the belt with the pulley system at the patient's back and wrap the belt arms around the patient, covering the thigh pad on the device, and securing in front. Make any adjustments necessary to keep the device in the correct position over the hip and aligned with the Greater Trochanter. Place Position marker on belt to identify optimal alignment of belt and frame.
6. Wrap the Adjustable Belt (F) around the waist and insert the Össur Fit Tool (G) (**Fig. 7**). Adjust the length of the belt arms to position the tool centrally. Secure the Belt Ends to desired sizing (**Fig. 8**).
Note: If the patient's waist exceeds 127 cm (50 in), an optional Belt Extension can be added to the belt arms to reach a maximum belt circumference of 178 cm (70 in).
7. Remove Össur Fit Tool from belt (**Fig. 9**).
Note: If an Össur Fit Tool is not available, lay the belt flat with Overlap Closures removed and compression system loosened. Fold each belt arm equally to achieve the desired waist circumference and secure the Belt Ends.
8. Place Belt Ends into the Overlap Closures (**Fig. 10**). Once Belt Ends are secured, the Overlap Closures will serve as Belt Ends.
Note: To obtain the lowest profile fit, the Belt Ends can be trimmed if later adjustments are not required.
9. Close belt securing the Overlap Closures (**Fig. 11**). Lay down the left-hand belt-arm first (1), followed by the right (2).
10. Slide thumbs through the holes in the two Compression System Handles and pull until the belt is at the appropriate tightness (**Fig. 12, 13**). If the cords are too long, remove device from patient and reset the pulley system by pulling outward on belt.

Device Adjustments

Range of Motion (ROM) Adjustment

Using the quick release pegs, adjust the Flexion/Extension Stops (H) to the desired range of motion and secure the pegs with a quarter turn twist until they lock back into position (**Fig. 14**).

For a left setting of the device:

- 90 degrees Flexion / 0 degrees Extension: place quick release pegs in hole marked "L" and hole marked "1" and "2" (as shipped).
- 45° Flexion: place quick release peg in holes L, 2 and 4.

For a right setting of the device:

- 90 degrees Flexion / 0 degrees Extension: place quick release pegs in hole marked "R" and hole marked "3" and "4" (as shipped).
- 45° Flexion place quick release peg in hole marked R, 1 and 3.

Orientation Adjustment

Note: While the device is available as a left or right orientation, it can be converted to the opposite orientation by reversing the Left (L) or Right (R) pegs at the base of the Hip Plate, and readjusting the Flexion/ Extension Stops for appropriate range of motion control. The device will be set at 90 degrees flexion and 0 degrees extension out of the box.

Device Removal

To remove the device,

1. Detach Compression System Handles slowly from the belt to loosen and re-attach them to the sides of the Adjustable Belt.
2. Detach the Thigh Strap.
3. Detach the Overlap Closures and remove the device.

Note: To ensure a proper fit, ensure the Belt Compression System is stretched to its full width before re-applying the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Device plastic can be washed with a damp cloth and mild soap.

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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