



Instructions for Use

RHEO KNEE[®]

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	14
FR	Notice d'utilisation	25
ES	Instrucciones para el uso	35
IT	Istruzioni per l'uso	46
NO	Bruksanvisning	57
DA	Brugsanvisning	67
SV	Bruksanvisning	77
EL	Οδηγίες Χρήσης	87
FI	Käyttöohjeet	98
NL	Gebruiksaanwijzing	108
PT	Instruções de Utilização	119
PL	Instrukcja użytkowania	130
TR	Kullanım Talimatları	141
RU	Инструкция по использованию	151
JA	取扱説明書	163
ZH	中文说明书	173
KO	사용 설명서	183

SYMBOLS

Medical Device

Caution Symbol

Consult instructions for use

Type B Applied Part

ON / OFF

Direct current

Serial number

Device contains electronic components and/or batteries that should not be disposed of in regular waste

IP34

Protected against harmful effects of objects over 2.5 mm and against the effects of splashing water

DESCRIPTION

The device is a microprocessor-controlled prosthetic knee. It provides resistance to movement by means of a magneto-rheological actuator.

1. Proximal pyramid adapter;
2. Magnetorheologic actuator;
3. Mechanical lock;
4. Power button;
5. Status indicator;
6. Charging port;
7. Battery;
8. Distal pyramid adapter.

Device labels (**Fig. 2**) can be found at the rear of the device, on the power supply and on the internal face of the battery pack.

This document is for patients (intended device operator) and healthcare professionals.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional. The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Unilateral transfemoral and knee disarticulation amputation;
- Bilateral transfemoral and knee disarticulation amputation.

No known contraindications.

The device is for low to moderate impact use, e.g., walking.

The weight limit for the device is 136kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury. The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Ensure the patient is able to recognize warning signals for power loss.

Warning: If device emits a warning signal, stop using device and refer to **Table 1** for more information. Note that certain errors can cause the device to shut down.

Warning: Avoid placing hands or fingers near moving joints.

Warning: A case of unexpected power failure or battery depletion during use will cause the device to stop functioning normally and may subject the patient to the risk of falling.

Note: In case of power failure, malfunction, or warnings it is recommended to engage the mechanical lock. When locked, the device can support full weight-bearing.

Warning: The device's misinterpretation of patient activity increases the risk of falling.

Note: High intensity uses or operation in high ambient temperature will cause the device actuator to heat up and become hot to the touch.

Note: Improper handling and/or adjustment of the device can cause malfunction. Patient should avoid:

- Direct impact to the device;
- Excessive shocks or vibrations;
- High-impact activities, sports, excessive loading, and heavy-duty use.

Warning: Not suitable for use in the presence of flammable anesthetic mixture with air, oxygen, or nitrous oxide.

Note: Do not attempt to modify device in any way.

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.
The device is for single patient use.

REQUIRED COMPONENTS

Össur Logic App

For the initial fitting, the healthcare professional must use the Össur Logic app. It is available on Apple App Store. See Össur Logic on the App Store for mobile device compatibility.

The healthcare professional can use the app to adjust the knee to the user's physiology (weight, strength, etc.), gait style, and personal preferences.

The patient can use the Össur Logic app to check the status of the device, monitor the battery charge level, count steps, access training programs, and view instructions for use.

See instructions on how to connect to the app in **Usage** section.

DEVICE SELECTION

See Össur catalog for a list of recommended components.

Use Össur titanium adapters to achieve weatherproof performance.

Use screws with anti-corrosive coating.

Note: Adapter installation and set-screw tightening shall be performed according to instructions provided in their respective accompanying documents.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through knee axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through knee axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i. e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Caution: At maximum flexion, ensure that a minimum distance of 3 mm (1/8") remains between the device and the socket (**Fig. 4**). If contact between device and socket cannot be avoided due to socket volume, make sure that the point of contact is on the device frame directly below the user interface panel. Keep the contact area flat and cushioned for pressure distribution.

Static Alignment

Turn the device ON.

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

1. Keep device powered ON.
2. Ensure the patient is familiar with how the device functions.
3. Tell the patient to walk between parallel bars with usual gait. This lets the patient get used to the behavior of the device.
4. Adjust dynamic alignment as required.
 - Make sure that the motion in swing and stance phases stays in the line of progression.
 - Make sure to achieve step-length symmetry.
5. The Össur Logic App must be used for setting device parameters during dynamic alignment.

Note: Use the “Initial Setup” function in the Össur Logic app to adjust the device for basic use. Use the “Advanced Settings” function to fine tune the device.

If stability modifications are desired after patient assessment, the reference line can be shifted up to 5 mm (3/16") anterior or 5 mm (3/16") posterior to the knee axis.

Caution: As the alignment reference line is positioned more posterior to the knee center, the patient will experience more stance flexion during loading response. Additional voluntary control will be required to maintain knee stability when the device is powered OFF.

Caution: Errors that can occur during alignment and adjustment may cause malfunction of the device and lead to a fall.

After Alignment

After alignment, the healthcare professional should:

- Have the patient load the heel and allow the knee to flex while it is loaded. Ensure the patient experiences the resistance provided by the knee.
- Have the patient experience the roll-over of the foot. Encourage the patient to apply hip extension.
- Ask the patient to sit down using the stance resistance of the device.
- Show the function of the manual lock as this could be helpful when the knee is turned OFF.
- Ask the patient to walk both short and slightly longer distances under supervision to evaluate the setup of the device fully.

The healthcare professional should instruct the patient on:

- How to operate the device.
- How to understand the warning signals that the device may issue (**see. Device Warning Signals section**).
- How to recognize and activate the different states of the device (**see. Supported Activities section**).

USAGE

Device operation

Turning the Device ON

1. Push and hold the power button for longer than 1 seconds.
2. The device will beep twice.
3. When the status indicator flashes green, the device is ready to use.

Turning the Device OFF

1. Push and hold the power button for longer than 1 seconds.
2. The device will emit 4 short and 1 long beeps and turn off.

Connecting to the App

1. Turn on the device.
2. Open the app and follow the instructions in the app.

Note: The required PIN code can be found on the device label at the back of the device.

Battery

Warning: A damaged power supply can lead to electric shock.

Warning: Always disconnect the power supply from the device before donning.

Caution: The Battery Pack charge connector is Nickel-plated which may cause allergic reactions. Fully charge device before first use.

Charging

1. Plug the power supply into the wall socket.
2. Plug the power supply into the charging port on the back of the device.
3. After charging, disconnect the power supply from the charging port.
4. Disconnect the power supply from the wall socket to save energy.

Note: It takes approximately 3 hours to charge the device battery pack when fully depleted.

Note: The battery pack cannot be overcharged.

The device can be charged regardless of whether it is turned ON or OFF.

Charge the battery daily for optimal performance.

Battery performance is affected by temperature below -10°C (14°F) or over 40°C (104°F), which may result in battery warnings or prevent device to perform as expected.

Note: A fully charged battery pack can last for 48 to 72 hours, depending on activity level.

To read the battery charge level:

Press the power button. The device will emit 1 to 5 beeps, where each beep indicates 20% charge level (i. e. 5 beeps indicate a fully charged device) (**Fig. 5**).

Device mechanical functions

Mechanical Lock

To engage mechanical lock, press mechanical lock button from the left (**Fig. 6**). To unlock, press the lock button from the right.

Supported activities

The device automatically adapts to activities such as walking, standing up and more. It has several different states that the patient must learn to recognize and activate.

Standing State

The device is in standing state by default and it will revert to standing state if it does not recognize a movement pattern.

The device will provide support when the patient applies weight to the prosthesis. The knee will release/swing freely when the prosthesis is unloaded/no weight is applied to it. The patient can take small steps and turns in a natural and physiological way.

Walking on Level Ground

The device provides adaptive support in stance phase to ensure consistent support and controlled swing motion. The device automatically adapts to the patient's walking speed and style.

To walk on level ground:

1. Before the first step, put weight onto the device to trigger the standing state.
 2. Take the first step with either the sound limb or prosthetic limb and continue walking at a preferred pace.
- If using a cane or a crutch, the patient must take care to put weight on the device during stance to maintain support.

Sitting Down

The device automatically detects sitting down motion and provides yielding resistance, allowing the patient to distribute weight equally to both legs and control speed while sitting down.

To sit down:

1. Stand in front of the chair.
2. Put weight on the device.
3. Lean back slightly and start flexing the knee until it yields.
4. Load the device and use the bodyweight to maintain a continuous knee flexion until seated. Once in a sitting position, the device is unlocked and moves freely.

Walking Down Stairs and Ramps

Caution: Always use the handrail when you walk down ramps and stairs.

The knee will give support while it flexes when walking down ramps and stairs.

To walk down stairs or ramps:

1. Take the first step down with the prosthesis.
2. Put weight on the prosthesis. Lean slightly backwards to flex the knee. The knee will flex and give support.
3. Lower the other foot onto the next step or ramp.
4. Continue to walk down the stairs or ramp. The knee will adapt the support to the walking speed.

Note: Leaning backwards keeps from falling forward if balance is lost.

Note: When walking down stairs, ensure to place half of the foot on the step to ensure adequate support.

Device warning signals

In case of warnings, stop using device and troubleshoot as follows:

1. Check battery status and recharge if necessary.
2. Check if the device is hot and allow to cool down.

If this does not help, contact your healthcare professional.

The Össur Logic App can also give information on device warnings.

Table 1 - Warning Signals			
Warning type	Device state	Device Feedback	Action
Low battery warning	Low battery (5% charge)	Audio and/or vibration feedback every 15 seconds (based on user settings). Device status indicator flashing orange. After 10-15 minutes: Audio AND vibration feedback every 15 seconds. Device status indicator flashing orange.	Charge battery
	Critical low battery (0% charge)	Repeated audio and vibration feedback for 20 seconds. Device status indicator flashing red. After 20 seconds: Device stops operating, gives shutdown audio feedback and shuts down.	Charge battery
Critical error	Critical sensor error	Audio and vibration feedback every 2 seconds. Device status indicator flashing red.	Stop using device. Send to Össur for service.
	Sensor calibration error	Audio and vibration feedback every 10 seconds. Device status indicator flashing red.	Calibrate the device
	Battery error	Audio and vibration feedback every 2 seconds on startup. Device status indicator flashing red.	Stop using device. Send to Össur for service.

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

It is recommended to maintain the device in good outward condition by cleaning the surface of the device regularly.

Caution: Do not use solvents stronger than isopropyl alcohol as it might degrade the material.

Caution: Do not use compressed air to clean device.

Note: Device is not delivered in a sterilized state or intended to be sterilized.

Environmental Conditions

The device is Weatherproof.

A Weatherproof device can be used in a wet or humid environment and can tolerate being splashed by fresh water (e. g., rain), no submersion is allowed.

No contact with salt water or chlorinated water is allowed.

Dry with a cloth after contact with fresh water or humidity.
Clean with fresh water in case of accidental exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

Warning: Do not use device in a dusty environment or sand.

Note: If sand or dirt gets into the device, range of motion may be affected.

Table 2 - Environmental Conditions				
	Use	Charging	Shipping	Extended storage
Temperature	-10°C to 40°C (14°F to 104°F)	0°C to 40°C (32°F to 104°F)	-25°C to 70°C (-13°F to 158°F)	0°C to 45°C (32°F to 113°F)
Relative humidity	0% to 100%, condensing	5% to 95%, non-condensing	10% to 100%, condensing	10% to 90%, non-condensing
Atmospheric pressure	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa	700 hPa to 1060 hPa

Note: Prolonged exposure to temperatures below 0°C can result in device malfunction.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.
Recommended interval is every 40 months.
Examine the device regularly for signs of damage. If you find signs of damage, consult your healthcare professional.

Warning: Do not attempt device repairs or remove covers. Risk of injury or damage to device. Send device to Össur for maintenance or repair.

Expected life duration under normal usage is 6 million steps. Depending on patient activity, this corresponds to a duration of use of 4 to 6 years.

SPECIFICATIONS

Table 3 - Device Specifications	
Device weight	1.6 kg (3.5 lbs)
Battery pack specification	Lithium-Ion / 2350 mAh / 33.84Wh
Energy to charge the battery pack	93.82 VAh at 230 VAC / 50 Hz 78.77 VAh at 110 VAC / 60 Hz
Primary packaging weight	1.4 kg (3.0 lbs)
Primary packaging materials	Polypropylene (PP 5) plastic and Polyurethane Foam

Table 4 - Power Supply Specifications	
Manufacturer	FRIWO
Model no.	FW8030M/24
Input	100–240 VAC, 50–60 Hz, 0.6–0.3 A
Output	24 VDC, 1.25A

ELECTROMAGNETIC COMPATIBILITY

Warning: Use of this device adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this device and the other equipment should be observed to verify that they are operating normally.

Warning: Use of accessories and cables other than those specified or provided by the manufacturer of this device could result in increased electromagnetic emissions or decreased electromagnetic immunity of this device and result in improper operation.

Warning: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the device, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The device is suitable for use in any environment except where immersion in water or any other fluid is possible, or in flammable / explosive atmospheres, or where exposure to highly electrical and/or magnetic fields can occur (e. g., electrical transformers, high-power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

The device may be susceptible to electromagnetic interference from portable and mobile RF communications equipment such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements.

Table 5 - Electromagnetic Emissions Compliance Information		
Emissions Test	Standard	Compliance Level
Conducted and Radiated RF Emissions	CISPR 11	Group 1 – Class B
Harmonics Emissions	IEC 61000-3-2	Class A
Voltage Fluctuations / Flicker Emissions	IEC61000-3-3	Complies

Table 6 - Electromagnetic Immunity Compliance Information		
Immunity Test	Standard	Compliance Level
Electrostatic discharge	IEC 61000-4-2	± 8 kV contact ± 15 kV air
Electrical fast transient/burst	IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines
Surge	IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth
Voltage dips	IEC 61000-4-11	0 % UT; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°
		0 % UT; 1 cycle and 70 % UT; 25/30 cycles Single phase: at 0°
Voltage interruptions	IEC 61000-4-11	0 % UT; 250/300 cycle
Power frequency (50/60 Hz) magnetic field	IEC 61000-4-8	30 A/m
Conducted disturbances induced by RF fields	IEC 61000-4-6	3 V 0,15 MHz – 80 MHz 6 V in ISM and amateur bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz
Radiated RF EM fields	IEC 61000-4-3	10 V/m 80 MHz – 2,7 GHz 80 % AM at 1 kHz
	EN ISO22523	12 V/m 26 MHz – 1 GHz

Table 7 - RF Wireless Communications Compliance Information				
Test frequency (MHz)	Band (MHz)	Service	Modulation	Immunity test level (V/m)
385	380-390	TETRA 400	Pulse modulation 18 Hz	27
450	430-470	GMRS 460, FRS 460	Pulse modulation 18 Hz	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	9
745				
780				
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	28
870				
930				
1720	1700-1990	GSM 1800, CDMA 1900, GSM 1900, DECT, UMTS LTE Band 1, 3, 4, 25	Pulse modulation 217 Hz	28
1845				
1970				
2450	2400-2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	9
5500				
5785				

Wireless modules regulated information

USA-Federal Communications Commission (FCC)

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy. If not installed and used in accordance with the instructions, it may cause harmful interference to radio communications. However, there is no ensured specification that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by tuning the equipment off and on, the user is encouraged to try and correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the distance between the equipment and the receiver.
- Connect the equipment to outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

Caution: Exposure to Radio Frequency Radiation – This device must not be co-located or operating in conjunction with any other antenna or transmitter.

Canada – Industry Canada (IC)

This device complies with RSS 210 of Industry Canada.

Operation is subject to the following two conditions:

- this device may not cause interference, and
- this device must accept any interference, including interference that may cause undesired operation of this device.

Caution: Exposure to Radio Frequency Radiation – The installer of this radio equipment must ensure that the antenna is located or pointed such that it does not emit RF field in excess of Health Canada limits for the general population; consult Safety Code 6, obtainable from Health Canada's website www.hc-sc.gc.ca/rpb

Table 8 - Wireless Modules			
Model	Regulatory certificates	Type and frequency characteristics	Effective radiated power
Bluetooth 2.0/4.1 Dual Mode Module: BT121-A-V2	FCC Contains FCC ID: QOQBT121	2402 - 2480 MHz (78 Channels)	0.27 mW/MHz
	Canada Contains transmitter module IC: 5123A-BGTBT121		
	Japan Contains transmitter with certificate number 209-J00171		
Bluetooth 4.1 Module: N550M8CC	FCC Contains FCC ID: O6R2398	F1D 2402-2480 MHz (1 MHz Interval, 79 Channels)	0.42 mW/MHz
	Canada Contains IC: 3797A-2398		
	Japan Contains transmitter with certificate number R203-JN6049		
	Australia & New Zealand ACMA supplier code number is N20233		

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

Note: The primary packaging is made from recyclable materials.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

COMPLIANCE

The device has been designed and tested according to the applicable international standards or in-house defined standards when no international standard exists or applies.

The device has been tested by an independent accredited laboratory and found to comply to:

- EN60601-1
- EN60601-1-2

This device has been tested according to ISO 10328 standard to three million load cycles.

The load level specification is ISO10328 - P6 - 136 kg

ISO 10328 - P6 - 136 kg *)

*) Body mass limit not to be exceeded!

For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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