




Instructions for Use









---

REBOUND<sup>®</sup> ACL  
REBOUND<sup>®</sup> ACL CUSTOM

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	13
ES	Instrucciones para el uso	17
IT	Istruzioni per l'uso	21
NO	Bruksanvisning	25
DA	Brugsanvisning	29
SV	Bruksanvisning	33
EL	Οδηγίες Χρήσης	37
FI	Käyttöohjeet	41
NL	Gebruiksaanwijzing	45
PT	Instruções de Utilização	49
PL	Instrukcja użytkowania	53
CS	Návod k použití	57
RU	Инструкция по использованию	61
JA	取扱説明書	65
ZH	中文说明书	69
KO	사용 설명서	72

Medical Device

## **INTENDED USE**

The device is intended for dynamic support and stabilization of ligaments of the knee.

The device must be fitted and adjusted by a healthcare professional.

### ***Indications for use***

Conditions benefiting from an agonistic force that varies throughout flexion in proportion to anterior cruciate ligament (ACL) tension applied to the thigh and calf in the sagittal plane. These may include:

Anterior cruciate ligament (ACL) injury

- Reconstruction
- Non-surgical treatment

No known contraindications.

### **Warnings and Cautions:**

- Care should be taken not to overtighten the device.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.
- The device is intended for post-surgical or post-injury rehabilitation. It is not intended to be used in moderate or high activity and/or impact sports. A healthcare professional must approve the types of activities that the patient may engage in while using the device.
- This device is a supportive brace only and is not intended or guaranteed to prevent injury.
- Any Cable tension adjustments must be made by a healthcare professional. In the event of an accidental adjustment to the Cable tension, the patient should immediately discontinue the use of the device and contact a healthcare professional.
- Cable tension should only be adjusted with the patient's leg in full extension. Adjusting cable tension in flexion will increase the risk of overloading the knee.
- The device must have an extension and flexion Stop installed in both Hinges during use.
- To avoid discomfort and potential impairment of circulation, the patient should avoid sitting on the frame for an extended period of time.
- To increase comfort while sitting, keep the braced leg extended.

### **GENERAL SAFETY INSTRUCTIONS**

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

## FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

The device can be ordered in a Custom Made variant and is a field adaptable knee brace. It should be modified to fit individual patient contour by a healthcare professional.

### Device Application

1. Unfasten all straps (**Fig. 2**). Ask patient to sit down and fully extend leg while fitting the device.

Ensure proper alignment of the device on the leg (**Fig. 3**).

- **Height positioning:** Align the center of the Hinges (A) slightly above the middle of the kneecap.
- **Side positioning:** Center of the Hinges should be at midline of leg.
- **Condyle Pad (B) positioning:** The Condyle Pads should lightly touch the medial and lateral condyles of the knee. The device is pre-assembled with standard thickness Condyle Pads and is included with pairs of thinner and thicker sized Condyle Pads, along with a hook and loop shim. Use any combination of the Condyle Pads and shim to widen or narrow the width of the device at the knee

**Note:** Ensure that the Cable Covers (C) are on the outside of the device and floating free of obstruction. Check both sides of the device.

2. All straps can be adjusted and/or can be trimmed to length. Fasten the straps in the following order:

- Fasten the Anti-Migration Strap (D) (**Fig. 4: a**). If needed, trim the Anti-Migration Strap (4b)
- While holding the device in place (**Fig. 5: a**), fasten the Upper (E) and then the Lower Shin Strap (F) by inserting them through their respective d-rings (5b, 5c). Position the Shin Pad (G) before fastening the Lower Shin Strap.
- Next, center the Thigh Pad (H) on the leg and then snugly fasten the Upper Thigh Strap (I) (**Fig. 6**).
- Fasten the Lower Thigh Strap Liner (J) (**Fig. 7**). If needed, shorten the strap by detaching the hook tab of the Lower Thigh Strap Liner and trim along its cut lines (**Fig. 8**).
- Center and affix the Hinges at the midline of the leg by:
- Fasten the Posterior Calf Strap (K) (**Fig. 9**).

**Note:** If the device is too anterior, loosen the Upper Shin Strap and tighten the Posterior Calf Strap. If the frame is too posterior, loosen the Posterior Calf Strap and tighten the Upper Shin Strap.

3. Ask the patient to stand up and fully extend leg. Place the required Shear Knob (**Fig. 10: a**) in the diamond shaped hole on the Spool (10b) of the Dynamic Tensioning System™ (DTS) (L).

Two different Shear Knobs are included with the device:

- White: standard load, and
- Black: low load.

Slowly turn the selected Shear Knob clockwise until it breaks and then discard the two pieces (**Fig. 10**). Ensure that the DTS is centered between the Cable Shrouds (M) with equal Cable (N) length on each side of the DTS. If not centered, please refer to the following section Re-tensioning the DTS. Place an included DTS sticker on the Spool to cover the diamond shaped hole and prevent DTS adjustments by anyone other than a healthcare professional (**Fig. 11**).

## ***Device Adjustments***

### **Range of Motion (ROM) adjustment**

The device comes with a 0° extension and 90° flexion Stops installed in the Hinges.

1. Peel Condyle Pads from Hinges to reveal Stop slots and Retaining Screws.
2. Using a screwdriver, remove the Retaining Screws and slide out the unfastened Stops.
3. Insert desired extension and/or flexion Stops.
4. Secure the Stops with Retaining Screws and replace the Condyle Pads on the Hinges.

**Note:** Do not overtighten the retaining screw.

### **Re-tensioning the DTS**

1. Remove the device from the leg and unfasten the DTS from the Lower Thigh Strap Liner.
2. Turn over the DTS to expose the underside.
3. Insert a pen/screwdriver into the hole in the Pawl (**Fig. 12: a**) and apply outward force in the same direction of the arrow to disengage the Pawl from the Spool (**Fig. 12**).
4. While maintaining the disengaged position of the Pawl, pull the Cable out of the DTS until the DTS slides freely on the Cable.
5. Pull the DTS away from the device and center it between the Cable Shrouds. Ensure that the Cable is of equal length on each side of the DTS (**Fig. 13**).
6. Place the desired Shear Knob into the diamond shaped hole on the Spool and slowly turn it clockwise approximately half turn to lock the DTS onto the Cable.
7. Turn over the DTS and re-affix it in its original position on the Lower Thigh Strap Liner.  
**Note:** The Cable should be loose but of equal length on each side of the DTS.
8. Repeat the Device Application instructions to fit the device and re-tension the DTS.

### ***Device Removal***

1. Ask patient to sit down with the leg extended and then unfasten all the straps on the front of the device (**Fig. 2**).
2. Attach the straps back to themselves, making it easier to put the device on next time and prolonging the life of the straps.

### ***Accessories and Replacement Parts***

Please refer to the Össur catalog for a list of available replacement parts or accessories.

## **USAGE**

### ***Cleaning and care***

Washing the device with the soft goods detached allows for more thorough cleaning.

### **Washing Instructions**

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

**Note:** Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

**Note:** Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

### **Hinge**

- Remove foreign materials (e. g., dirt or grass) and clean using fresh water.

### **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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