












Instructions for Use

REBOUND[®] POST-OP ELBOW

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	13
NO	Bruksanvisning	16
DA	Brugsanvisning	18
SV	Bruksanvisning	20
EL	Οδηγίες Χρήσης	22
FI	Käyttöohjeet	24
NL	Gebruiksaanwijzing	26
PT	Instruções de Utilização	28
PL	Instrukcja użytkownika	31
CS	Návod k použití	33
TR	Kullanım Talimatları	35
HR	Upute za uporabu	37
RU	Инструкция по использованию	39
JA	取扱説明書	41
ZH	中文说明书	43
KO	사용 설명서	45
RO	Indicații de utilizare	48
HU	Használati javallat	50
BG	Указания за употреба	52

Medical Device

INTENDED USE

The device is intended to provide controlled range-of-motion for post-traumatic, post-surgical repair of the elbow

The device must be fitted and adjusted by a healthcare professional.

Indications for use

Indications requiring elbow immobilization and/or range of motion restriction. This may include:

- Stable fractures of the elbow, distal humerus, proximal radius, or ulna
- Tendon and ligament injuries or repairs
- Collateral ligament reconstructions
- Elbow hyperextension
- Triceps tendon repair
- Chronic elbow injuries
- Tennis elbow release

No known contraindications.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

1. Unfasten all buckles and straps of the device.
2. Locate the hinge pivot center over the lateral humeral epicondyle (outer elbow) (**Fig. 2**).
3. Device length should be adjusted to provide optimal elbow stabilization. To adjust the position of the Cuffs:
 - Push down the length adjustment button to unlock and slide telescoping paddle until desired length is achieved (**Fig. 3**).

Positioning Guidelines:

- Forearm Cuff (A) extension should extend up to the styloid process (bump on the medial side of wrist), avoid impinging on the styloid process.
- Upper Arm Cuff (B) should extend up to the most comfortable position, avoid impinging the axilla (armpit).
- Adjust the Inner Cuffs (C) if needed.

Note: Ensure all length adjustment buttons are in locked position.

4. The device struts may be bent to achieve a better fit. Place the device on a tabletop with one strut protruding off the edge. Apply downward pressure to bend the strut.
5. Form the cuffs to optimize the fit around the upper and forearm (**Fig. 4**).
6. Adjusting Cuff strap length:
 - Beginning with buckles closest to the elbow, secure buckles and fasten and tighten straps (**Fig. 5**).
 - Straps may be trimmed if necessary (**Fig. 6**).
7. To adjust Range-Of-Motion Hinge (D), pull the flexion (white) or extension (black) button and rotate around hinge until center of button is pointing at desired angle and release (**Fig. 7**).
8. (Optional) Place Range-Of-Motion (ROM) Stop-Locks (E) over flexion and extension buttons to limit ROM adjustment (**Fig. 8**).

Shoulder Strap Application

9. Position the Shoulder Strap (F) over the patient's shoulder and attach the strap clasp with D-Ring to the loop in the forearm paddle (**Fig. 9**).
10. Position the shoulder strap across the back and under the injured arm. Attach the clasp to the D-Ring (**Fig. 10**).
11. Adjust the length of the Shoulder Strap and reposition the strap pad until the forearm is held at the proper elevation and the neck and forearm are held at a comfortable position (**Fig. 11**). Trim the strap to the appropriate length, if needed.

Device Adjustments

Drop-Lock Adjustment

To prevent Range-of-Motion, lock the device by sliding the Blue Drop-Lock (G) toward the center of the hinge after positioning the arm at the appropriate angle as marked by blue lines on the hinge (**Fig. 12**).

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

- Remove pads and straps and wash separately.
- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Note: Do not immerse hinge and uprights in water.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

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