

Instructions for Use

REBOUND[®] DIABETIC WALKER


		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	18
NO	Bruksanvisning	21
DA	Brugsanvisning	24
SV	Bruksanvisning	27
FI	Käyttöohjeet	30
NL	Gebruiksaanwijzing	33
PT	Instruções de Utilização	37
PL	Instrukcja użytkownika	40
TR	Kullanım Talimatları	43
RU	Инструкция по использованию	46
JA	取扱説明書	49
ZH	中文说明书	52
KO	사용 설명서	55
RO	Indicații de utilizare	58
BG	Указания за употреба	61

9

10

ENGLISH

Medical Device

Caution: This device contains natural rubber **latex** which may cause allergic reactions.

INTENDED USE

The device is intended for plantar surface offloading and immobilization. The device must be fitted and adjusted by a healthcare professional.

Indications for use

- Plantar foot ulcers and wound prevention.

Contraindications

This device is not intended for use on patients with chronic Charcot feet with a rocker bottom deformity.

Warnings and Cautions:

The device must be prescribed by a healthcare professional.

- Weight limit: The patient's maximum weight is 136 kg (300 lbs).
- Hexplugs should only be removed or altered by a healthcare professional.
- Wounds under treatment should be monitored closely for maintenance of pressure relief and healing progress.
- The insole should not be used in conjunction with any other footwear device.
- In case of emergency, the compliance strap can be cut with scissors for the removal of the device.
- Before reapplying the device inspect:
 - The foot for any rash, sores, new ulcers, or other indications of skin irritation.
 - The boot, bladders, and liners for any damage.
- If the device has become loose or uncomfortable to wear, readjust the straps as necessary and contact your healthcare professional if discomfort is not relieved.
- This device is only intended for treating or preventing plantar ulcers of the foot.
- Patient should use caution when walking on slippery or wet surfaces to avoid injury.

Do not remove the device unless instructed by a healthcare professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Device Application

1. Open the liner by first unfastening straps and removing the front shell.
2. Cover the patient's wound with an appropriate dressing.

Customization of the Insole:

3. Apply the supplied sock over the patient's foot/dressing, prior to placing the foot on the heated insole.
4. Do not remove hexplugs prior to heat molding.
5. Pre-heat the oven to 110 °C (230 °F) and then heat the insole for three minutes.
6. Have the patient stand, wearing the sock, on the heated insole for one minute (**Fig. 1**).
7. Have the patient remove the foot from the insole.
8. Mark the wound location on the wound dressing with a transferable mark using iodine, lipstick, indelible ink pencil, etc. (**Fig. 2**).
9. Reposition the insole within the liner, with the heel properly located at the back of the device.
10. While seated, place the foot inside the device. Ensure the proper location of the patient's foot before transferring the mark to the insole (**Fig. 3**).
11. Wrap liner securely around foot and leg. Have user bear weight to ensure transfer of relief zone mark.
12. With the patient seated, remove the device. Remove liner and insole. Lightly press on the transferred mark on the insole and remove elevated hexplugs (**Fig. 4**).

Caution: Do not make the relieving area too large, which may allow the ulcer to contact the bottom of the device while weight-bearing.

13. Re-insert the insole into the liner with the hexplugs facing away from the patient's foot. Re-insert the liner and insole into the device.
14. While seated, place the leg into the device.
15. Wrap liner securely around the foot and leg and place the front shell over the foot and leg.
16. In the following order, fasten the ankle strap, the foot strap, and then the calf strap (**Fig. 5**). Tighten straps until snug and comfortable. Strap adjustments may be necessary with changes in lower limb volume.
17. The lower leg should be kept at a 90° angle with the foot. Use the positioning portals at the back of the device to confirm that the ankle is centered and that the heel is not pressed against the back of the device (**Fig. 6**).
18. Inflate air chambers by pressing "+" on the air bulb, and deflate by pressing "-". Inflate or deflate until snug (**Fig. 7,8**).

Note: The Air System is designed with a feature to prevent over-inflation.

19. A disposable compliance strap may be placed around the ankle of the device to inhibit removal. Place the strap through the strap guide on the front shell before securing the compliance strap (**Fig. 9**).

Note: Use of the compliance strap is recommended. The compliance strap is for single use only and should be replaced by the clinician when the device is removed. Cut the compliance strap with scissors to remove the device.

Device Adjustments

The device circumference around the calf can be widened by cutting Flex Edge® (**Fig. 10**) to accommodate patients with a large calf.

Note: Ensure to remove the calf strap before cutting the Flex Edge®.

Device Removal

Before removing the device, deflate air chambers by pressing “-“. Undo the straps in any order and attach the straps back to themselves. This will make it easier when fitting the device next time and will prolong the life of the straps. Open the Liner to lift the foot and lower leg out of the device.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing Instructions

Device plastic can be washed with a damp cloth and mild soap. Remove liner from the device before washing.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Insole:

- Remove the entire insole from the device. Do not remove any additional hexplugs.
- Once the insole is fully dried it can be re-inserted into the device.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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