







Instructions for Use

REBOUND® CARTILAGE
REBOUND® CARTILAGE CUSTOM

		3
EN	Instructions for Use	4
DE	Gebrauchsanweisung	6
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	12
IT	Istruzioni per l'uso	14
NO	Bruksanvisning	17
DA	Brugsanvisning	19
SV	Bruksanvisning	22
EL	Οδηγίες Χρήσης	24
FI	Käyttöohjeet	27
NL	Gebruiksaanwijzing	30
PT	Instruções de Utilização	32
PL	Instrukcja użytkowania	35

ENGLISH

Medical Device

Caution: This device contains natural rubber **latex** which may cause allergic reactions.

INTENDED USE

The device is intended for unicompartmental unloading of the knee. The device must be fitted and adjusted by a healthcare professional.

Indications for use

Cartilage knee trauma that may benefit from unloading or movement restriction.

For knee conditions that may benefit from unicompartmental load reduction or range of movement restriction, such as:

- Articular cartilage defect repair
- Meniscal repair
- Avascular necrosis
- Condylar bone marrow lesions (bone bruises)
- Tibial plateau fractures
- Mild to severe unicompartmental knee osteoarthritis

No known contraindications.

Warnings and Cautions:

- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- Make sure the device fits properly to minimize the possibility of skin irritation. Increase usage time gradually as the skin adapts to the device. If redness appears, temporarily decrease the usage time until it has subsided.
- Care should be taken not to overtighten the device.
- Make sure the device fits properly to achieve effective pain relief.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

1. Open Upper (A) and Lower (B) Buckles. Ask patient to sit down and extend leg while fitting the device. Ensure that the Upper (C) and Lower (D) SmartDosing® Dials are set to the “0” position. Place the device on the patient’s leg with the Hinge (E) on the affected side of the knee.
Ensure proper alignment of the device on the leg (**Fig. 2**).
 - Height positioning (**2a**): Align Össur logo on the Hinge with upper edge of the patella.
 - Side positioning (**2b**): Front edge of the Hinge cap should be at midline of leg.
2. Fasten the Buckle buttons to their color-matched Keyholes (F, G).
 - Place the blue Lower Buckle button into the blue Calf Shell Keyhole (F) above Buckle Stability Shelf (H) and use the palm of the hand to snap the Lower Buckle closed (**Fig. 3**). Adjust Calf Strap (I) to appropriate length by tensioning around the calf and folding into the Alligator Clip (J) so that it keeps the device securely and correctly positioned on leg.
 - Bend the patient’s knee to 80°. Place the yellow Upper Buckle button into the yellow Thigh Shell Keyhole (G) and use the palm of the hand to snap the Upper Buckle closed (**Fig. 4**). Adjust Thigh Strap (K) to appropriate length by tensioning around leg and folding into Alligator Clip.
3. Adjust the Lower Cartilage Protection Strap™ (CPS) (L) length until it sits firmly against the leg, and then fold it into Alligator Clip (**Fig. 5**). At this point, patient should not experience any tension or unloading. Adjust the Upper CPS (M) in the same way.
4. Ask patient to bend knee to 90° (**Fig. 6**). Turn upper and lower SmartDosing Dials clockwise (**Fig. 7a**) until indicators are in “5” position (**Fig. 7b**). Have the patient stand up and take a few steps to verify correct positioning of the device and tightness of the straps.
 - Determine optimal CPS tension based on patient’s pain relief feedback.
 - If patient requires more or less tension with indicator at the “5” position, adjust the length of the CPS accordingly.
 - Aim for final SmartDosing Dial setting at the “5” position since this will give patient the ability to adjust dosing during daily life activities.
5. When final fit has been confirmed, trim the straps to the appropriate length starting with the Calf Strap so that device sits correctly on leg while trimming other straps.
 - Trim the straps sufficiently so that the alligator clips are positioned away from the popliteal area (**Fig. 8**). This reduces the bulk behind the knee.

Device Removal

1. Ask patient to sit down with the leg extended.
2. Turn both Dosing Dials counterclockwise to release tension on the CP Straps and adjust back to the “0” position (**Fig. 6**).
3. Bend the patient’s knee to 90° and open both Lower and Upper Buckles.
4. Pull the Buckles buttons out of Keyholes.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Hinge

- Remove foreign materials (e. g., dirt or grass) and clean using fresh water.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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