The patient should stop using the device and contact a healthcare professional:
 • If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
 • If any pain, skin irritation, or unusual reaction occurs with the use of the device. The device is for single patient – multiple use.

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Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener. Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:
 • Device not maintained as instructed by the instructions for use.
 • Device assembled with components from other manufacturers.
 • Device used outside of recommended use condition, application, or environment.

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



MD Medical Device

INTENDED USE

The device is intended for external support or stabilization of the knee.

Indications for use

For knee conditions that may benefit from increased stability around the knee, such as:

- Mild medial/lateral knee instabilities.

No known contraindications.

Warnings and Cautions:

If you have a medical condition that decreases blood circulation in your limbs, then consult your healthcare professional before you use the device.

Such medical conditions include, but are not limited to, diabetes or peripheral vascular disease.

Inspect the skin under the device's contact area before and after each use.

The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

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- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device. The device is for single patient – multiple use.

USAGE

Device Application

Note: For best results when contouring hinge arms, remove hinges and bend hinge arms as desired (Fig. 2). Re-apply hinges before application, ensuring that the hinge arms are completely seated and secured within the hinge pockets.

Note: Universal buttress can be adjusted to multiple positions to optimize patella support (sleeve version).

1. Unfasten the straps.
 2. If necessary, add the buttress accessory and adjust it to the patella (sleeve version only).

3. Depending upon the version of the device chosen, with the leg in full extension, grasp the device with two hands on the side and pull the device over the knee (sleeve version) or wrap the device around the knee and fasten the contact closures over the front of the knee (wrap version). Position the device, such that the kneecap opening and the middle of the hinge are centered on the kneecap (Fig. 3).

4. Insert the straps through the corresponding D-ring (Fig. 4) and fasten the strap at the calf and then at the thigh (Fig. 5).

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