




Instructions for Use



REBOUND[®] KNEE

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	9
ES	Instrucciones para el uso	11
IT	Istruzioni per l'uso	13
NO	Bruksanvisning	15
SV	Bruksanvisning	17
EL	Οδηγίες Χρήσης	19
FI	Käyttöohjeet	21
NL	Gebruiksaanwijzing	23
PT	Instruções de Utilização	25
PL	Instrukcja użytkowania	27
JA	取扱説明書	29
ZH	中文说明书	31
KO	사용 설명서	33

5

ENGLISH

Medical Device

INTENDED USE

The device is intended for external support or stabilization of the knee.

Indications for use

For knee conditions that may benefit from increased stability around the knee, such as:

- Mild to moderate ACL or PCL, MCL and LCL instabilities, sprains, and strains.

No known contraindications.

Warnings and Cautions:

If you have a medical condition that decreases blood circulation in your limbs, then consult your healthcare professional before you use the device.

Such medical conditions include, but are not limited to, diabetes or peripheral vascular disease.

Inspect the skin under the device's contact area before and after each use.

- If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this device, discontinue the use immediately and consult with a healthcare professional.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

Read these instructions carefully before use. Keep them for future reference.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

USAGE

Device Application

Note: For best results when contouring the hinge arms, remove hinges and bend hinge arms as desired. Re-apply hinges before application, ensuring that the hinge arms are completely seated within the hinge pockets and the loop at the back of the hinge is secured to the hook on the device.

Note: Do not remove the yellow pull tab until after the initial fitting.

1. Unfasten the straps.
2. Depending upon the version of the device chosen, with the leg in full extension, either pull the device over the knee using the pull-tabs (sleeve version) or wrap the device around the knee (lower portion first) and fasten the contact closures over the front of the knee (wrap version).

Position the device, such that the kneecap opening and the middle of the hinge are centered on the kneecap (**Fig. 1**).

3. Fasten the straps at the calf and thigh in the order shown by inserting them through the corresponding D-ring (**Fig. 2**). Straps should wrap inside the pull-tabs at the top of the device (**Fig. 3**).

Note: When the final fit is confirmed, remove the yellow pull tab to secure the strap in position (**Fig. 3**).

Device Adjustments

Range of Motion (ROM) Adjustment

The device comes with 10° extension stops installed and includes the following additional stops:

- Extension stops are labeled 0°, 10°, 20°, 30° and 40°.
 - Flexion stops are labeled 0°, 45°, 60°, 75°, and 90°.
1. To remove the extension stop, use a Phillips head screwdriver to remove the fixing screws (**Fig. 4**).
 2. Use a screwdriver or pen to slide the 10° stop out of both hinges. Re-apply the desired extension stop and gently secure them with the retaining screws (do not overtighten) (**Fig. 5**).
 3. Follow the same procedure for the addition of a flexion stop on the other side of each hinge.

Device Removal

- Undo the straps.
- Grip the device at the level of the calf and pull it off in a downward direction or unwrap the device.

Cleaning and care

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Note: Remove hinges and secure the hook and loop fasteners before washing to avoid damaging the textiles.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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