











Instructions for Use

TOTAL KNEE[®] 2000
TOTAL KNEE[®] 2100


| | | |
|----|-----------------------------|-----|
| | | 3 |
| EN | Instructions for Use | 6 |
| DE | Gebrauchsanweisung | 13 |
| FR | Notice d'utilisation | 22 |
| ES | Instrucciones para el uso | 30 |
| IT | Istruzioni per l'uso | 38 |
| NO | Bruksanvisning | 46 |
| DA | Brugsanvisning | 54 |
| SV | Bruksanvisning | 62 |
| EL | Οδηγίες Χρήσης | 70 |
| FI | Käyttöohjeet | 79 |
| NL | Gebruiksaanwijzing | 87 |
| PT | Instruções de Utilização | 95 |
| PL | Instrukcja użytkowania | 103 |
| TR | Kullanım Talimatları | 111 |
| RU | Инструкция по использованию | 119 |
| JA | 取扱説明書 | 127 |
| ZH | 中文说明书 | 134 |
| KO | 사용 설명서 | 141 |


12


A

B

C

ENGLISH

Medical Device

DESCRIPTION

The device is a polycentric knee joint with a Geometric Locking System and a 3-phase hydraulic swing control that accommodates changes in walking speed.

Proximal adapter options are available, and the device has a Distal Tube Connector.

Interchangeable Bumpers (**Fig. 7A,7B**) allow for knee flexion to be adjusted to suit the patient's preference and activity levels. The adjustable Extension Promoter (7E) assists in limiting heel rise and moving the device into a fully extended position.

INTENDED USE

The device is intended as a part of a prosthetic system that replaces knee function of a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for moderate to high impact use, e.g., walking and occasional running.

Total Knee 2000

The weight limit for the device is 100 kg.

Total Knee 2100

The weight limit for the device is 125 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Warning: Avoid placing hands or fingers near moving joints.

Caution: Do not adjust other screws than described in these instructions.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Warning: Risk of structural failure. Components from other manufacturers have not been tested and may cause excessive load on the device.

Proximal Adapter Options

The following proximal adapter options are available for the device (**Fig. 1**):

- 4-Prong Socket Adapter (1A),
- Female Pyramid Adapter (1B),

- Euro 4-Hole Adapter (1C), or
- Male Pyramid Adapter (1D).

Caution: Only use the specified Total Knee proximal adapters, please refer to Össur catalogue for more information.

1. Fully engage the threads of the selected proximal adapter.
Caution: Adapters from other manufacturers may not fully engage all the threads and compromise the strength of the prosthesis.
2. Align the Pinch Bolt (1E) of the Male Pyramid Adapter (1D), 4-Prong Socket Adapter (1A), and the Euro 4-hole Adapter (1C) perpendicular to the anterior/posterior axis of the device, except for the Female Pyramid Adapter (1B), which should be angled 45°.
Caution: Maximum connection strength will not be established unless the adapter is correctly aligned.
3. For secure fastening to the device, tighten the adapter's Pinch Bolt according to its instructions for use.
Caution: Do not tighten the proximal adapter against the housing of the device.

Caution: When prong adapters are used, resin may get under the adapter and hinder full thread engagement. Carefully remove resin without damaging the adapter so that adapter threads can be fully engaged when mounted on the device.

Caution: Contact between socket/adapter and proximal 2/3 portion of the device's back link must be avoided when fully flexing the knee (**Fig. 2**). Contact in this area can damage the device. Acceptable contact is between the socket and the distal 1/3 portion of the device's back link.

ALIGNMENT INSTRUCTIONS

Bench Alignment (Fig. 3)

Alignment Goal

Alignment reference line (B) should:

- pass through midpoint of socket at the ischial tuberosity level (D)
- pass through pivot axis (A)
- fall at the 1/3 mark on the inside of the foot cover.

Note: Prioritize knee alignment over foot alignment if there is a mismatch.

Alignment Instructions

1. Position the foot so that the alignment reference line (B) falls at the 1/3 mark on the inside of the foot cover (with the foot cover and shoe on). Consider the external rotation of the foot.
2. Use the applicable adapters to connect the knee to the foot and establish the correct knee center height.
3. Position the knee so that the alignment reference line passes through pivot axis (A)
4. On the lateral side of the socket, make a first mark at the midpoint of the socket at the ischial tuberosity level (D). Make a second mark at the midpoint of the socket distally (E). Draw a line through both marks.
5. Position the socket so the alignment reference line (B) goes through first mark at the midpoint of the socket at the ischial tuberosity level (D).
6. Adjust the socket flexion to 5° in addition to the existing position (i.e., hip flexion contracture) and set the height of the full prosthesis.
7. Use the applicable adapters to connect the knee to the socket.

Warning: After adjustments, all screws must be secured with a medium strength threadlocker and tightened with the correct torque.

Caution: Adapters used on the distal connection must be cut straight and inserted down to the end stop of the device's tube receiver. No spacer should be used (**Fig. 4**).

- Proximal Adapter Pinch Bolt (**Fig. 1E**): Please refer to the respective proximal adapter's instructions for use.
- Tube Clamp Screw (**Fig. 5A**): 15 Nm.

Bumper Selection (Fig. 6,7)

The selection of a bumper determines the amount of stance flexion.

1. Carefully insert a flat head screwdriver or a similar flat tool under the bottom of the Bumper Cover (7D) and remove it (**Fig. 6**).
2. Remove the Bumper (7A,7B) by inserting a small and flat hand tool underneath it and lift upwards. Leave the underlying Shim (7C) in place.
3. Select a Bumper to install based on which Total Knee device is being used, and on the patient's weight:

| Total Knee 2000 (Fig. 7: A) - The device is shipped with the Red Bumper installed. | | | |
|---|-------------------|--------------|-------------|
| Weight kg | Weight lbs | Color | Type |
| 35 - 50 | 80 - 110 | Green | Soft |
| 50 - 70 | 110 - 155 | Black | Regular |
| 70 - 90 | 155 - 200 | Red | Firm |
| 90 - 100 | 200 - 220 | Yellow | Extra Firm |

| Total Knee 2100 (Fig. 7: B) - The device is shipped with the Yellow Bumper installed. | | | |
|--|-------------------|--------------|-------------|
| Weight kg | Weight lbs | Color | Type |
| 50 - 70 | 110 - 155 | Black | Soft |
| 70 - 90 | 155 - 200 | Red | Regular |
| 90 - 110 | 200 - 240 | Yellow | Firm |
| 110 - 125 | 240 - 275 | Orange | Extra Firm |

4. Re-install the Bumper Cover (7D) by hooking the groove under the chassis front and pushing the two pins into the chassis holes.

Influencing Factors on Bumper Selection:

- Active patients may require a stiffer bumper.
- Patients who are accustomed to wearing a prosthesis without a stance-flexion feature may initially prefer a stiffer Bumper. After a trial period of two weeks, selecting a softer Bumper is recommended. Increasing the softness of the Bumper allows the patient to get familiarized with the device and facilitates the transition to using its stance-flexion feature.

Static Alignment

- Make sure the patient stands with equal weight on both legs.
- Check for correct prosthesis length.
- Check internal / external rotation.
- Check for correct load on toe and heel.

Dynamic Alignment

Ensure the patient is familiar with the functioning of the device.

Ensure the patient can activate and deactivate the geometric lock.

Swing Control (Fig. 8)

The hydraulic system is supplied in its basic setting.

In the basic setting, Valves "F" and "H" are both opened by $\frac{3}{4}$ turn. Valve "E" is opened by $1\frac{1}{2}$ turns.

Caution: For safety, please make initial assessments and adjustments with the patient standing between parallel bars.

Adjusting the Flexion and Extension valves

The following valves can be adjusted:

- Valve "F": Affects swing flexion resistance above 60°.
- Valve "H": Affects swing flexion resistance from 0° - 60°.
- Valve "E": Affects swing extension resistance.

Start with adjusting the flexion with the patient walking slowly and at faster walking speeds.

For each adjustment, turn the valve in small increments of about a $\frac{1}{4}$ of a turn. Always check the result right after each adjustment step.

Possible observations and actions to take are outlined in the following:

- If the patient is walking fast and you observe excessive heel rise, increase the swing flexion resistance by turning Valve "F" to the right until heel rise is normalized.
- To smooth flexion resistance changes between 0° and 60°, turn Valve "H" to the right or left to increase or decrease the resistance, respectively.

Afterwards, adjust the extension to harmonize the gait.

- Increase the swing extension resistance by turning Valve "E" to the right until terminal impact at full extension is reduced. Patient should feel a slight bump at full extension.
- Fine tune Valves "F" and "E" until a smooth and secure gait pattern is achieved for slow and fast walking speeds.

Caution: The valves should never be fully closed. Over tightening valves or flexing knee with all valves completely closed may damage valves.

Warning: Valves adjustment must always allow for complete extension of the knee. Excessive swing extension resistance will prevent knee from moving into full extension and may compromise stance security.

Note: Generally, it is not required to adjust Valve "H" unless the patient's walking speed is slow and excessive heel rise is observed. Valve "H" is covered by an adhesive label on the right side of the device's hydraulic housing (8H).

Bumper adjustment

As needed, try different Bumpers, and select the one that best suits the need of the patient.

The final bumper selection might differ from the initial selection based on the patient's weight (please refer to the previous section: Bumper Selection).

The torque required to release the Geometric Lock and the stance flexion motion:

- increases when using a soft Bumper, and
- decreases when using a firm Bumper.

Depending on which foot is selected, combine:

- softer bumper with soft heels, and
- firmer bumper with firm heels.

Bumper height adjustment

A thin plastic plate, or a Shim (**Fig. 7C**), is placed underneath the Bumper (7A,7B) to adjust its height. There are two types of Shims:

- Grey: 0.38 mm thick, and
- Clear: 0.76 mm thick (the device is shipped with the Clear Shim installed).

Adding Shims changes the release point of the Geometric Lock that allows for knee flexion during the swing phase. Adjusting the height of the Bumper affects knee stability, along with ease of walking and sitting:

- Higher Bumper: facilitates the release into swing flexion and decreases stance flexion.
- Lower Bumper: increases the device's stability and stance flexion, while delaying the release into swing flexion.

Gradually increase the height of the bumper until the patient feels premature release of the knee. Reduce the number of Shims until the patient finds a suitable balance between device stability and at which point the knee releases into swing flexion.

Warning: Do not install more Shims than the two different types combined. Adding more Shims can deactivate the Geometric Lock, and negatively affect the stability of the knee during stance.

Note: If no Shims are used, the device will exhibit some flexion/extension play during stance.

Extension Promoter

The Extension Promoter reduces excessive heel rise and promotes faster knee extension. It is delivered in the minimum setting.

- Flex the knee to access the Extension Promoter Adjustment Screw (7E).
- To reduce excessive heel rise, especially when walking fast, turn the Extension Promoter Adjustment Screw clockwise.

Extension Assist

The Extension Assist is an optional feature that facilitates full knee extension. Use of the Extension Assist is only recommended when the patient is unable to extend the knee to engage the Geometric Lock.

Three springs of different stiffness are included with the Extension assist:

- Yellow: regular (comes pre-installed with the unit),
- Blue: firm, and
- Red: extra firm.

Note: Before installing the Extension Assist, try reducing the swing extension resistance by turning Valve "E" (**Fig. 8: E**) to the left to facilitate full knee extension.

Installation of the Extension Assist (Fig. 9)

1. Remove the Bumper Cover (A), Bumper (B), and Shims (C).
2. Install the Spring Cup (D) into the Hole (E) inside the Distal Tube Connector.
3. Firmly seat the Spring Cup by gently tapping it with a rubber mallet (or a similar hand tool).
4. Place the Ferrule (F) over the Cable (G).
5. Place the Spring (H) over the Cable and Ferrule (try the yellow spring first).
6. Insert the Cable and Spring into the Spring Cup, with the ball end of the Cable first.

- Slide ball end of the Cable through the Hole (E) in the Tube Connector and insert the ball into the Slot on the device's Bottom Link (I) by using needle-nose pliers or a small screwdriver.
- Re-install the selected Shims, Bumper, and Bumper Cover.
- Let the patient walk and adjust the swing extension resistance with Valve "E" (**Fig. 8: E**). If sufficient extension is not obtained with a fully open Valve "E", try a stiffer spring.
- Readjust the swing flexion resistance after installing the Extension Assist unit with Valves "F" and "H" (**Fig. 8: F, H**).

Follow-up check

- Reassess adjustments after a couple of weeks.
- It may be necessary to change the Shim and Bumper adjustments after the patient has gained some experience with the device.
- Over time, the patient's increased device experience, confidence during use, and proficiency may warrant a change to a firmer Bumper and/or an increase in Shim thickness.

USAGE

Cleaning and care

Wipe the device with a soft cloth. Do not use solvents.

Environmental Conditions

The device should not come in contact with fresh water, salt water or chlorinated water.

Caution: Device should not be used in dusty environment. Exposure to sand, talcum, or similar should be avoided.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

The stance flexion Bumper must be replaced if cracked or worn, depending on the patient's activity level.

Check for damage, excessive wear and dirt.

If the device is subjected to excessive moisture or a corrosive environment, it is recommended to clean and lubricate the device frequently.

In a slightly flexed position, three sets of roller bearings will be visible (**Fig. 10A**). Place few drops of general purpose oil or sewing machine oil on the roller bearings. Move the knee several times and wipe off excess oil with a soft cloth. No other parts need external lubrication.

Wipe the device with a soft cloth moistened with small amount of general purpose oil or sewing machine oil.

Caution: Do not use compressed air to clean device. Air forces pollutants into bearings and may cause malfunctions and wear.

Caution: Do not use talcum powder for lubrication.

SUPPLEMENTARY DEVICES

Cosmetic Options (Fig. 11,12)

The device is intended to protect the prosthesis and to restore the appearance of a missing limb.

- A discontinuous cosmetic solution that does not interfere with knee function is recommended.
- Continuous cosmetic solutions are also a possibility.

Following cosmetic solutions are available, please refer to Össur catalogue for further information:

Total Knee 2000

- Knee Cap, Small, Black (11A),
- Knee Cap, Large, Beige (11B),
- Shin Fairing, Black (11C),
- Cosmetic Foam Cover, Regular (12A), and
- Cosmetic Foam Cover, Large (12B).

Total Knee 2100

- Knee Cap, Small, Black (11D),
- Knee Cap, Large, Beige (11E),
- Shin Fairing, Black (11F), and
- Cosmetic Foam Cover, High Activity, Large (12C).

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance

This device has been tested according to ISO 10328 standard to three million load cycles.

Depending on patient activity, this may correspond to 3-5 years of use.

Total Knee 2000

The load level specification is ISO10328 – P5 – 100 kg.

Total Knee 2100

The load level specification is ISO10328 – P6 – 125 kg.

| | |
|---|--|
| ISO 10328 - "P" - "m"kg *)  | |
|  | <p>*) Body mass limit not to be exceeded!</p> <p>For specific conditions and limitations of use see manufacturer's written instructions on intended use!</p> |

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