






Instructions for Use

UNITY[®] VALVE TT
UNITY[®] VALVE TF
UNITY[®] TF KIT


| | | |
|----|-----------------------------|----|
| | | 3 |
| EN | Instructions for Use | 6 |
| DE | Gebrauchsanweisung | 10 |
| FR | Notice d'utilisation | 14 |
| ES | Instrucciones para el uso | 18 |
| IT | Istruzioni per l'uso | 22 |
| NO | Bruksanvisning | 26 |
| DA | Brugsanvisning | 30 |
| SV | Bruksanvisning | 34 |
| EL | Οδηγίες Χρήσης | 38 |
| FI | Käyttöohjeet | 42 |
| NL | Gebruiksaanwijzing | 46 |
| PT | Instruções de Utilização | 50 |
| PL | Instrukcja użytkowania | 54 |
| TR | Kullanım Talimatları | 58 |
| RU | Инструкция по использованию | 62 |
| JA | 取扱説明書 | 66 |
| ZH | 中文说明书 | 70 |
| KO | 사용 설명서 | 73 |

8

ENGLISH

Medical Device

DESCRIPTION

The Unity Valve Kit is an auto expulsion valve for use in a vacuum socket system.

Components: (Fig. 1 & 2)

- A. Cap Key
- B. Cap
- C. Rubber Membrane
- D. Air Filter
- E. Valve Inner Housing
- F. O-ring
- G. Valve Outer Housing
- H. Check Valve
- I. Tube Connector
- J. Tube
- K. Foam Cap
- L. Set Screw
- M. Socket Insert

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

Directly to the socket

1. Fabricate the socket as usual.
2. Mark the position of the valve onto the socket.
3. Using an 8.5 mm drill bit, make a hole perpendicular to the socket wall.
4. Tap the hole using M10x1.5 tap.
5. Prior to shortening the threads, make sure that the Inner Housing has been separated from the Outer Housing.

6. Ensure the Cap, Rubber Membrane and Air Filter have been removed prior to shortening the threads to prevent debris lodging in the airways.
 7. Measure socket wall thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.
- Warning:** Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off (**Fig. 3**).

Using the Socket Insert (Fig. 2)

1. The Socket Insert is used as a fixation device for the valve, and should always be used when fabricating a flexible inner socket.
2. Prepare the plaster cast as usual.
3. Fix the Socket Insert to the cast using the Set Screw. Ensure it sits flush to the cast.
4. Plaster may have to be removed or added.
5. Apply the Foam Cap to the Set Screw.
6. Fabricate the socket as usual.
7. Grind material away over the socket insert until Foam Cap is fully exposed.
8. Unscrew the Set Screw to remove socket from cast.
9. Ensure the Cap, Rubber Membrane and Air Filter have been removed prior to shortening the threads to prevent debris lodging in the airways.
10. Measure socket wall and Socket Insert thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.

Warning: Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off (**Fig. 3**).

Valve Installation

1. Assemble the Inner Housing to the Outer Housing. Ensure that the Inner Housing snaps in place.
2. Apply silicone adhesive to the threads and housing, where it comes in contact with the socket.
Note: It is not recommended to use epoxy or similar as it can compromise the seal.
3. Install valve and adjust the orientation of the Tube Connector. Tighten with a 5 mm hex key.
Note: Be careful not to over-tighten as it may damage the threads.
Note: After installation distribute silicone adhesive around the valve / socket contact point to ensure an air tight seal.
4. Ensure the valve parts and airways are free from debris before final assembly.
5. Replace Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity. (**Fig. 4**)
Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.
6. Insert one Air Filter to Valve Inner Housing and Rubber Membrane to Cap (**Fig. 1**).
7. Screw the Cap on using the Cap Key.

For Unity TF Kit:

Note: The Tube should be placed on the medial side of the prosthesis.

1. Thread the Tube through the tube fixations.
2. Place the large tube fixation on the prosthetic knee (**Fig. 5a**).
3. Place the small tube fixation on the pylon above knee (**Fig. 5b**).
4. Fix the tube clamp to the pylon between foot and knee (**Fig. 5c**).
5. Make sure the tube w/ protective sleeve is placed next to the knee joint and not behind it (**Fig. 5d**).

For Unity TF Kit used with NAVii® knee and cover:

Note: The Tube should be placed on the medial side of the prosthesis.

1. Stamp two 8mm holes on the medial side of the knee cover (**Fig. 6**). Use the small circular indentations on the knee cover as placement guides.
2. Attach the proximal and distal tube clips to the knee cover by inserting the tube clips through the holes and securing them with the rivets (**Fig. 7**).
3. Thread the tube through, then place the small tube fixation on the pylon above the knee (**Fig. 8a**).
4. Place the tube clamp on the pylon between foot and knee (**Fig. 8b**).
5. Thread the tube through the tube clamp (**Fig. 8b**) and the two tube clips (**Fig. 8c**).
6. Make sure the tube w/protective sleeve is placed next to the knee joint and not behind it (**Fig. 8d**).

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Disassembly

1. Unscrew the Cap using the Cap Key.
2. Remove the Rubber Membrane carefully.
3. Remove filter material.
4. Remove Check Valve.
5. Clean all parts with water and mild detergent.
6. Hold down the release button while cleaning airways.

Reassembly

1. Re place Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity **(Fig. 4)**.
2. Replace the filter material to Valve Inner Housing **(Fig. 1)**.
3. Replace the Check Valve, note its orientation **(Fig. 1)**.
4. Screw the Cap on using the Cap Key.

Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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