






Instructions for Use





---

UNITY<sup>®</sup> FOR PRO-FLEX<sup>®</sup> XC &  
PRO-FLEX<sup>®</sup> XC TORSION

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	9
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	18
IT	Istruzioni per l'uso	22
NO	Bruksanvisning	27
DA	Brugsanvisning	31
SV	Bruksanvisning	35
EL	Οδηγίες Χρήσης	39
FI	Käyttöohjeet	44
NL	Gebruiksaanwijzing	48
PT	Instruções de Utilização	52
PL	Instrukcja użytkowania	56
TR	Kullanım Talimatları	60
RU	Инструкция по использованию	64
JA	取扱説明書	69
ZH	中文说明书	73
KO	사용 설명서	76
ET	Kasutusjuhend	80

### DESCRIPTION

The device is a noninvasive mechanical vacuum pump intended for vacuum assisted prosthetic suspension.

Suspension is obtained through the act of a pump which draws air out from the socket in each step creating an elevated vacuum inside the socket interface.

#### Components (Fig. 1)

1. Membrane Housing
2. Inlet Port
3. Outlet Port
4. Duckbill
5. Bolt
6. Protection Pad
7. Tube
8. Tube Clamp

### INTENDED USE

The device is intended as a part of a prosthetic system that provides suspension via subatmospheric pressure.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

It is intended to be used with the Pro-Flex XC or Pro-Flex XC Torsion, respective foot cover and a Unity Valve Kit for connection to a suspension socket with an Iceross® Seal-In Liner.

#### *Indications for use*

- Lower limb amputation and/or congenital deficiency
- No known contraindications

### GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

**Warning:** If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

### ASSEMBLY INSTRUCTIONS

1. Clean the surface of the carbon blade and the Protection Pad with acetone or isopropyl alcohol (Fig. 2).
2. Apply Loctite 495 to the bottom surface of the Protection Pad, labelled with size (Fig. 3)  
**Warning:** The Loctite 495 glue cures easily to skin and fabric. Be aware of residues when attaching the part.
3. Place the Protection Pad on top of the carbon blade so its indents are adjacent to the shims on the foot (Fig. 4)
4. Apply pressure on the Protection Pad for 2 minutes or until the glue has dried (Fig. 5)

5. Cut the middle of the Protection Pad along the blade gap (**Fig. 6**)
6. The leftover from the Protection Pad needs to be removed. Final appearance after the cutting can be seen in (**Fig. 7**)
7. Attach the device to the prosthetic foot by applying a medium strength threadlocker to the Bolt and tighten to 5 Nm (**Fig. 8**)
8. Device placement depends on the foot size and category. To achieve optimum device placement the device must rest on the Protection Pad or with a distance up to 0.5 mm (**Fig. 9**)

Tube connection:

- i. Connect the Tube between the Inlet Port and the device valve on the close fitted socket.
- ii. Cut off any excess tube, if necessary.
- iii. Secure the Tube to the prosthetic foot medially using the Tube Clamp to avoid damage during ambulation.

**Caution:** Ensure that the device Tubing is not compressed or folded by the application of cosmetic finishing (**Fig. 10**)

**Warning:** A loose or dangling Tube can create a tripping hazard. Secure the Tube to the prosthetic foot in a permanent manner.

## USAGE

### *Cleaning and care*

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

**Caution:** Do not use solvents stronger than isopropyl alcohol as it might degrade the material or get through the valve and into the socket.

### *Environmental Conditions*

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

The device can be used in temperatures between -15°C to 50°C.

## MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity. Recommended interval is every 6 months.

The following components might need cleaning or replacing during the lifetime of the device.

- Duckbill valves
- Outlet Port
- Tube

**Note:** The Duckbill valves must be oriented correctly as instructed in (**Fig. 11**) for correct function. Torque the Outlet Port and Inlet Port to 4 Nm.

### *Disassembly*

Remove the Tube, Outlet Port and Inlet Port to access the Duckbill valves.

## Reassembly

Assemble the parts in reverse order.

## Troubleshooting Guide

In case the device is not functioning properly check the following:

Problem	Possible Cause	What to do
Noise from the Unity product. (clicking, tapping, squeaking, cratching, whistling)	The unity has become loose.	Align the Unity accordingly and reattach. Apply medium strength threadlocker and torque the bolt to 5 Nm (Fig. 8).
	Sand is between bumper and the blade.	Clean according to cleaning instructions.
	Sand and/or debris is between the blades.	Clean according to cleaning instructions.
	Duckbill is worn, torn or dirty.	Change both of the Duckbills (Inlet and Outlet Port).
	The filter is dirty.	Change outlet port.
Lost vacuum	Seal-In is leaking.	Visually check the seal and/or liner for damage, change if necessary. Use syringe to check if the vacuum is holding.
	Unity Socket Valve leaking.	To check for this, seal the valve on the inside of the socket with a silicone pad and use a syringe to check if the valve works as intended.
	The Tube is not fully inserted.	Make sure the tube is fully inserted onto the Inlet Port and the Unity Socket Valve.
	Duckbills are incorrectly installed.	Inspect if the Duckbills are oriented in the correct direction and change it if needed (Fig. 11).
	The socket is leaking.	Use the provided syringe and the associated equipment to check the vacuum in the socket.
Noise and/or lost vacuum	The yellow warning label has not been removed from the pyramid of the foot.	Remove the yellow warning label, torque according to specifications (see the adapter instructions for use).
	The female pyramid is incorrectly connected to the male foot pyramid.	Readjust, use medium strength threadlocker and torque correctly (see adapter instructions for torque value).
	The Tube is leaking.	Change the Tube and fit according to instructions.
	Inlet and/or Outlet Port have become loose.	Apply medium strength threadlocker and torque the Inlet and/or Outlet Port to 4 Nm. Make sure the Tube is correctly installed upon the Inlet Ports (Fig. 11).
	Spectra sock restricts the functional performance of the product.	Refit the Spectra sock and make sure it is fitted accordingly (Fig. 10)

## **REPORT OF SERIOUS INCIDENT**

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

## **DISPOSAL**

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

## **LIABILITY**

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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