






Instructions for Use

UNITY[®] FOR PRO-FLEX[®] ST

		3
EN	Instructions for Use	4
CS	Návod k použití	8
BG	Указания за употреба	12
HR	Upute za uporabu	17
HU	Használati javallat	21
RO	Indicații de utilizare	25
SK	Indikácie použitia	29
SL	Indikácie pre použitie	33

DESCRIPTION

The device is a noninvasive mechanical vacuum pump intended for vacuum assisted prosthetic suspension.

Suspension is obtained through the act of a pump which draws air out from the socket in each step creating an elevated vacuum inside the socket interface.

Components (Fig. 1)

1. Membrane Housing
2. Inlet Port
3. Outlet Port
4. Duckbill
5. Bolt
6. Elastic Tube Fixation
7. Tube
8. Tube Clamp

INTENDED USE

The device is intended as a part of a prosthetic system that provides suspension via subatmospheric pressure.

Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

It is intended to be used with the Pro-Flex ST, a respective foot cover, and a Unity Valve Kit for connection to a suspension socket with an Iceross® Seal-In Liner.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS

1. Remove the plastic cover from the foot pyramid (**Fig. 2**). There is a small recess distally to insert a flat tool to ease removal of the plastic cover.
2. Attach the device to the prosthetic foot by applying a medium strength threadlocker to the Bolt and tighten to 5 Nm (**Fig. 3**)
3. Place the elastic tube fixation on the device so it covers the bolt (**Fig. 4**).
4. Insert the tube through the elastic tube fixation.

Tube connection:

- i. Connect the tube to the Inlet Port on the Unity. Connect the other end of the tube to the device valve on the close fitted socket.
- ii. Cut off any excess tube, if necessary.
- iii. Secure the tube to the prosthetic foot medially using the Tube Clamp to avoid damage during ambulation.

Caution: Ensure that the device tubing is not compressed or kinked (Fig. 5).

Caution: Ensure that the device tubing is not compressed or folded by the application of cosmetic finishing (Fig. 6).

Warning: A loose or dangling tube can create a tripping hazard. Secure the tube to the prosthetic foot in a permanent manner.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Caution: Do not use solvents stronger than isopropyl alcohol as it might degrade the material or get through the valve and into the socket.

Environmental Conditions

The device is Waterproof.

A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.

It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.

It can also tolerate occasional exposure to sand, dust, and dirt.

Continuous exposure is not allowed.

Dry with a cloth after contact with fresh water or humidity.

Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

The device can be used in temperatures between -15°C to 50°C.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Recommended interval is every 6 months.

The following components might need cleaning or replacing during the lifetime of the device.

- Duckbill valves
- Outlet Port
- Tube

Note: The Duckbill valves must be oriented correctly as instructed in (Fig. 7) for correct function. Torque the Outlet Port and Inlet Port to 4 Nm.

Disassembly

Remove the Tube, Outlet Port and Inlet Port to access the Duckbill valves.

Reassembly

Assemble the parts in reverse order.

Troubleshooting Guide

In case the device is not functioning properly check the following:

Problem	Possible Cause	What to do
Noise from the Unity product. (clicking, tapping, squeaking, cratching, whistling)	The unity has become loose.	Align the Unity accordingly and reattach. Apply medium strength threadlocker and torque the bolt to 5 Nm (Fig. 3).
	Sand is between bumper and the blade.	Clean according to cleaning instructions.
	Sand and/or debris is between the blades.	Clean according to cleaning instructions.
	Duckbill is worn, torn or dirty.	Change both of the Duckbills (Inlet and Outlet Port). Before inserting the Duckbills, add isopropanol in the inlet and outlet holes of the housing for cleaning.
	The filter is dirty.	Change outlet port.
Lost vacuum	Seal-In is leaking.	Visually check the seal and/or liner for damage, change if necessary. Use syringe to check if the vacuum is holding.
	Unity Socket Valve leaking.	To check for this, seal the valve on the inside of the socket with a silicone pad and use a syringe to check if the valve works as intended.
	The Tube is not fully inserted.	Make sure the tube is fully inserted onto the Inlet Port and the Unity Socket Valve.
	Duckbills are incorrectly installed.	Inspect if the Duckbills are oriented in the correct direction and change it if needed (Fig. 7).
	The socket is leaking.	Use the provided syringe and the associated equipment to check the vacuum in the socket.

Problem	Possible Cause	What to do
Noise and/or lost vacuum	The yellow warning label has not been removed from the pyramid of the foot.	Remove the yellow warning label, torque according to specifications (see the adapter instructions for use).
	The female pyramid is incorrectly connected to the male foot pyramid.	Readjust, use medium strength threadlocker and torque correctly(see adapter instructions for torque value).
	The Tube is leaking.	Change the Tube and fit according to instructions.
	Inlet and/or Outlet Port have become loose.	Torque the Inlet and/or Outlet Port to 4 Nm. Make sure the Tube is correctly installed upon the Inlet Port.
	Spectra sock restricts the functional performance of the product.	Refit the Spectra sock and make sure it is fitted accordingly (Fig. 6)

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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