







Instructions for Use

UNLOADER ONE® RATCHET
UNLOADER ONE® RATCHET
CUSTOM

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	7
FR	Notice d'utilisation	10
ES	Instrucciones para el uso	13
IT	Istruzioni per l'uso	15
NO	Bruksanvisning	18
DA	Brugsanvisning	21
SV	Bruksanvisning	23
EL	Οδηγίες Χρήσης	26
FI	Käyttöohjeet	28
NL	Gebruiksaanwijzing	31
PT	Instruções de Utilização	34
JA	取扱説明書	36
ZH	中文说明书	39
KO	사용 설명서	41

INTENDED USE

The device is intended for unicompartmental unloading of the knee
The device must be fitted and adjusted by a healthcare professional.

Indications for use

- Mild to severe unicompartmental knee osteoarthritis
- Degenerative meniscal tears
- Other unicompartmental knee conditions that may benefit from unloading such as:
 - Articular cartilage defect repair
 - Avascular necrosis
 - Tibial plateau fracture
 - Bone marrow lesions (bone bruises)

No known contraindications.

Warnings and Cautions:

- Regular healthcare professional supervision is recommended for patients with peripheral vascular disease, neuropathy, and sensitive skin.
- Make sure the device fits properly to minimize the possibility of skin irritation. Increase usage time gradually as the skin adapts to the device. If redness appears, temporarily decrease the usage time until it has subsided.
- Care should be taken not to overtighten the device.
- Make sure the device fits properly to achieve effective pain relief.
- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, excessive pressure, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

Device Application

1. Open Upper (A) and Lower (B) Buckles. Ask patient to sit down and extend leg while fitting the device. Ensure that the Dosing Ratchet (C) set to the “1” position. Place the device on the patient’s leg with the Hinge (D) on the affected side of the knee.

Ensure proper alignment of the device on the leg (**Fig. 2**)

- Height positioning (2a): Align Össur logo on the Hinge with upper edge of the patella.
 - Side positioning (2b): Front edge of the Hinge cap should be at midline of leg.
2. Fasten the Buckle buttons to their color-matched Keyholes (E, F).
 - Place the blue Lower Buckle button into the blue Calf Shell Keyhole (E) above Buckle Stability Shelf (G) and use the palm of the hand to snap the Lower Buckle closed (**Fig. 3**). Adjust Calf Strap (H) to appropriate length by tensioning around the calf and folding into the Alligator Clip (I) so that it keeps the device securely and correctly positioned on leg.
 - Bend the patient's knee to 80°. Place the yellow Upper Buckle button into the yellow Thigh Shell Keyhole (F) and use the palm of the hand to snap the Upper Buckle closed (**Fig. 4**). Adjust Thigh Strap (J) to appropriate length by tensioning around leg and folding into Alligator Clip.
 3. Adjust the length of the Dynamic Force System™ (DFS) Straps (K, L).
 - With the patient's knee fully extended (**Fig. 5**), adjust the Upper DFS Strap (K) length until it sits firmly against the leg, and then fold it into Alligator Clip. At this point, patient should not experience any tension or unloading.
 - Bend the patient's knee to 80° (**Fig. 6**). Adjust the Lower DFS Strap (L) length until it sits firmly against the leg, and then fold it into Alligator Clip.
 4. Ask patient to bend the knee with the foot flat on the floor. Pull the Dosing Ratchet to the "5" position (**Fig. 7**). Have the patient stand up and take a few steps to verify correct positioning of the device and tightness of the straps.
 - Determine optimal DFS Strap tension based on patient's pain relief feedback.
 - If patient requires more or less tension with Dosing Ratchet at the "5" position, adjust the length of the DFS Straps accordingly.
 - Aim for final Dosing Ratchet setting at the "5" position since this will give patient the ability to adjust dosing during daily life activities.
 5. When final fit has been confirmed, trim the straps to the appropriate length starting with the Calf Strap so that device sits correctly on leg while trimming other straps.
 - Ensure that the Strap Pad (M) is not wrinkled and positioned where the DFS straps cross in the popliteal fossa (**Fig. 8**).
 - Trim back the straps sufficiently so that the alligator clips are positioned away from the popliteal area. This reduces the bulk behind the knee.

Device Removal

1. Ask patient to sit down with the leg extended.
2. Press button on Dosing Ratchet and place it to the "0" position to release tension on the Upper DFS Strap (**Fig. 7**).
3. Bend the patient's knee to 90° and open both Lower and Upper Buckles.
4. Pull Buckle buttons out of Keyholes.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

Hinge

- Remove foreign materials (e.g., dirt or grass) and clean using fresh water.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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