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Quellen

1. Deutsche Gesellschaft für Orthopädie und Unfallchirurgie, Nachbehandlungsempfehlungen 2023, 9. Aufl., Sektion Physikalische Therapie und Rehabilitation der DGOU.
2. Bahr R, Bahr IA (1997) Incidence of acute volleyball injuries: a prospective cohort study of injury mechanisms and risk factors. *Scand J Med Sci Sports* 7:166–171
3. Babins EM (2012) Lace-up ankle braces reduced acute ankle injuries in high school basketball players. *Clin J Sport Med* 22(4):379–380
4. Balduini FC, Tetzlaff J (1982) Historical perspectives on injuries of the ligaments of the ankle. *Clin Sports Med* 1:3–12
5. Best R, Brüggemann P, Petersen W, Rembitzki I, Ellermann A, Gösele-Koppenburg A, Liebau C (2011) Aktuelle und neue Konzepte in der Behandlung akuter Außenbandverletzungen des Sprunggelenkes, *Deutsche Zeitschrift für Sportmedizin* Jahrgang 62, Nr. 3:57–62
6. Beynon BD, Renström PA, Haugh L, Uh BS, Barker H (2006) A prospective, randomized clinical investigation of the treatment of first-time ankle sprains. *Am J Sports Med* 34(9):1401–1412
7. Bleakley CM, O'Connor SR, Tully MA, Roche LG, Macauley DC, Bradbury I, Keegan S, McDonough SM (2010) Effect of accelerated rehabilitation on function after ankle sprain: randomised controlled trial. *BMJ* 10(340):c1964
8. Boyce SH, Quigley MA, Campbell S (2005) Management of ankle sprains: a randomised controlled trial of the treatment of inversion injuries using an elastic support bandage or an Aircast ankle brace. *Br J Sports Med* 39(2):91–96
9. Cass JR, Settles H (1994) Ankle instability: in vitro kinematics in response to axial load. *Foot Ankle Int* 15:134–140
10. Cooke MW, Marsh JL, Clark M, Nakash R, Jarvis RM, Hutton JL, Szczepura A, Wilson S, Lamb SE; CAST trial group (2009) Treatment of severe ankle sprain: a pragmatic randomised controlled trial comparing the clinical effectiveness and cost-effectiveness of three types of mechanical ankle support with tubular bandage. The CAST trial. *Health Technol Assess* 13(13):iii, ix, x, 1–121
11. Chan KW, Ding BC, Mroczek KJ (2011) Acute and chronic lateral ankle instability in the athlete. *Bull NYU Hosp Jt Dis* 69(1):17–26
12. Costa ML, Achten J, Wagland S, Marian IR, Maredza M, Schlüssel MM, Liew AS, Parsons NR, Dutton SJ, Kearney RS, Lamb SE, Ollivere B, Petrou S (2020) Plaster cast versus functional bracing for Achilles tendon rupture: the UKSTAR RCT. *Health Technol Assess*;24(8):1-86
13. Dubin JC, Comeau D, McClelland RI, Dubin RA, Ferrel E (2011) Lateral and syndesmotic ankle sprain injuries: a narrative literature review. *J Chiropr Med* 10(3):204–219
14. Ekstrand J, Gillquist J (1983) Soccer injuries and their mechanism: a prospective study. *Med Sci Sports Exerc* 15:367–370
15. Engebretsen AH, Myklebust G, Holme I, Engebretsen L, Bahr R (2010) Intrinsic risk factors for acute ankle injuries among male soccer players: a prospective cohort study. *Scand J Med Sci Sports* 20(3):403–410
16. Fong DT, Chan YY, Mok KM, Yung PSh, Chan KM (2009) Understanding acute ankle ligamentous sprain injury in sports. *Sports Med Arthrosc Rehabil Ther Technol* 30(1):14
17. Freeman MA (1965) Instability of the foot after injuries to the lateral ligament of the ankle. *J Bone Joint Surg Br* 47(4):669–677
18. Freeman MA, Dean MR, Hanham IW (1965) The etiology and prevention of functional instability of the foot. *J Bone Joint Surg Br* 47(4):678–685
19. Gerber JP, Williams GN, Scoville CR, Arciero RA, Taylor DC (1998) Persistent disability associated with ankle sprains: a prospective examination of an athletic population. *Foot Ankle Int* 19(10):653–660
20. Handoll HH, Rowe BH, Quinn KM, de Bie R (2001) Interventions for preventing ankle ligament injuries. *Cochrane Database of Systematic Reviews*. Issue 3, Art. No. CD000018
21. Hertel J (2000) Functional instability following lateral ankle sprain. *Sport Med* 29(5):361–371
22. Holme E, Magnusson SP, Becher K (1999) The effect of supervised rehabilitation on strength, postural sway, position sense and re-injury risk after acute ankle ligament sprain. *Scand J Med Sci Sports* 9(2):104–109
23. Houglum PA (1992) Soft tissue healing and its impact on rehabilitation. *J Sport Rehabil* 1:19–23
24. Hupperets MD, Verhagen EA, van Mechelen W (2009) Effect of unsupervised home based proprioceptive training on recurrences of ankle sprain: randomised controlled trial. *BMJ* 9(339):b2684. doi:10.1136/bmj.B2684
25. Junge A, Engebretsen L, Mountjoy ML, Alonso JM, Renström PA, Aubry MJ, Dvorak J (2009) Sports injuries during the Summer Olympic Games 2008. *Am J Sports Med* 37(11):2165–2167
26. Kerkhoffs GMMJ, Rowe BH, Assendelft AJJ, Kelly K, Struis PAA, van Dijk CN (2002) Immobilisation and functional treatment for acute lateral ligament injuries in adults (Review) *Cochrane Database Syst Rev* Issue 3, Art. No. CD003762. doi:10.1002/14651858.CD003762
27. Kerkhoffs GM, Struijs PA, Marti RK, Assendelft WJ, Blankevoort L, van Dijk CN (2002) Different functional treatment strategies for acute lateral ankle ligament injuries in adults. *Cochrane Database Syst Rev* 3:CD002938
28. Kerkhoffs GMMJ, Handoll HHG, de Bie R, Rowe BH, Struis PAA (2007) Surgical vs. conservative treatment for acute injuries of the lateral ligament complex in adults (Review). *Cochrane Database Syst Rev* Issue 2, Art. No. CD000380. doi:10.1002/14651858.CD000380
29. Collaborative, B. (2019). Weight-bearing in ankle fractures: An audit of UK practice. *Foot (Edinburgh, Scotland)*, 39, 28–36.
30. Bruce, J., & Sutherland, A. (2013). Surgical versus conservative interventions for displaced intra-articular calcaneal fractures. *The Cochrane Database of Systematic Reviews*, 25(1), CD008628. <http://doi.org/10.1002/14651858.CD008628.pub2>
31. De Boer, A. S., Van Lieshout, E. M. M., Van Moolenbroek, G., Hartog, Den, D., & Verhofstad, M. H. J. (2018). The effect of time to post-operative weightbearing on functional and clinical outcomes in adults with a displaced intra-articular calcaneal fracture; A systematic review and pooled analysis. *Injury*, 49(4), 743–752. <http://doi.org/10.1016/j.injury.2018.02.021>.
32. Pituckanotai, K., Arirachakaran, A., Piyapittayanun, P., Tuchinda, H., Peradhammanon, E., & Kongtharvonskul, J. (2018). Comparative Outcomes of Cast and Removable Support in Fracture Fifth Metatarsal Bone: Systematic Review and Meta-Analysis. *The Journal of Foot and Ankle Surgery: Official Publication of the American College of Foot and Ankle Surgeons*, 57(5), 982–986.
33. El-Akkawi, A. I., Joanroy, R., Barfod, K. W., Kallelose, T., Kristensen, S. S., & Viberg, B. (2018). Effect of Early Versus Late Weightbearing in Conservatively Treated Acute Achilles Tendon Rupture: A Meta-Analysis. *The Journal of Foot and Ankle Surgery: Official Publication of the American College of Foot and Ankle Surgeons*, 57(2), 346–352. <http://doi.org/10.1053/j.jfas.2017.06.006>.
34. Kastoft, R., Bencke, J., Speedtsberg, M. B., Penny, J. Ø., & Barfod, K. (2018). Early weight-bearing in nonoperative treatment of acute Achilles tendon rupture did not influence mid-term outcome: a blinded, randomised controlled trial. *Knee Surgery, Sports Traumatology, Arthroscopy: Official Journal of the ESSKA*, 96(6), 1497–8. <http://doi.org/10.1007/S00167-018-5058-4>.
35. AWMF-Leitlinien Unfallchirurgie; AWMF: 012-022, ICD S-93.40, S-93.41, S-93.42, S-93.43.
36. Erfahrungswerte aus der klinischen Praxis.