










Instructions for Use




AeroFit® Seal-In

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	15
IT	Istruzioni per l'uso	19
NO	Bruksanvisning	23
DA	Brugsanvisning	26
SV	Bruksanvisning	30
EL	Οδηγίες Χρήσης	33
FI	Käyttöohjeet	37
NL	Gebruiksaanwijzing	40
PT	Instruções de Utilização	44
PL	Instrukcja użytkowania	48
TR	Kullanım Talimatları	51
RU	Инструкция по использованию	55
JA	取扱説明書	59
ZH	中文说明书	62
KO	사용 설명서	65
CS	Návod k použití	69
SK	Indikácie použitia	72
RO	Indicații de utilizare	76
HR	Upute za uporabu	79
HU	Használati javallat	83
BG	Указания за употреба	86
SL	Indikácie pre použitie	90
UA	Інструкція з використання	94
AR	مرادختس ال ا ت امريلعت	97

DESCRIPTION

The device is a transfemoral silicone liner that covers the residual limb and provides an interface between the residual limb and the socket of the prosthesis. The device forms a suspension system where a distal vacuum chamber is created below the level of the seal that reliably suspends the residual limb in the prosthesis.

The device features a breathable silicone structure designed to reduce the build-up of moisture within the device during use and is used solely with the AeroFit® Socket.

INTENDED USE

The device is a prosthetic interface with suspension properties intended to be used as part of a system that replaces a missing lower limb. When used in combination with the socket, the breathable nature of the device reduces moisture on the residual limb. Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.

The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
 - No known contraindications
- The device is for low to high impact use, e.g., walking and occasional running.

The device should not be used by patients with the following conditions:

- Poor cognitive abilities that affect the proper use of the device.
- Drastic volume change within a short period of time, exceeding 6 ply.
- An extremely short residual limb (See. **Table 1**).

Caution: The device might not be suitable for patients with deep invaginations (scars) that pass the seal area (between the distal end and the breathable section).

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.

Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

Misuse of the device may result in loss of suspension.

Potential side effects are skin irritation/reaction problems.

The device must not come in contact with glass, carbon fibers, or other abrasive particles. Washing alone may not be sufficient to eliminate the problem. If the device is exposed to abrasive substances or chemicals, it should be returned to a healthcare professional.

The device is for single patient use.

DEVICE SELECTION

To determine the correct size:

- Measure the circumference of the residual limb 4 cm from the distal end (**Fig. 1**).
 - Choose the corresponding device size. If the measurement taken does not have a corresponding device size, round the measurement down.
- A too small/tight device can cause excessive pressure on the residual limb and increased pistoning.

A too large/loose device can influence the obtained suspension and cause excessive movement between the device and the residual limb.

Caution: If increased movement, pistoning, numbness, blistering, or rashes are experienced, the patient should stop using the device and contact a healthcare professional immediately.

Refer to **Table 1** for the minimum residual limb length required for each device size.

Table 1 - Minimum Residual Limb Length Relative to Device Size

Device Size	Minimum Residual Limb Length (mm)
23.5-28	200

Device Size	Minimum Residual Limb Length (mm)
30-32	210
34-36	220
38-40	230
45	250
50	260
55	280

Residual limb length is measured from the ischial tuberosity to the distal end, with the tissue hanging down.

FITTING INSTRUCTIONS

Trimming Instructions

Ensure that the area of interior silicone with no micro-hole pattern is aligned to the groin or any area subject to higher pressure or force before trimming (**Fig. 2**).

The device may be trimmed to increase range of motion, or for patient preference (**Fig. 3**). The Iceross® Curvemaster is a recommended tool, specifically designed to produce a smooth, rounded edge when trimming.

If not available, use sharp scissors to trim the device to the desired shape and length.

When trimming the device, ensure that no sharp edges or notches are created, as these may cause tears in the material or lead to skin irritation (**Fig. 4**).

Do not trim the device below the proximal edge of the socket brim.

Ensure the device is not overlength to prevent the proximal edge from rubbing on the groin area.

Socket Fabrication Guidelines

Refer to the AeroFit Socket Instructions for Use for the fabrication of the socket.

The device is Unity® compatible.

Socket Fitting

Achieving the intended function of the device is dependent on the fit of the socket. Appropriate socket design is essential to achieve suspension, control, and comfort.

The following checklist should be used when checking the socket fit:

- Patient can don the socket easily without the use of lubricant aids.
- Full contact between the seal lips and the socket inner wall.
- No air pockets develop when the patient loads the prosthesis.
- No pistoning in the socket
- No visible movement when pulling on the prosthesis.
- Full movement possible without compromising the suspension.

If feedback indicates any of the above conditions are not met, the socket will need to be adjusted or remade accordingly.

USAGE

Donning

Do not apply lotion to the residual limb immediately before donning the device. The residual limb should be clean and dry to prevent device damage.

Caution: Open wounds or damaged skin should be covered with a bandage or other appropriate covering to prevent direct contact between the wound and device.

Liner Donning

The proximal area of the device includes a region without a micro-hole pattern on the interior silicone surface. This region should be aligned with the groin or with any area subject to higher pressure or force to help reduce the risk of skin irritation (**Fig. 2**).

Caution: Follow the liner donning instructions provided carefully. Incorrect donning can result in skin reactions, such as blisters and rashes.

1. Grip the top of the device from the inside and slide it over the hand until the inner surface is fully exposed (**Fig. 5**). Take care not to damage the device with fingernails.
2. After exposing as much of the distal end of the device as possible, position it against the distal end of the residual limb (**Fig. 6**).
3. With light compression, roll the device all the way upwards onto the limb (**Fig. 7**).
Do not tug or pull the upper end of the device when rolling it upwards onto the limb.
4. Check that no air pockets are present. If air pockets are found, reapply the device.
5. Ensure the seal sits correctly and horizontally around its full circumference (**Fig. 8**).

6. Ensure the seal is correctly positioned by running fingers around the seal and assessing the tension.

Socket Donning

Do not use lubricant spray as a donning aid for the device, as this can destroy the low-friction coating, making socket donning difficult. With the device (and seal ring) correctly configured, step into the socket, expelling all air until fully engaged into the distal end (**Fig. 9**).

Cleaning and Care

Daily cleaning of the residual limb is essential. Use of a mild, pH balance, 100% fragrance- and dye-free liquid soap is recommended. If the skin is dry, apply a pH balance, 100% fragrance- and dye-free lotion to nourish and soften the skin.

Warning: The device and residual limb must be cleaned and inspected daily. Before first use, the device must be cleaned and thoroughly inspected. The device must be fully dried before use.

Machine Wash

- Place the device in the provided wash bag with the textile facing outwards.
- Machine wash the device at 40°C with a mild detergent and at a low spin speed.
- Allow the device to air dry using the liner stand.

The device must be machine washed after a maximum of 3 days of use. Do not tumble dry, iron, bleach, or wash with fabric softener.

Daily Cleaning After Use

- Turn the device inside out (**Fig. 5**).
- Use a damp cloth or a wipe with water or a mild soap solution to clean the exposed interior surface.
- Either allow the device to air dry using the liner stand or pat dry with a lint-free cloth.
- Always return the device to its neutral state immediately after cleaning. Do not store inside-out as device damage may result.
- The device should not be exposed to excessive heat or sunlight.
- The seal must be clean and free from dirt before use.

Skin health and the device should be monitored daily. In case of abnormal symptoms, discontinue use of the device and contact the healthcare professional immediately.

Be aware that common household or bath products may cause or contribute to skin irritation, e.g., soaps, deodorants, perfumes, abrasive cleaners, aerosol- or alcohol sprays.

Caution: When the device is used in circumstances or environments where excessive perspiration can be expected, such that stability or suspension are compromised, it is recommended to remove the device and dry off both the device and the residual limb.

Warning: After contact with salt and/or chlorinated water, wash the device with tap water.

Environmental Conditions

The device is safe to use with caution with waterproof components.

Warning: If the device is submerged in water, remove it and ensure that the residual limb, the device's skin contact layer, and the socket are dry before ambulation.

Caution: The suspension is compromised if water gets between the skin and the device or the seal and the socket.

Caution: Frequent exposure to salt and/or chlorinated water might affect the durability of the device.

SUPPLEMENTARY DEVICES

The following accessories can be used in conjunction with Össur seal liners. Refer to Össur catalog for more information, e.g., sizing.

Iceross® Pads

The device is a prosthetic interface intended to be used as part of a system that replaces a missing lower limb.

The pad can be used for prosthetic volume management to temporarily accommodate lost distal end volume as well as to assess distal contact while fitting a check socket.

Place the pad in the distal end of the socket. If placed in the distal end of the liner, it may compromise suspension.

Note: For daily care, refer to Cleaning and Care instructions.

Iceross Seal-In® Socks

The device is intended to be used for accommodating residual limb volume reduction.

If the patient experiences insecurity or decreased contact with the socket due to volume fluctuations, one or more socks can be pulled over the liner to restore socket fit.

1. Invert the seal completely in the distal direction (**Fig. 10**).
2. Pull the sock over the inverted seal until the distal end of the sock is level with the proximal border of the seal (**Fig. 11**).
3. Reflect the seal over the end of the sock (**Fig. 12**).
4. Don the prosthesis as normal.

Caution: Loss of suspension may be encountered if a sock is placed over the seal.

Machine wash at max 40°C with light-colored fabric. Do not tumble dry.

Caution: If volume fluctuation cannot be adequately compensated for by this means, the socket will need to be remade to suitably reduced dimensions.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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