## ENGLISH

**MD** Medical Device

### INDICATIONS FOR USE

For conditions that may benefit from restricting motion of the wrist, such as:

- Mild to severe sprains and strains
- Tendonitis
- Post-cast healing
- Other selected soft tissue injuries

### WARNINGS AND PRECAUTIONS

- This product is designed to be worn directly on the skin. Inspect your hand, wrist, and forearm before and after each use.
- Never wear the brace directly over an open wound.
- This device is for single patient use.
- Do not overtighten the brace, risking loss of blood flow. Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease need to be extra careful and are advised to consult with a healthcare professional before using the brace. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this brace, discontinue the use of it immediately and consult with a healthcare professional.
- If there is a change or loss in device functionality or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

### FITTING INSTRUCTIONS

**Positioning of the brace:**  
Align the palmer stay to the palm of the hand (**Figure 1**).

### Putting on the brace:

- Undo all contact closure straps (**Figure 2**)

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- Slide hand into the brace and place thumb into thumb hole. Ensure stockinette is positioned on top of the hand (**Figure 3**)
- Fasten all contact closure straps and secure to a snug fit (**Figure 4 and 5**).

- Adjustments may be necessary as hand, wrist, and thumb swelling patterns change.
- Straps may be cut as desired (**Figure 6**)

### WARNINGS AND PRECAUTIONS

- This product is designed to be worn directly on the skin. Inspect your hand, wrist, and forearm before and after each use.
- Never wear the brace directly over an open wound.
- This device is for single patient use.
- Do not overtighten the brace, risking loss of blood flow. Patients with already compromised blood circulation in their extremities, like those with diabetes or peripheral vascular disease need to be extra careful and are advised to consult with a healthcare professional before using the brace. If you experience pain, swelling, sensation changes, or if your extremity shows signs of insufficient blood flow (turns blue, white, or cold) while using this brace, discontinue the use of it immediately and consult with a healthcare professional.
- If there is a change or loss in device functionality or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

### REPORT A SERIOUS INCIDENT

**Important notice to users and/or patients established in Europe:**

The user and/or patient must report any serious incident that has occurred in relation to the device to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

### DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

### LIABILITY

Ossur does not assume liability if the device is not maintained as instructed by the instructions for use, if the device is assembled with components from other manufacturers, or if the device is used outside of recommended use condition, application, or environment.

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