




Instructions for Use



MCPDriver®
BPF 202 | Gen 3


		3
EN	Instructions for Use	6
DE	Gebrauchsanweisung	10
FR	Notice d'utilisation	14
ES	Instrucciones para el uso	19
IT	Istruzioni per l'uso	23
NO	Bruksanvisning	28
DA	Brugsanvisning	32
SV	Bruksanvisning	37
EL	Οδηγίες Χρήσης	42
FI	Käyttöohjeet	47
NL	Gebruiksaanwijzing	51
PT	Instruções de Utilização	55
PT-BR	Instruções de Uso	60
PL	Instrukcja użytkowania	64
TR	Kullanım Talimatları	69
JA	取扱説明書	73
ZH	中文说明书	78
KO	사용 설명서	81
CS	Návod k použití	85
SK	Indikácie použitia	90
RO	Indicații de utilizare	94
HR	Upute za uporabu	98
HU	Használati javallat	103
BG	Указания за употреба	107
SL	Indikácie pre použitie	112
ET	Kasutusjuhend	116
UA	Інструкція з використання	120
HE	הוראות שימוש	125
AR	تعليمات الاستخدام	129


5


6a


6b


6c


6d


6e


6f


6g

6h

LIST OF PARTS (Fig. 1&2)

- a. Wrist Alligator Tabs
 - b. Backplate
 - c. Abduction/Adduction Screw
 - d. Metal Linkages
 - e. Ring Adjustment Screw & Slider
 - f. P2 Silicone Pad
 - g. Device Tip
 - h. Hinge Joint
 - i. Backplate Pad
 - j. Suspension Ring For more information: visit ossur.com
 - k. P2 Fairing *
 - l. Silicone Tip Pad
 - m. Shim Holes (optional)
 - n. Single Articulating Tip (optional)
- * not present in Single Articulating Tips

Click Strap Parts (Fig. 3)

- a. Clasp Strap
- b. Cinch Strap
- c. Clasp Plate
- d. Cinch Plate
- e. Strap Retainer Clip

Classic Strap Parts (Fig. 4)

- a. Wrist O-Ring
- b. Wrist Strap

Palmar Strap Parts (Fig. 5)

- a. 15mm Strap
- b. Y Alligators
- c. APS Steel Guide
 - Single-digit
 - Multi-digit
- d. APS Washer

Components Adjustable by a Certified Prosthetist (Fig. 6)

PRODUCT DESCRIPTION

The device is a body-driven mechanical prosthesis made of titanium linkages, a nylon backplate and finger fairings with a silicone interface. It is anchored to the dorsal surface of the hand and is attached with a wrist strap and nylon ring on the residual finger(s).

Indications for Use

The MCPDriver is designed for a partial-finger limb difference distal to the metacarpophalangeal (MCP) joint when there is digit length to engage the Suspension Ring (Fig. 2j). The MCPDriver is designed to be used in conjunction with a provided Backplate system or a custom socket fabricated by a Certified Prosthetist. MCPDrivers with Single Articulating Tip (SAT) design are used to compensate for excess digit length or limited range of motion at the MCP joint.

Intended Use

The device is a prosthesis that aims to restore hand function to patients with partial-finger limb

difference.

The device is intended to be fit by a certified prosthetist.

Contraindications for Use

The residual digit must be properly healed before using device. Do not use device if open wounds or lesions are present.

Device Application

1. Place the device on the hand by adding the ring to the appropriate finger(s) and positioning the backplate on the back of hand.
 - Align the Hinge Joint (**Fig. 2h**) directly above, or slightly proximal to the MCP joint of the patient's hand.
 - The Abduction/Adduction Screw (**Fig. 1c**) should be $\frac{1}{2}$ " – $\frac{3}{4}$ " proximal to the center of the MCP joint.
2. Check articulation of device.
 - The Suspension Ring (**Fig. 2j**) should fit loosely, allowing the ring to slide along the digit(s) during active range of motion.
 - Shorter residual fingers may benefit from a smaller sized ring or silicone shim to help capture range of motion.
 - The suspension ring should never bind on the tissue of the treated digit.
 - If volume change is experienced, use shims to achieve optimal fit.
3. Select appropriate wrist strap option for the wearer.
 - Determine length of strap required by measuring the wrist circumference between the alligator tabs.
 - Position the Clasp Plate (**Fig. 3c**) or O-Ring (**Fig. 4a**) at a comfortable position on the wrist and trim ends of the strap to length. Ensure to only cut the loose ends of the straps. The classic strap O-Ring can be positioned on either side for ease of donning.
 - Re-install the strap into the alligator tabs (**Fig. 1a**).
 - For Click Strap Only: Determine a comfortable position for the Strap Retainer Clip (**Fig. 3e**). if you need to remove length, do not cut the retainer clip but remove excess length from the loose end.
4. Fasten straps.
 - The wrist strap should be distal to the wrist crease in a position which does not block wrist motion.
 - Some devices will come with an additional Palmar Strap (**Fig. 5**) to help with suspension. The Palmar strap should go across the palm in a comfortable position.
 - If the MCPDriver migrates proximally or significantly impedes wrist extension, reposition the backplate distally on the dorsum of the hand and reattach the wrist strap at an angled or offset position (**Fig. 6g**).
 - Straps should not be so tight that motions of the fingers or wrist are restricted or painful.
5. (For MCPDrivers with an SAT Tip)
 - Adjust the tip angle (**Fig. 6b**). A flexed tip may be used for smaller items or to compensate for lack of range of motion.

Usage Recommendations

- Occupational or physical therapy is highly recommended to optimize functional use after receiving the device.
- Follow a healthcare provider's direction regarding length of time and use, particularly during the adoption phase.

Caution

- For device wearers with sensitive skin conditions, use the MCPDriver as recommended by a healthcare provider.
- If redness appears on the skin, decrease usage of the device until skin status returns to normal.
- If redness persists or discomfort arises, stop usage of the device immediately and contact a healthcare provider.

- Any concerns regarding fit or function should be addressed by a prosthetist.

Important note: No part of the metal should touch or cause pressure on the skin of the hand or digit, nor should the device make contact with the distal end of the residual digit at any position – if either occurs, immediately stop usage of the device and contact a prosthetist.

Serious incidents should be reported to Össur Iceland ehf. and appropriate government authorities in which the device wearer is established.

Warning: Do not use device in environments or situations that would cause harm to a natural hand.

Included With The Device

- Additional Wrist Strap
- Alternative Suspension Ring Sizes
- Palmar Strap Kit (if requested)
- Silicone Shims (Small, Medium, and Large) for rings with shim holes
- (1) Torx 8-Wing Driver
- (1) .050" Hex Wrench

Device Maintenance (See Fig. 6 for Components Adjustable by a Certified Prosthetist)

- The MCPDriver is water-, dirt-, chemical- and heat-resistant.
- Working regularly with degreasers may strip lubrication from hinges. If necessary, apply a food-safe mineral oil to the hinge locations.
- Contact your prosthetist to facilitate servicing or replacement components.

Replacement Components Available Separately:

- Device Tip
- Suspension Ring
- Palmar Strap
- P2 Fairing
- Strap

Instructions:

- To replace the Device Tip, use a 0.050" Hex Wrench. Access the screw through the hidden port between the Silicone Tip Pad and Device Tip at the tip leading-edge centerline (**Fig. 6a**).
- To replace the tips of an MCPDriver with an SAT design, use a flat-head screw driver to loosen and remove the knurled nut. Remove and replace the new tip, then hand tighten the nut. (**Fig. 6b**).
- To replace the P2 Fairing, use a 0.050" Hex Wrench. Access the screw through the side port (**Fig. 6c**).

Caution: Do not over-tighten. Doing so will increase stiffness of the device.

- Rings and Backplates can be replaced using a Torx 8-Wing Driver (**Fig. 6d and 6e**).
- Shims can be added or removed by inserting them into or pulling them out of the shim holes (**Fig. 6f**).
- Wrist Straps can be replaced by attaching the new strap to the Wrist Alligator Tabs using the Velcro closure (**Fig. 6g**).
- To install the Palmar Strap, remove the metal linkage and abduction/adduction washer. Place APS Steel Guide between backplate and linkage and insert APS Abduction/Adduction washer (**Fig. 6h**). For multi-digit configuration see (**Fig. 6i**). Straps can be trimmed to patient's preference.

Note: LOCTITE® 243 should be applied each time the ring adjustment and abduction/adduction screws are removed. Ensure that the Loctite has appropriate setting time each use (**Fig. 5d and 5e**); refer to Loctite usage instructions.

Important note: For safety reasons, hinges should only be disassembled by the device manufacturer.

Warning: Do not modify device in any way or attempt maintenance other than intended by the manufacturer. Unintended modifications or maintenance may cause device failure.

Care & Cleaning

- Prolonged exposure to strong solvents and degreasers may damage surfaces and remove hinge lubrication.
- Rinse with fresh water after exposure to salt water.
- Clean device with warm water and mild detergent as necessary; pat- or air-dry.
- Scrub dirt or grime using a soft bristle brush, remove dust with compressed air if necessary.
- Do NOT machine wash; Do NOT place in dryer.

Note: Strap components will absorb water and take longer to dry than the rest of the device or system.

Warranty

Össur Iceland ehf. includes a standard 2-year warranty (from ship date) on the mechanical function of the MCPDriver. Please contact Össur Iceland ehf. for additional, extended warranty options.

Change in Performance

In the event of change in performance of a device, always cease using device and contact your clinic.

Warning: Do not use device if contact between device and patient's skin is causing pain, excessive abrasion, or bleeding.

Report of Serious Incident

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

Liability

Össur Iceland ehf. does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device disassembled or amended with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

PRODUCT TECHNICAL SPECIFICATIONS

Materials

- Nylon
- Silicone
- Stainless Steel
- Titanium
- Chemical Resistance
- Bleach
- Ammonia
- Acetone
- Rubbing Alcohol (Isopropanol)
- Gasoline
- Antifreeze
- Ethanol

Operational Temperature

- 0 – 120° F (-17 – 48° C)

Expected Service Life

- 3 years

Force Output

The MCPDriver is designed such that any force the patient can generate at their MCP joint can be transmitted by the devices. Force output is dependent on the patient's strength.

Silicone

Silicone will absorb gasoline or diesel if exposed and will experience a volume change. It will return to a normal shape if cleaned after exposure.

For more information: visit ossur.com

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