
Instructions for Use





ThumbDriver®







BPF203 | Gen 3

		3
EN	Instructions for Use	5
DE	Gebrauchsanweisung	8
FR	Notice d'utilisation	12
ES	Instrucciones para el uso	16
IT	Istruzioni per l'uso	20
NO	Bruksanvisning	24
DA	Brugsanvisning	28
SV	Bruksanvisning	31
EL	Οδηγίες Χρήσης	35
FI	Käyttöohjeet	39
NL	Gebruiksaanwijzing	43
PT	Instruções de Utilização	46
PT-BR	Instruções de Uso	50
PL	Instrukcja użytkowania	54
TR	Kullanım Talimatları	58
JA	取扱説明書	62
ZH	中文说明书	65
KO	사용 설명서	69
CS	Návod k použití	72
SK	Indikácie použitia	76
RO	Indicații de utilizare	80
HR	Upute za uporabu	83
HU	Használati javallat	87
BG	Указания за употреба	91
SL	Indikácie pre použitie	95
ET	Kasutusjuhend	99
UA	Інструкція з використання	102
HE	הוראות שימוש	106
AR	تعليمات الاستخدام	110

5

LIST OF PARTS (Fig. 1, 2)

- a. Anchor Plate
- b. Abduction/Adduction Screw
- c. Thumb Anchor
- d. Metal Linkages
- e. Ring Adjustment Screw & Slider
- f. Device Tip
- g. Anchor Plate Pad
- h. Hinge Joint
- i. Knuckle Pad
- j. Suspension Ring
- k. Tip Angle Adjustment Screw
- l. Silicone Tip Pad
- m. Shim Holes (optional)

Fig. 3: Click Strap Parts: (a) Clasp Strap, (b) Cinch Strap, (c) Clasp Plate, (d) Cinch Plate, (e) Strap Retainer Clip, (f) Dual-Ended Alligator,

Fig. 4: Classic Strap Parts: (a) Wrist O-Ring, (b) Wrist Strap

Fig. 5: Components Adjustable by a Certified Prosthetist

PRODUCT DESCRIPTION

The device is a body-driven mechanical prosthesis made of titanium linkages, a nylon backplate and finger fairings with a silicone interface. It is anchored to the dorsal surface of the hand and is attached with a wrist strap and nylon ring on the residual finger(s).

The device allows for user adjustment of thumb tip flexion or extension.

Indications for Use

The ThumbDriver is designed for a partial-finger limb difference distal at the MP joint or phalanges of the thumb when there is digit length to engage the Suspension Ring (**Fig. 2j**) and motion at the CMC joint. The ThumbDriver is driven off two joints: the metacarpophalangeal (MP) and the carpometacarpal (CMC) joints of the thumb.

Note: The device will only achieve active articulation if the MP joint is intact and functional. However, the device will provide most of its functional opposition via CMC tracking, therefore full ROM at the MP joint is not required for successful device fitting.

Intended Use

The device is a prosthesis that aims to restore hand function to patients with partial-finger limb difference.

The device is intended to be fit by a certified prosthetist.

Contraindications for Use

The residual digit must be properly healed before using device. Do not use device if open wounds or lesions are present.

Device Application

1. Place the device on the hand by adding the ring to the thumb.
 - Align the Hinge Joint (**Fig. 2h**) directly above or slightly proximal to the MCP joint.
 - Align the Anchor Plate (**Fig. 1a**) directly above the intact CMC joint.
2. Check articulation of device.
 - The Suspension Ring (**Fig. 2j**) should fit loosely, allowing the ring to slide along the thumb during active range of motion.

- Shorter residual thumbs may benefit from a smaller sized ring or silicone shim to help capture range of motion.
 - The Suspension Ring should never bind on the tissue of the treated digit.
 - If volume change is experienced, use shims to achieve optimal fit.
3. Select appropriate strap option for the wearer.
 - Determine length of strap required by measuring the wrist circumference accounting for the anchor plate.
 - Position the Clasp Plate (**Fig. 3c**) or O-Ring (**Fig. 4a**) at a comfortable position on the wrist. Remove the Dual-Ended Alligator Tab and trim the loose ends as required.
 - Re-install the strap into the Dual-Ended Alligator (**Fig. 3f**).
 - For Click Strap Only: Determine a comfortable position for the Strap Retainer Clip (**Fig. 3e**). if you need to remove length, do not the cut the retainer clip but remove excess length from the loose end.
 4. Fasten the strap at the base of the hand.
 - The strap should be distal to the wrist crease in a position which does not block wrist motion.
 - If the ThumbDriver migrates proximally or significantly impedes wrist extension, reposition the Anchor Plate distally on the dorsum of the hand and reattach the strap at an angled or offset position.

Usage Recommendations

- Occupational or physical therapy is highly recommended to optimize functional use after receiving this device.
- Follow a healthcare provider's direction regarding length of time and use, particularly during the adoption phase.

Caution

- For device wearers with sensitive skin conditions, use the ThumbDriver as recommended by a healthcare provider.
- If redness appears on the skin, decrease usage of the device until skin status returns to normal.
- If there is persistent redness or discomfort, stop usage of the device immediately and contact a healthcare provider.
- Any concerns regarding fit or function should be addressed by a prosthetist.

Important Note: No part of the metal should touch or cause pressure on the skin of the hand or digit, nor should the device make contact with the distal end of the residual digit at any position – if either occurs, immediately stop usage of the device and contact a prosthetist.

Serious incidents should be reported to Össur Iceland ehf. and appropriate government authorities in which the device wearer is established.

Warning: Do not use device in environments or situations that would cause harm to a natural hand.

Included With The Device

- Additional Wrist Strap
- Alternative Suspension Ring Sizes
- Silicone Shims (Small, Medium, and Large)
- (1) Torx 8-Wing Driver
- (1) .050" Hex Wrench

Device Maintenance (See **Fig. 5** for Components Adjustable by a Certified Prosthetist)

- The ThumbDriver is water-, dirt-, chemical- and heat-resistant.
- Working regularly with degreasers may strip lubrication from hinges. If necessary, apply a food-safe mineral oil to the hinge locations.
- Contact your prosthetist to facilitate servicing or replacement components.

Replacement Components Available Separately:

- Device Tip
- Knuckle Pad
- Strap
- Suspension Ring
- Twisted Tip (optional)

Instructions:

- To replace the Device Tip, use a flat-head screw driver to loosen and remove the Tip Angle Adjustment Screw. Remove and replace the tip. Re-tighten the new Device Tip until the teeth embed into the plastic (**Fig. 5a**).

Caution: Do not over-tighten.

- Rings, Anchor Plates, and Knuckle Pads can be replaced using a torx 8-wing driver (**Fig. 5b, 5c, and 5f**).
- Shims can be added or removed by inserting them into or pulling them out of the shim holes (**Fig. 5d**).
- Straps can be replaced by attaching new strap to the Wrist Alligator Tabs using the Velcro closure (**Fig. 5e**).

Note: LOCTITE® 243 should be applied each time the ring adjustment and abduction/adduction screws are removed.

Ensure that the Loctite has appropriate setting time each use (**Fig. 3b and 3c**); refer to Loctite usage instructions.

Important Note: For safety reasons, hinges should only be disassembled by the device manufacturer.

Warning: Do not modify device in any way or attempt maintenance other than intended by the manufacturer. Unintended modifications or maintenance may cause device failure.

Care & Cleaning

- Prolonged exposure to strong solvents and degreasers may damage surfaces and remove hinge lubrication.
- Rinse with fresh water after exposure to salt water.
- Clean device with warm water and mild detergent as necessary; pat- or air-dry.
- Scrub dirt or grime using a soft bristle brush, remove dust with compressed air if necessary.
- Do NOT machine wash; Do NOT place in dryer.

Note: Strap components will absorb water and take longer to dry than the rest of the device or system.

Warranty

Össur Iceland ehf. includes a standard 2-year warranty (from ship date) on the mechanical function of the ThumbDriver. Please contact Össur Iceland ehf. for additional, extended warranty options.

Change in Performance

In the event of change in performance of a device, always cease using device and contact your clinic.

Warning: Do not use device if contact between device and patient's skin is causing pain, excessive abrasion, or bleeding.

REPORT OF SERIOUS INCIDENT

Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

LIABILITY

Össur Iceland ehf. does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.

- Device disassembled or amended with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

PRODUCT TECHNICAL SPECIFICATIONS

Materials

- Nylon
- Silicone
- Stainless Steel
- Titanium

Chemical Resistance

- Bleach
- Ammonia
- Acetone
- Rubbing Alcohol (Isopropanol)
- Gasoline
- Antifreeze
- Ethanol

Operational Temperature

- 0 – 120° F (-17 – 48° C)

Expected Service Life

- 3 years

Force Output

The ThumbDriver is designed such that any force the patient can generate at the MP and/or CMC joint can be transmitted by the devices. Force output is dependent on the patient's strength.

Silicone

Silicone will absorb gasoline or diesel if exposed and will experience a volume change. It will return to a normal shape if cleaned after exposure.

For more information: visit ossur.com

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