
ÖSSUR®
LIFE WITHOUT LIMITATIONS

Instructions for Use

C180™ ROCKET

3

EN | Instructions for Use 4

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NO | Bruksanvisning 17

DA | Brugsanvisning 19

SV | Bruksanvisning 22

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


1
F
E
B
A
C
D


2


3

4

5

6

7

8
0°
0°
20°- 30°
20°- 30°
30°- 40°
30°- 40°
10°- 20°
10°- 20°
5°- 10°
5°- 10°

ENGLISH

MD Medical Device

INTENDED USE

The device is intended for external support, stabilization and protection of the knee.
The device must be fitted and adjusted by a healthcare professional.

Indications for use

ACL, MCL, LCL, PCL, rotary and combined instabilities.
No known contraindications.

Warnings and Cautions:

- The use of the device may increase the risk of deep vein thrombosis and pulmonary embolism.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.
The patient should stop using the device and contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

This device is a supportive brace only and is not intended or guaranteed to prevent injury.

FITTING INSTRUCTIONS

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text **(Fig. 1)**.

Device Application

1. Unfasten all straps from the outer side of the device.
2. Bend the knee at a 90° angle, with the foot flat on the floor **(Fig. 2)**.
3. Place the device on the leg such that the patella is centered between the Condyle Pads (A) **(Fig. 3)**.

Ensure proper alignment of the device on the leg:

4. Height positioning: Align the center of the Hinge (B) slightly above the middle of the patella **(Fig. 3)**.
5. The device can be adjusted for individual knee width. To widen or narrow the width at the knee, the 6 mm (1/4") padding located under each cover can be removed and replaced with thinner pads [3 mm (1/8") or 1.5 mm (1/16")] or you can install a shim behind the 6 mm (1/4") pad to make the fit more snug. Any combination of pads can be used to achieve a more comfortable fit.

Fasten the straps in the following order. All straps can be adjusted to length.

6. Fit the Lower Cruciate Strap (C) first, located directly below the knee **(Fig. 4)**. This strap is located above the calf muscle and aids in the correct positioning of the device as it helps create the device-to-bone contact required to help keep the device in place.
7. Next, snugly fasten the Bottom Strap (D) by inserting it through the D-Ring.
8. Continue with the Upper Cruciate Strap (E), located just above the knee **(Fig. 5)**.
9. Fasten the Top Strap (F). Make sure that the strap hook tab does not overlap the elastic insert on the strap.

Note: Overtightening the top two straps can cause the device to migrate.

Device Adjustments

Range of Motion (ROM) Adjustment

The device comes with a 0° extension stop installed and includes the following additional stops:
5–10°, 10–20°, 20–30°, and 30–40° **(Fig. 8)**

1. Remove the Condyle Pad from the hinge.
2. Use a short size 2 Philips head screwdriver to remove the screw securing the extension stop **(Fig. 6)**.
3. Remove the extension stop out of the hinge.
4. Select the desired extension stop and place the stop into the hinge from the back with the small lip facing down **(Fig. 7)**.
5. Position the hole in the extension stop with the screw hole in the hinge plate. Re-insert the screw and fasten. Making sure to not overtighten.
6. Replace the condyle pad.
7. Repeat this procedure for the opposite hinge.
 Note: Ensure the same extension stop is used for both hinges.

Heat Molding Adjustment

For optimal fit and comfort, the shells can be heat molded.

1. Remove the straps and pads from the knee brace.
2. When adjusting a painted frame, remove the hook-and-loop from underneath the frame by heating it for a few seconds with a heat gun to soften the glue. Turn the buckles out of the way, if necessary.
3. Holding an industrial-style heat gun approximately 7.6-12.7 cm (3-5") from the surface of the knee brace, heat the desired area for 30-45 seconds. Always keep the gun moving to avoid concentrating the heat in one spot, as excessive heat may cause the material to delaminate.

Note: When heating a painted frame, heat only the underside of the frame as it will distort and, possibly, crack the paint.

Device Removal

1. Ask the patient to sit down with the leg bent 80° with the foot flat on the floor.
2. Unfasten all straps from the outer side of the device and attach the straps back to themselves. This will make it easier to put the device on next time and will also prolong the life of the straps.

Accessories and Replacement Parts

Please refer to the Össur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing the device with the soft goods detached allows for more thorough cleaning.

Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.

Note: If the device is used in salt water or chlorinated water, rinse with fresh water and air dry.

Hinge

- Remove foreign materials (e.g., dirt or grass) and clean using fresh water.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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©Copyright Össur 2022-09-30 IFU0324 1254_001 Rev. 6