ENGLISH

SYMBOLS

MD Medical Device

INTENDED USE

The device is intended for foot protection post injury and plantar offloading
The device must be fitted and adjusted by a healthcare professional.

Indications for use

- Post-operative use
- Plantar foot ulcers

No known contraindications.

Warnings and Cautions:

- Hot and/or cold therapy should not be used with this device.
- This device should not be used on patients with bony deformities.
- This device has been designed for the treatment of ulcers that are located on the plantar aspect of the foot.
- Hexes should be only removed and altered by a healthcare professional.
- The insole should not be used in conjunction with any other footwear device.
- If the relieved zone is too large and inhibits offloading and pressure relief by allowing the ulcer to bottom out, this device should not be used.
- In order to prevent individual hexes from collapsing into the relieved area, trim (scythe) the corners of the hexes lying on the perimeter of the relieved area. This method should only be performed when warranted by a healthcare professional and should only be performed by a healthcare professional.
- It is the healthcare professional's responsibility to monitor the patient's progress. Be sure the patient performs circulation checks.
- If a loss of circulation is felt, or if the device feels too tight, release contact closure straps and adjust to a comfortable level. If discomfort continues, the patient should contact a healthcare professional immediately.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device. Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.
The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Device Application

While carrying out the following instructions, please refer to the overview figure for locating components mentioned in the text (**Fig. 1**).

1. Remove the Hexagonal Insole (A) and Toe-Piece (B) from the shoe and then remove the Toe-Piece from the insole (**Fig. 2**).
2. Re-insert the Insole into the shoe. Make sure the Insole is centered evenly and there is no gap between the Insole and the heel counter (**Fig. 3**).

Note: The Insole must be positioned correctly for accurate marking.

3. Cover the ulcer with tape or gauze. Mark material with a transferable marker using iodine, lipstick, indelible ink pencil, etc. (**Fig. 4**).
4. Place foot inside the shoe, heel first (**Fig. 5**). Fasten Ankle Strap (C) and Contact Closure (D). To center the mark to the hexes, have the patient stand. The heel should be positioned comfortably against the heel counter.
5. Remove foot from the shoe and remove marked hexes. Scythe hex edges around the relief zone with a sharp cutting tool (**Fig. 6**).
Note: Scything hex edges prevents the hexes from "falling" into the relief zone.
6. Remove tape or gauze from the foot.

Open-toe option:

7. Place the foot inside the shoe, heel first. Position Ankle Strap Pad (E) such that it is centered over the foot (**Fig. 7**). Fasten the Ankle Strap and Contact Closure.
Ankle Strap Pad slides along the length of the Ankle Strap. The strap needs to be centered correctly to protect the foot from the Ankle Strap. Trim excess strap length.

Closed-toe option:

7. Lift the heel section of the Insole from the outside of the shoe. Insert the Toe-Piece onto the front of the Insole (**Fig. 8**). Make sure there is no gap between the front of the Insole and inside of the Toe-Piece.
Note: Do not remove the entire Insole from the shoe to keep the orientation of the offloading.
8. Insert the foot into the shoe, toe first. Fasten Ankle Straps and Contact Closure. Center the Ankle Strap Pad (E) for extra cushioning. The Insole should be positioned under the Ankle Strap to provide cushioning (**Fig. 9**).

Accessories and Replacement Parts

Please refer to the Ossur catalog for a list of available replacement parts or accessories.

USAGE

Cleaning and care

Washing Instructions

- Hand-wash using mild detergent and rinse thoroughly.
- Air dry, do not expose to direct heat.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.
Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Ossur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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© Copyright Ossur 2022
Ossur Africa (Pty) Ltd