
ENGLISH

MD Medical Device

DESCRIPTION
The device is an auto expulsion valve for use in a vacuum socket system.

Components (Fig. 1)
1. Bottom Attachment with Expulsion Valve
2. Fabrication Ring with O-Ring, Silicone Ring, and Silicone Seal (Included in L-641000)
3. Cross Driver (Included in L-680000)

INTENDED USE
The device is intended to connect and release a prosthetic system that replaces a missing lower limb.
The device must be fitted and adjusted by a healthcare professional.
Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.
It is intended to be used with a cushion liner and a sleeve.

Indications For Use and Target Patient Population
- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS
The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.
The device is for single patient use.

ASSEMBLY INSTRUCTIONS
Following the socket manufacture
1. Make sure that the seal between the Fabrication Ring and Attachment is uncontaminated, without debris.
2. Screw the Bottom Attachment with Expulsion Valve into the bottom of the socket.
3. Torque to 8 Nm using the Cross Driver and an 8mm torque wrench.

Warning: Failure to do this step properly can lead to loss of vacuum suspension.

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