
ENGLISH

MD Medical Device

DESCRIPTION
The Icelock 651 Smooth is a mechanical lock that offers a smooth transition and locks in every position. It is easily interchangeable and removable from the inside of the socket.

Components (Fig. 1)
1. Attachment Pin Smooth, 53 mm.
2. Screw Funnel (Included in L-651000).
3. Smooth Lock.
4. Release Unit.
5. Protective Sleeve.
6. Push Button.
7. Dust Cap.
8. Cross Driver (Included in L-680000).

INTENDED USE
The device is intended to connect and release a prosthetic system that replaces a missing lower limb.
Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.
The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population
- Lower limb amputation and/or congenital deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS
The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS
Note: Use with L-680000 Fabrication Tooling Kit (reusable)
1. Place the Smooth Lock into the socket, align the release unit hole with the hole in the socket.
2. Screw the Release Unit firmly and completely through the hole on the socket, into the Smooth Lock.
 Note: Ensure that sufficient material is removed around the release hole to allow the Release Unit to be fully screwed in position without restriction from the socket wall, else the function of the lock may be compromised.
3. Trim the Push Button to a suitable length and replace the Dust Cap.
4. If used with a prosthetic cosmesis, trim the protective sleeve to a suitable length.
5. Screw the Funnel on the threads above Smooth Lock and torque to 8 Nm using the Cross Driver.

USAGE
Cleaning and care
Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions
The device is Waterproof.
A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour. It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.
It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.
Dry with a cloth after contact with fresh water or humidity.
Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE
Unscrew the Release Unit, blow with compressed air alternately between the Pin hole and the hole for the Release Unit.
Caution: Do not use lubricant or any cleaning material. Only use air. Regular cleaning is necessary to prevent compromise on function of locking mechanism.

REPORT OF SERIOUS INCIDENT
Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL
The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY
Össur does not assume liability for the following:
- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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