
Instructions for Use



ÖSSUR®
LIFE WITHOUT LIMITATIONS


Icelock® 551
Expulsion Valve


			3
EN	\|	Instructions for Use	4
DE	\|	Gebrauchsanweisung	6
FR	\|	Notice d'utilisation	8
ES	\|	Instrucciones para el uso	11
IT	\|	Istruzioni per l'uso	13
NO	\|	Bruksanvisning	16
DA	\|	Brugsanvisning	18
SV	\|	Bruksanvisning	20
EL	\|	Οδηγίες Χρήσης	23
FI	\|	Käyttöohjeet	25
NL	\|	Gebruiksaanwijzing	28
PT	\|	Instruções de Utilização	30
PL	\|	Instrukcja użytkowania	33
TR	\|	Kullanım Talimatları	35
RU	\|	Инструкция по использованию	37
JA	\|	取扱説明書	40
ZH	\|	中文说明书	42
KO	\|	사용 설명서	44
CS	\|	Návod k použití	46
SK	\|	Indikácie použitia	49
RO	\|	Indicații de utilizare	51
HR	\|	Upute za uporabu	54
HU	\|	Használati javallat	56
BG	\|	Указания за употреба	58
SL	\|	Indikácie pre použitie	61
UA	\|	Інструкція з використання	63
AR	\|	تعليمات الاستخدام	65

1
4
3
2
1

2
7
6
5

3
min 4mm (3/16")

4

MD Medical Device

DESCRIPTION
The device is an auto-expulsion valve with push button for suction release.

Components: (Fig. 1 & 2)
1. Valve Housing
2. Rubber Membrane
3. Cap
4. Cap Key
5. Socket Insert
6. Set Screw
7. Foam Cap

INTENDED USE
The device is intended to connect and release a prosthetic system that replaces a missing lower limb.
Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.
The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population
- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for high impact use, e.g., walking and occasional running.

GENERAL SAFETY INSTRUCTIONS
The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.

The device is for single patient use.

ASSEMBLY INSTRUCTIONS
Directly to the socket
1. Fabricate the socket as usual.
2. Mark the position of the valve onto the socket.
3. Using an 8.5 mm drill bit, make a hole perpendicular to the socket wall.
4. Tap the hole using M10x1.5 tap.
5. Ensure the Cap and Rubber Membrane have been removed prior to shortening the threads to prevent debris lodging in the airways.
6. Measure socket wall thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.

Warning: Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off **(Fig. 3)**.

Using the Socket Insert (Fig. 2)
1. The Socket Insert is used as a fixation device for the valve, and should always be used when fabricating a flexible inner socket.
2. Prepare the plaster cast as usual.
3. Fix the Socket Insert to the cast using the Set Screw. Ensure it sits flush to the cast.

4. Plaster may have to be removed or added.
5. Apply the Foam Cap to the Set Screw.
6. Fabricate the socket as usual.
7. Grind material away over the socket insert until Foam Cap is fully exposed.
8. Unscrew the Set Screw to remove socket from cast.
9. Ensure the Cap and Rubber Membrane have been removed prior to shortening the threads to prevent debris lodging in the airways.
10. Measure socket wall and Socket Insert thickness at insertion point with a caliper, and shorten the threads accordingly using a hacksaw and smooth edges.

Warning: Do not shorten the threads beyond 4 mm, as it may affect the seal and increase the risk of the valve breaking off **(Fig. 3)**.

Valve Installation

1. Apply silicone adhesive to the threads and housing, where it comes in contact with the socket.
 Note: It is not recommended to use epoxy or similar as it can compromise the seal.
2. Install valve and hand tighten with a 5 mm hex key.
 Note: Be careful not to over-tighten as it may damage the threads.
 Note: After installation distribute silicone adhesive around the valve / socket contact point to ensure an air tight seal.
3. Ensure the valve parts and airways are free from debris before final assembly.
4. Replace Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity. **(Fig. 4)**
 Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.
5. Screw the Cap on using the Cap Key.

USAGE

Cleaning and care

Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions

The device is Waterproof.
A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.
It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.
It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.
Dry with a cloth after contact with fresh water or humidity.
Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE

The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

Disassembly

1. Unscrew the Cap using the Cap Key.
2. Remove the Rubber Membrane carefully.
3. Clean all parts with water and mild detergent.
4. Hold down the release button while cleaning airways.

Reassembly

1. Replace Rubber Membrane to Cap. Ensure the concave surface is facing the cap cavity **(Fig. 4)**.
2. Screw the Cap on using the Cap Key.

Caution: When reassembling ensure the Membrane is on the Cap and not in the Housing.

REPORT OF SERIOUS INCIDENT
Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL
The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY
Össur does not assume liability for the following:
- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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©Copyright Össur 2026-03-13