
Instructions for Use



ÖSSUR®
LIFE WITHOUT LIMITATIONS


Icelock® 850 Hybrid for Direct Socket



			3
EN	\|	Instructions for Use	6
DE	\|	Gebrauchsanweisung	10
FR	\|	Notice d'utilisation	14
ES	\|	Instrucciones para el uso	18
IT	\|	Istruzioni per l'uso	22
NO	\|	Bruksanvisning	26
DA	\|	Brugsanvisning	30
SV	\|	Bruksanvisning	34
EL	\|	Οδηγίες Χρήσης	38
FI	\|	Käyttöohjeet	42
NL	\|	Gebruiksaanwijzing	46
PT	\|	Instruções de Utilização	51
PL	\|	Instrukcja użytkowania	55
TR	\|	Kullanım Talimatları	59
RU	\|	Инструкция по использованию	63
JA	\|	取扱説明書	67
ZH	\|	中文说明书	71
KO	\|	사용 설명서	75
CS	\|	Návod k použití	78
SK	\|	Indikácie použitia	82
RO	\|	Indicații de utilizare	86
HR	\|	Upute za uporabu	91
HU	\|	Használati javallat	95
BG	\|	Указания за употреба	99
SL	\|	Indikácie pre použitie	103
UA	\|	Інструкція з використання	107
AR	\|	تعليمات الاستخدام	111


A
B
C
D
E
F
G
H
I
J
K
L
M


1


2


3
2-3x

4
2 cm


5

6


7

8
2
1


9


10


11


12



13



14



15
b
c
3Nm
a


16
4-6 cm
b
c
a


17


18

19

20

21

22

23
b
a

MD Medical Device

DESCRIPTION

The device offers a combination of mechanical and vacuum suspension. The mechanical suspension uses a lanyard, allowing the patient to pull themselves into the socket. The device can be used exclusively as mechanical suspension or in combination with an expulsion valve to provide either passive or active suspension.

Components

A. Attachment Pin for Icelock 850
B. Flat-head Screws for Direct Socket
C. Sealing Insert
D. Lock Housing for 4-hole
E. Quick Connect
F. Strap
G. Routing Insert
H. Latch
I. Slider Button with screw
J. Strap Lock with screws
K. Latch Housing Frame
L. Latch Housing
M. Quick Connect retrieval tool

INTENDED USE

The device is intended to connect and release a prosthetic system that replaces a missing lower limb.
Suitability of the device for the prosthesis and the patient must be evaluated by a healthcare professional.
The device must be fitted and adjusted by a healthcare professional.

Indications For Use and Target Patient Population

- Lower limb loss, amputation, or deficiency
- No known contraindications

The device is for low to high impact use, e.g., walking and occasional running.
The weight limit for the device is 166 kg.

GENERAL SAFETY INSTRUCTIONS

Warning: Use of a lower limb prosthetic device carries an inherent risk of falling which may lead to injury.
The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Warning: If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions, the patient should stop using the device and contact a healthcare professional.
Warning: Only use Össur replacement parts for the device for ensured durability.
The device is for single patient use.

REQUIRED COMPONENTS

The device should be used in combination with an Iceross Seal-In® X Locking liner, Iceross Seal-In® X Seal and Icelock® 552 Expulsion Valve or Unity® Valve.

ASSEMBLY INSTRUCTIONS

Socket Construction

1. Follow the instructions for use for the Direct Socket TF fabrication process until the fibers have been cut to length, before putting on the second silicone sheath.

If using the Icelock 552 Expulsion Valve, perform the following two steps.

2. Put a piece of double-sided adhesive tape onto the back of the red Dummy.
3. Stick the Dummy onto the inner silicone insulation sheath on the residual limb **(Fig. 1)**.

If using the Unity Valve, perform the following steps.

4. Put a piece of double-sided adhesive tape onto the back of the Socket Insert.
5. Place the Socket Insert in between the middle fiber layers.

The Latch lock can be placed either anteriorly or laterally on the socket, according to the patient's preference. Consider the patient's dominant hand and hand dexterity when choosing where to place the Latch lock. The selected side cannot be changed after fabrication.

6. On the selected side, mark the distal attachment at the half-way point between the threaded holes for guidance during fabrication. The Latch Housing will later be placed in line with the mark **(Fig. 2)**.
 Note: The injection port should always be placed anteriorly, regardless of which side the latch will be placed on.
 Note: As the central hole will be used for the mechanical fastening of the lock, the offset on the large distal attachment cannot be utilized.
7. Attach the fiber braids with distal attachment to the residual limb.
8. Turn the distal attachment nut onto the distal Attachment Pin and tighten nut firmly by hand.
 Note: Make sure no fiber braids are caught between the silicone cap and the Distal Connector.
9. Fold the top 2-3 layers of fiber off the residual limb **(Fig. 3)**.
10. Place the Latch Housing Frame so that it is at least 2 cm below the expected proximal trimline and in a direct line from the middle of the 4-hole connection (directly in line with the mark made in step 4), so the Strap will be routed straight **(Fig. 4)**. Orient the Frame so the round end points towards the proximal end. If it is placed incorrectly, the socket will be unusable.
 Note: It is recommended to place the Latch Housing as proximally as possible for easier donning.
11. Fold the layers of fiber over the Latch Housing Frame **(Fig. 5)**.
12. Push the Latch Housing into the Frame, with the fiber layers in between, until it snaps securely into place **(Fig. 6)**
 Note: Make sure that the arrow on the Latch Housing points towards the proximal end.
13. Prepare resin injection as described in instructions for use for the Direct Socket TF.
14. When the resin is distributed, maneuver resin by hand to make sure no air pockets form around the Latch Housing and resin fully saturates the fibers underneath the Latch Housing **(Fig. 7)**.
15. Finish the socket as described in the instructions for use for the Direct Socket TF.

If using the Icelock 552 Expulsion Valve.

16. Grind to expose the valve Dummy.

If using the Unity Valve.

17. Using an 8.5 mm drill bit, make a hole through the Socket Insert.
18. Tap the hole using M10x1.5 tap.

Lock Assembly

1. Remove all dust, dirt and other contamination from the socket and Lock Housing.
2. Remove plastic cap and apply pressurized air into the hole in the soft Dummy to loosen it from the socket and pull the soft Dummy out of the Latch Housing.
3. Apply medium strength threadlocker to the screw on the Latch lock and fasten the Latch lock by sliding it into the slot on the Latch Housing and hand-tightening the screw **(Fig. 8)**.
4. Route the Quick Connect into the Lock Housing and the Sealing Insert **(Fig. 9)**.
5. Slide the Routing Insert into the Lock Housing.
6. Assemble the Distal Connector and the lock to the socket, along with a 4-hole adapter of

choice. Apply medium strength threadlocker to the Flat-Head Screws and torque to 10 Nm **(Fig. 10)**.

7. Refer to the instructions for use for the Icelock 552 Expulsion Valve or Unity Valve for the assembly process.

Note: If the Strap end gets pulled out of the Routing Insert during assembly, it can be reassembled and threaded **(Fig. 11)**.

FITTING INSTRUCTIONS

Warning: Failure to adhere to the following instructions may cause the pin to detach from the liner resulting in loss of suspension, which may lead to injury.
Note: Setting the length of the Strap should be done before donning the system to ensure correct engagement.

1. Apply Loctite 410/411 or an equivalent compound onto the threads on the Attachment Pin and screw it into the liner. Torque the Attachment Pin to 4 Nm.
2. Pull the Quick Connect out of the socket and fasten to the Attachment Pin on the liner **(Fig. 12)**.
 Caution: Make sure to fully engage the pin into the Quick Connect.
3. Pull hard on the outermost Strap (a) so the liner is fully donned into the socket **(Fig. 13)**.
 Note: This ensures correct final position of the residual limb inside the socket once the lock has been fastened.
4. Adjust the length by using the inner two Straps. Pull the middle Strap (b) to elongate the system and the inner Strap (c) to shorten it **(Fig. 13)**. Adjust the Straps so they become tight when the Latch is nearly closed **(Fig. 14)**.
5. Fasten the two inner Straps (b) and (c) with the Strap Lock and hand-tighten the screws **(Fig. 15)**.
 Note: The screws have a pre-applied threadlocking patch.
6. Fit the device to the patient using the Donning and Doffing instructions described in the next section.
 Note: Only minor length adjustments can be made on the Strap after fitting the device. The Strap should not be lengthened by more than a few millimeters to ensure the Quick Connect will still enter the Lock Housing and seal with vacuum.
7. When the fitting is finished, cut the excess Strap (c) 4-6 cm above the Routing Insert and heat seal the edge to prevent fraying **(Fig. 16)**.

Note: If necessary, during the lifetime of the device, the strap should be retightened by a healthcare professional. Correct strap length ensures optimal function of device.
Caution: Make sure to fasten the two inner Straps with the Strap Lock, not only one, as it will not provide a secure enough fastening.

USAGE

Donning Instructions

1. Don the liner and the seal.
2. Secure the Quick Connect onto the liner pin **(Fig. 12)**.
 Caution: Make sure to fully engage the pin into the Quick Connect.
3. While sitting or standing, pull the latch while donning the socket to keep tension on the strap. This will help guide the Quick Connect into the distal end **(Fig. 17)**.
4. Stand up to push remaining air out of the socket before securing the latch **(Fig. 18)**.
 Caution: The latch won't be secured if the socket fit is too tight. Do not put excessive force into pulling the residual limb into the socket as this may cause too much tension on the distal end.

Doffing Instructions

1. Slide the latch lock down to release the latch **(Fig. 19)**. Remove the latch from the socket.
2. Push the valve to allow airflow inside the socket while pulling the residual limb out **(Fig. 20)**.
3. Push the black button on the Quick Connect and remove from the liner pin **(Fig. 21)**.

Store the device so the Quick Connect is pulled out of the socket and the latch is not locked **(Fig. 22)**.

Note: If the latch is locked while no liner is attached, the Quick Connect will be unreachable by hand at the distal end. To retrieve the Quick Connect, insert the liner **(Fig. 23a)** or the Quick Connect retrieval tool **(Fig. 23b)** into the socket so that the pin engages the Quick Connect, then pull it out of the distal end.

Cleaning and Care
Clean with a damp cloth and a mild soap. Dry with a cloth after cleaning.

Environmental Conditions
The device is Waterproof.
A Waterproof device can be used in a wet or humid environment and submerged in up to 3-meter-deep water for a maximum of 1 hour.
It can tolerate contact with: Salt water, chlorinated water, perspiration, urine, and mild soaps.
It can also tolerate occasional exposure to sand, dust, and dirt. Continuous exposure is not allowed.
Dry with a cloth after contact with fresh water or humidity.
Clean with fresh water after exposure to other liquids, chemicals, sand, dust, or dirt and dry with a cloth.

MAINTENANCE
The device and the overall prosthesis should be examined by a healthcare professional. Interval should be determined based on patient activity.

REPORT OF SERIOUS INCIDENT
Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.

DISPOSAL
The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY
Össur prosthetic devices are designed and verified to be safe and compatible in combination with each other and custom-made prosthetic sockets with Össur adapters, and when used in accordance with their intended use.

Össur does not assume liability for the following:
- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

Compliance
This device has been tested according to ISO 10328 standard to three million load cycles. Depending on patient activity, this may correspond to 3-5 years of use.

ISO 10328 - P7 - 166 kg *) ⚠

*) Body mass limit not to be exceeded!

⚠ For specific conditions and limitations of use see manufacturer's written instructions on intended use!

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©Copyright Össur 2025-10-31