
ENGLISH

SYMBOLS

MD Medical Device

MR Magnetic Resonance (MR) safe

INTENDED USE

The device is intended to provide gross immobilization to the cervical spine
The device must be fitted and adjusted by a healthcare professional.
X-ray and CT lucent.

Indications for use

Conditions requiring gross immobilization of the cervical spine. This may include:

- C-Spine precaution for trauma patients
- Immobilization for pre and post c-spine surgery
- Other conditions requiring gross immobilization of the mid-cervical spine

Contraindications

- Patients with a compromised airway or known spinal deformities such as ankylosing spondylitis.
- Patients with penetrating trauma injuries.

Warnings and Cautions:

Warning: If an unstable fracture is suspected or unknown, with or without a sustained trauma, ensure additional spinal precautions are implemented to immobilize the spine.
Warning: Use of a cervical collar may increase intracranial pressure (ICP) through jugular venous compression.
Warning: Use of a cervical collar may increase complexity of airway management.
Warning: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Impaired respiratory effort and forced expiratory volume
- Pneumonia
- Aspiration
- Worsening of existing cervical spine injury
- Severe neurological deterioration in patients with ankylosing spondylitis
- Triggering of non-compliance or agitation
- Delay to definitive treatment
- Impaired physical examination/secondary survey

Caution: Cervical spine immobilization, including use of a cervical collar, has been associated with:

- Increased pain and discomfort, which may lead to increased spinal movement
- Impaired swallowing

Caution: Perfumes and harsh cleaning agents should not be used on or under the collar, as they may compromise the integrity of materials.
Caution: Use of a cervical collar for cervical spine immobilization is not recommended in situations where patients are awake, alert, not intoxicated, without neck tenderness or pain and without abnormal sensory or motor findings on examination.
Caution: Regular cleaning of the collar, pads and the skin beneath along with inspection of the skin for any signs of irritation is required to reduce the risk of skin ulceration. Bedridden patients are at an increased risk for skin ulceration.

Precautions:

- At least two persons are needed for initial collar placement: one to maintain the patient's head and neck in proper alignment, the second to fit the collar.
- Collar removal or any adjustments to the collar should be done with physician permission only and according to the physician's instructions.
- Unless otherwise specified by the physician, the patient should not remove the collar except to wash under it and change the pads.
- The patient will need a second person to help when cleaning the collar.

These instructions do not supersede hospital protocol and/or direct orders of the patient's healthcare professional.

GENERAL SAFETY INSTRUCTIONS

The healthcare professional should inform the patient about everything in this document that is required for safe use of this device.
Any serious incident in relation to the device must be reported to the manufacturer and relevant authorities.
The patient should immediately contact a healthcare professional:

- If there is a change or loss in device functionality, or if the device shows signs of damage or wear hindering its normal functions.
- If any pain, skin irritation, or unusual reaction occurs with the use of the device.

The device is for single patient – multiple use.

FITTING INSTRUCTIONS

Size Selection

Proper sizing is critical to achieve comfort and proper immobilization.
Adjusting the size of the device is best achieved before the device is applied to the patient.
Ensure that correct spinal protocol is followed.
Height: Place the chin-cup of the front to the patients chin and slide the lower front down to the patient´s sternum **(Fig. 1)**. Lift the front away from the patient and push the two tabs to lock the size **(Fig. 2)**.
Alternatively, measure the vertical distance from the top of the patient's shoulder to the tip of the chin **(Fig. 3)** to identify the height setting of the collar. Remove the *"PEEL"* sticker and push the two tabs to lock the size **(Fig. 2)**.
Note: If the height falls between two sizes, select the smaller size.
Note: This device is designed to be sized and locked into a single position only once; it is not intended to be resized.
Circumference: Identify your patient´s profile type **(Fig. 3)** and refer to the size guide below. Alternatively, measure the patient's neck circumference. Adjust the back piece to the correct circumference by tearing away sections of the foam to change to a smaller circumference **(Fig. 4)**.

Patient Type	Neck Height	Circumference	Circumference Adjustment
Stout	2 (2.25" / 5.7 cm)	Large (16" - 19" / 41 - 48 cm)	No tear off
Average Male	3 (3.25" / 8.3 cm)	Medium (13" - 16" / 33 - 41 cm)	Tear off 1 section on each side
Thin Female	4 (4.25" / 10.8 cm)	Small (10" - 13" / 25 - 33 cm)	Tear off 2 sections on each side

Supine Application

Ensure that correct spinal protocol is followed.

- Slide the back piece behind the patient´s neck and center it, noting arrow on panel pointing up **(Fig. 5)**.

Note: Long hair should be placed outside of the back piece.

- Apply the front piece with the chin secured in the recess and center it to secure neutral alignment. The front piece overlaps the back piece to ensure effective immobilization and comfort **(Fig. 6)**. The front arrow should point upward.
- Apply the Hook and Loop strap on both sides. Tighten straps alternately to an equal length on both sides **(Fig. 7)**. Excess may be trimmed.
- The Plastazote material can be trimmed or heat moulded as necessary depending on the patient's anatomy. Ensure that all edges are rounded off smoothly.

Note: Device application must be secure to ensure proper fit, maintain alignment, and prevent the chin from slipping inside.

Device Removal

Ensure that correct spinal protocol is followed.
Undo the straps, remove the front piece and slide the back piece out from behind the head.

USAGE

Cleaning and care

- Wipe with a damp cloth and mild soap.
- Air dry.

Note: Do not machine-wash, tumble dry, iron, bleach, or wash with fabric softener.
Note: Avoid contact with salt water or chlorinated water. In case of contact, rinse with fresh water and air dry.

DISPOSAL

The device and packaging must be disposed of in accordance with respective local or national environmental regulations.

LIABILITY

Össur does not assume liability for the following:

- Device not maintained as instructed by the instructions for use.
- Device assembled with components from other manufacturers.
- Device used outside of recommended use condition, application, or environment.

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WWW.OSSUR.COM ©Copyright Össur 2025-01-02 IFU 1386_003 Rev. 7