




### System Height Example



		SPORT SOLUTIONS										JUNIOR SOLUTIONS		
		TT & TF										TT & TF		
		CHEETAH® XPLORE	FLEX-RUN™ WITH NIKE SOLE	CHEETAH® XCEED™	FLEX-FOOT CHEETAH®	CHEETAH® XTEND WITH NIKE SPIKE PAD	CHEETAH® XTREME WITH NIKE SPIKE PAD	CHEETAH® XCEL™ (SMALL)	CHEETAH® XCEL™ (LARGE)	CHEETAH® XPANSE™ (SMALL)	CHEETAH® XPANSE™ (LARGE)	CHEETAH® XPLORE JUNIOR	FLEX-RUN™ JUNIOR WITH NIKE SOLE	CHEETAH® JUNIOR WITH NIKE SOLE
ACTIVITIES	Long Jumping													
	Sprinting			•		• M	•	•	•	•	•			
	Distance running		•	•	•	•		•	•			•		•
	Recreational activities	•	•	•								•	•	•
SYSTEM HEIGHT		13 11/16" - 20 1/2" H 348 - 521mm H	10 1/4"† 261mm†	13 9/16"‡ 345mm‡	16 3/8" 411mm	13 - 17 3/4" 330 - 451mm	14 3/4 - 20 7/8" 375 - 530mm	17 7/16 - 21 3/16" 443 - 538mm	20 3/16 - 22 3/8" 513 - 569mm	14 15/16 - 18 3/8" 379 - 466mm	18 15/16 - 22 1/16" 481 - 576mm	7 5/16 - 12 1/4" 185 - 311mm	7 7/8" 200mm	7 3/16 - 12 1/16" 182 - 307mm
PATIENT WEIGHT LIMIT		325lbs (147kg)	287lbs (130kg)			325lbs (147kg)					256lbs (116kg)			
CONNECTORS*	Lamination connector	•			•	•	•	•	•	•	•	•		•
	L-bracket connector				•	•	•	•	•					
	Male bolt-on pyramid		•	•								•		
	Female bolt-on pyramid		•	•								•		
ALIGNMENT ADAPTERS*	Dynamic L-Adapter <sup>1</sup> 13°				•	•	•	•	•					
	Bolted L-Adapter <sup>2</sup> 13°				•	•	•	•	•					
	Dynamic Adapter <sup>1</sup> 10°				•	•	•	•	•					
	Bolted Adapter <sup>3</sup> 10°				•	•	•	•	•					
	Alignment Adapter Xpanse <sup>1,3</sup>									•		•		
SOLES*	Nike Traction SoleX			• St				•	•					
	Nike Spike SoleX			•				•	•					
	Nike Spike Pad 2.0			•	•	•	•	•	•	•	•			•
	Nike Spike Pad			•	•	• St	• St	•	•	•	•			•
	Nike Sole		• St											
	Nike Sole Junior											• St		• St

\*Not included unless stated otherwise, †Measured without sole and male/female pyramid, ‡Measured with Traction SoleX and without male/female pyramid, St Standard (comes equipped with), M Mid-distance sprinting, H Height of posterior brim, 1 Alignment only, 2 Long term, 3 Short term.