

INNOVATOR X®

Benefiting from our proven, patented technologies, Innovator X, with extending struts, can help improve patient outcomes. With superior comfort, clean aesthetics and the most user-friendly hinge in post-op bracing, Innovator X is designed to increase patient satisfaction and compliance.



INDICATIONS FOR USE

- Controlled elbow flexion/extension following surgery or trauma

FEATURES AND BENEFITS

- Easy-to-use Innovator hinge
- Slide-to-Size™ struts lead to a more customized fit
- Sensil® silicone padding on the strapping system reduces slipping and malleable aluminium cuffs allow for customization
- Patented Formfit® 3-dimensional padding and lycra lining provide superior comfort and compliance
- Comfort pad at the neck and shoulder
- Flexion range can be quickly adjusted to 0, 10, 20, 30, 45, 60, 75, 90, 105 and 120. Extension range can be adjusted to 0, 10, 20, 30, 45, 60, 75 and 90.
- Optional anti-supination/pronation arm bar kit is available

REQUIRED MEASUREMENTS

- Universal, one size

ORDERING INFORMATION

INNOVATOR X

Side	Part#
Left	W-823101
Right	W-823100

ACCESSORIES

Description	Part#
Innovator X Bar Kit	W-A82310
Replacement Softgoods Left	W-823111
Replacement Softgoods Right	W-823110