

MIAMI LSO™

Miami LSO is a spinal orthosis designed for pain relief, protection of injured ligaments or muscles and post-surgical immobilization. It provides low-profile, anatomical immobilization, easy fit and wear, options for universal or modular support and an advanced pulley system with a length-adjustable pulley cord.



INDICATIONS FOR USE

Conditions requiring gross immobilization of the trunk in the lumbar region. Such conditions may include:

- Stable, non-displaced spinal fractures (L1-L5)
- Spinal stenosis
- Herniated discs
- Degenerative spinal pathologies
- Spondylolisthesis
- Spondylolysis

Contra-indication: Unstable, displaced spinal fractures

FEATURES AND BENEFITS

- Slim, breathable material
- Proven Miami double-pulley system with length-adjustable pulley cord
- Intuitive fitting without trimming with fold-over belt arms
- Optional fit tool for perfect fitting directly on the patient
- Universal belt sizing accommodates waists up to 178 cm/70"
- New Universal Posterior panel to accommodate more patients with less inventory
- Contoured posterior panels sized S-XL for a precise, anatomical fit
- Quick-Fit pocket enclosure and blue highlights for patient ease of use
- Immobilization from L1 – L5 with coverage from S1 – T9


REQUIRED MEASUREMENTS



- Waist circumference






MIAMI LSO™ ORDERING INFORMATION

OPTION 1: Order braces for inventory

MIAMI LSO BELT		
		
SIZE	WAIST CIRC.	PART #
6" Universal (Slim)	26-50" (66-127 cm)	MLSO-1
8" Universal (Standard)	26-50" (66-127 cm)	MLSO-2
6" (Slim) Belt Extension	26-70" (66-177 cm)	MLSO-1EXT
8" (Standard) Belt Extension	26-70" (66-177 cm)	MLSO-2EXT

CONTOURED POSTERIOR PANEL			CONTOURED POSTERIOR PANEL WITH LATERAL SUPPORTS		
					
SIZE	WAIST CIRC.	PART #	SIZE	WAIST CIRC.	PART #
Small (S)	26-34" (66-86 cm)	MLSO-P11	Small (S)	26-34" (66-86 cm)	MLSO-P12
Medium (M)	34-38" (86-96 cm)	MLSO-P21	Medium (M)	34-38" (86-96 cm)	MLSO-P22
Large (L)	38-50" (96-127 cm)	MLSO-P31	Large (L)	38-50" (96-127 cm)	MLSO-P32
Extra Large (XL)	50-66" (127-167 cm)	MLSO-P41	Extra Large (XL)	50-66" (127-167 cm)	MLSO-P42

STANDARD POSTERIOR PANEL		UNIVERSAL POSTERIOR PANEL		ANTERIOR PANEL	
					
SIZE	PART #	TYPE	PART #	SIZE	PART #
Standard	MLSO-P56	with Lateral	MLSO-P55	6" Panel (Slim)	MLSO-PA6
		without Lateral	MLSO-P57	8" Panel (Standard)	MLSO-PA8

OPTION 2: Order a brace in a kit.

KIT	POSTERIOR PANEL	WITHOUT LATERAL SUPPORTS	WITH LATERAL SUPPORTS
	WITH ANTERIOR PANEL		
Slim Belt	S	MLSO-111-1	MLSO-112-1
	M	MLSO-121-1	MLSO-122-1
	L	MLSO-131-1	MLSO-132-1
	XL	MLSO-141-1	MLSO-142-1
	UNIVERSAL	MLSO-157-1	MLSO-155-1
		WITHOUT ANTERIOR PANEL	WITH ANTERIOR PANEL
	STANDARD	MLSO-156-2	MLSO-156-1

KIT	POSTERIOR PANEL	WITHOUT LATERAL SUPPORTS	WITH LATERAL SUPPORTS
	WITH ANTERIOR PANEL		
Standard Belt	S	MLSO-211-1	MLSO-212-1
	M	MLSO-221-1	MLSO-222-1
	L	MLSO-231-1	MLSO-232-1
	XL	MLSO-241-1	MLSO-242-1
	UNIVERSAL	MLSO-257-1	MLSO-255-1
		WITHOUT ANTERIOR PANEL	WITH ANTERIOR PANEL
	STANDARD	MLSO-256-2	MLSO-256-1