

ÖSSUR LEGS TRANSTIBIAL PROSTHESIS ORDER FORM

FACILITY INFORMATION

Company Name:

Clinician Name:

Contact Number:

Email:

Purchase Order#

Ship to Address:

REQUIRED PATIENT INFORMATION

Patient Name: Impact Level: ☐ High ☐ Med ☐ Low
 Weight: kg. K-Level: ☐ K1 ☐ K2 ☐ K3 ☐ K4 ☐ Left ☐ Right







1. TT LINERS (check one then fill in desired sizes and options)

Limit 2 liners of same type, 2 Seal-Rings, 2 Distal Cups, 1 Sleeve - per order.

LINER OPTIONS

Locking Liners







☐ Iceform Locking ☐ Original Locking ☐ Comfort Locking ☐ Dermo Locking ☐ Synergy Locking ☐ Sport Locking

☐ Wave ☐ Conical ☐ Wave ☐ Wave

Cushion Liners

☐ Iceform Cushion ☐ Comfort Cushion ☐ Dermo Cushion ☐ Synergy Cushion ☐ Relax 3C Cushion ☐ Activa Cushion

☐ Conical ☐ Wave ☐ Wave

Seal-In Liners

☐ Seal-In V with Wave ☐ Seal-In V HP with Wave ☐ Dermo Seal-In ☐ Seal-In X5 with Wave ☐ Seal-In X with Rings ☐ Seal-In X Locking

1. TT LINERS (check one then fill in desired sizes and options)

LINER SPECIFICATIONS

Liner Size	Optional: Add additional liner of same type
<div> <div>LINER SIZE</div> <div> <input type="text"/> </div> </div> <input type="checkbox"/> 3mm <input type="checkbox"/> 6mm <input type="checkbox"/> 9mm (Iceform Only)	<div> <div>LINER SIZE</div> <div> <input type="text"/> </div> </div> <input type="checkbox"/> 3mm <input type="checkbox"/> 6mm <input type="checkbox"/> 9mm (Iceform Only)










OPTIONAL ADD-ONS

Optional: Add Distal Cups (Limit 2)	Optional: Add Sleeve (Limit 1)	Optional: Add Seal-Rings (Limit 2)
<div> <div>DISTAL CUP SIZE</div> <div> <input type="text"/> </div> </div> <input type="checkbox"/> Distal Cup 1	<div> <div>SLEEVE SIZE</div> <div> <input type="text"/> </div> </div> <input type="checkbox"/> Iceform Sleeve <input type="checkbox"/> Iceross Sleeve <input type="checkbox"/> Genu Sleeve	<div> <div>Seal-Ring 1</div> <div> <input type="checkbox"/> Classic <input type="checkbox"/> Grip <input type="checkbox"/> Volume </div> </div> <div> <div>Seal-Ring 2</div> <div> <input type="checkbox"/> Classic <input type="checkbox"/> Grip <input type="checkbox"/> Volume </div> </div>
<div> <div>DISTAL CUP SIZE</div> <div> <input type="text"/> </div> </div> <input type="checkbox"/> Distal Cup 2	<div> <div>RING SIZE</div> <div> <input type="text"/> </div> </div>	<div> <div>RING SIZE</div> <div> <input type="text"/> </div> </div>

NOTES:






2. SUSPENSION OPTIONS (Choose for conventional socket fabrication) Refer to section 3 for Direct Socket Suspension Options Limit 1 lock/valve kit per order. All socket adapters and modular components are discounted at 10%.

LOCKING OPTIONS

<input type="checkbox"/> 680 Fabrication Adaption Kit 	<input type="checkbox"/> 621 Ratchet Pin 	<input type="checkbox"/> 651 Smooth Pin 	<input type="checkbox"/> 562 Hybrid 	<input type="checkbox"/> 214 Clutch 	<input type="checkbox"/> 200 Series Lanyard 
Works With					
<input type="checkbox"/> 600 Series Pyramid Socket Adapter 	<input type="checkbox"/> 600 Series 4-Hole Socket Adapter 				<input type="checkbox"/> 211 Clutch 

☐ Add Extra Pin: ☐ Long ☐ Standard ☐ Short

EXPULSION OPTIONS

EXPULSION OPTIONS				UNITY
<input type="checkbox"/> 544 Expulsion Plate Kit with Unity 	<input type="checkbox"/> 544 Expulsion Plate Kit for Seal-In 	<input type="checkbox"/> 544 Expulsion Valve Kit for Sleeves 	<input type="checkbox"/> 551 TT Expulsion Valve 	<input type="checkbox"/> Unity Valve Kit 
				Works With
			<input type="checkbox"/> 544 Socket Adapter 	<input type="checkbox"/> Male Pyramid Socket Adapter 

3. DIRECT SOCKET SUSPENSION OPTIONS (Choose when fabricating your socket using a Direct Socket Kit)

Limit 1 lock/valve kit per order

LOCKING

☐ 621 Ratchet Pin



☐ 651 Smooth Pin



☐ 214 Clutch Lock



☐ 562 Hybrid Lock with Unity (for TT only)



☐ 562 Hybrid Lock with Passive Vacuum (for TT only)



Requires: Modular Spacer Plate



☐ Add Extra Pin: ☐ Long ☐ Standard ☐ Short

EXPULSION

☐ Expulsion Valve



Requires: Bottom Seal



☐ Expulsion Plate



Requires: DS Adapter Plate



UNITY

☐ Unity TT Valve



Requires: Bottom Seal



☐ Unity Plate



Requires: DS Adapter Plate



4. OPTIONAL: DIRECT SOCKET TT KITS (Select material and size)

Limit 1 Direct Socket Kit per order

GLASS

Resin: ☐ Black ☐ Neutral

Braid: ☐ 4 inch ☐ 5 inch ☐ 7 inch



CARBON

☐ 5 inch with black resin

☐ 7 inch with black resin



BASALT

☐ 4 inch with black resin

☐ 5 inch with black resin

☐ 7 inch with black resin


























































5. ENDOSKELETAL COMPONENTS (check all that apply)

Limits apply on some components. Require a minimum of 2 units; Max 5 total components. Refer to below items;

Socket components: Limit 2; Compatible types
Pylons: Limit 2

Tube Clamps: Quantity limitations vary depending on choice of pylons or pylon kits
Torsion Shock Adapter: Limit 1 (when foot does not include torsion/shock adapter)

4-HOLE ADAPTERS	DOUBLE ADAPTERS	SINGLE ADAPTERS
 <input type="checkbox"/> 4-Hole Male 	 <input type="checkbox"/> Female Double Adapter 32mm 	 <input type="checkbox"/> Male Single Adapter Short 
 <input type="checkbox"/> 4-Hole Female 	 <input type="checkbox"/> Female Double Adapter 45mm 	 <input type="checkbox"/> Male Single Adapter Long 
 <input type="checkbox"/> 4-Hole Male w/ Rotation 	 <input type="checkbox"/> Female Double Adapter 60mm 	 <input type="checkbox"/> Female Single Adapter Short 
 <input type="checkbox"/> 4-Hole Male w/ 8mm Ax Offset 	 <input type="checkbox"/> Female Double Adapter 75mm 	 <input type="checkbox"/> Female Single Adapter Long 
 <input type="checkbox"/> 4-Hole Female w/ Rotation 	TORSION SHOCK ADAPTER	
 <input type="checkbox"/> 4-Hole Female Pyramid 20/20 Offset 10° Right 	 <input type="checkbox"/> Cat 1  <input type="checkbox"/> Cat 2  <input type="checkbox"/> Cat 3 	
 <input type="checkbox"/> 4-Hole Female Pyramid 20/20 Offset 10° Left 	TUBE CLAMPS	
PRONG ADAPTERS		
 <input type="checkbox"/> 4-Prong Socket Adapter 	PYLONS	
 <input type="checkbox"/> 3-Prong Socket Adapter 	 <input type="checkbox"/> Female Short   <input type="checkbox"/> Female Kit Short   <input type="checkbox"/> Female Long   <input type="checkbox"/> Female Kit Long 	
 <input type="checkbox"/> Male Pyramid Insert for Prong 	HEIGHT ADJUSTABLE PYLONS	
 <input type="checkbox"/> Female Pyramid Insert for Prong 	 <input type="checkbox"/> Standard  <input type="checkbox"/> Long 	
 <input type="checkbox"/> Female Pyramid Insert for Prong LP 	<input type="checkbox"/> Please tick this box if you would prefer Titanium adaptors. * Patients under 100kg will be supplied with Aluminium adaptors unless specified. * Patients between 100-166kg will be provided with Titanium adaptors. * Titanium pylon kits are waterproof.	
SPECIFY OTHER COMPONENTS AND/OR QUANTITIES DESIRED		

NOTES:

6. FOOT: WITH OR WITHOUT UNITY (check one and fill out the options below)

Limit 1 per order. Trial Proprio Foot units are available for use prior to purchase. Bionic units may not be returned or exchanged following purchase. Mechanical feet may be exchanged during the 60-day trial period.

FOOT OPTIONS: Foot Size: cm.

Category (1-9)

☐ Left ☐ Right

Foot Cover: ☐ Beige ☐ Brown

LOW ACTIVITY



- ☐ Flex-Foot® Balance
☐ with Unity®



- ☐ Flex-Foot® Balance
with D/P Flexion
☐ with Unity®



- ☐ Balance™ Foot J
☐ with Unity®



- ☐ Flex-Foot® Assure
☐ with Unity®



- ☐ Balance™ Foot S
☐ with Unity®



- ☐ Balance™ Foot S
Torsion
☐ with Unity®



BIONIC SOLUTION

LOW TO MODERATE ACTIVITY



- ☐ PROPRIO FOOT®
☐ with Unity®



MODERATE ACTIVITY



- ☐ LP Vari-Flex®
☐ with Unity®



- ☐ Vari-Flex®
☐ with Unity®



- ☐ Vari-Flex® Modular



- ☐ Talux®



- ☐ Pro-Flex® Pivot
☐ with Unity®



- ☐ Pro-Flex® LP
☐ with Unity®



- ☐ Pro-Flex® LP Align



HIGH ACTIVITY



- ☐ Cheetah® Xplore
☐ Foot Cover Kit
☐ Lamination
Connector Kit



- ☐ Re-flex Shock™
☐ with Unity®
☐ Standard
☐ Tall



- ☐ Re-flex Rotate™
☐ with Unity®
☐ Standard Male
☐ Tall Male
☐ Standard Female
☐ Tall Female



- ☐ Pro-Flex® XC Torsion
☐ with Unity®



- ☐ Pro-Flex® XC
☐ with Unity®



- ☐ Pro-Flex® LP Torsion
☐ with Unity®



☐ Other (Please Specify)

FOOT AND ADAPTER CHARACTERISTICS



WATERPROOF

Provides permanent protection against harmful ingress of water. Even after submersion in water, the foot's structural parts maintain their integrity and continue to provide full support, function and durability.



WEATHERPROOF

Allows for use in a wet and/or humid environment but does not allow submersion. Fresh water splashing against the enclosure from any direction shall have no harmful effect. Dry thoroughly after contact with fresh water. Fresh water: Includes tap water. Excludes salt and chlorinated water.

NOTES: