



Nothing measures like *SmartMeasure*™

An app for iPhone® & iPad® designed to order pre-fabricated knee braces or measure custom knee braces. *SmartMeasure* is Össur's most intuitive, efficient and accurate way to measure and order Össur knee braces, including Rebound® DUAL, Rebound® PCL and Rebound® Cartilage.



INTUITIVE

SmartMeasure utilises a section-based, step-by-step approach to guide you through the entire brace selection, measuring and order process.

With built-in tutorials and completion indicators for each section, every aspect was designed with the user experience in mind.



EFFICIENT

SmartMeasure makes every step of the measuring and ordering process quicker. Favorites and smart fields make the data entry process faster and, once completed, the order is sent directly from your device.



ACCURATE

SmartMeasure uses the iPhone and iPad's gyroscopic sensor to reduce lens angle distortion on the anterior and lateral photos. This, in conjunction with a series of specific anatomical landmarks and Össur's proprietary software, provides the precise data needed to build an intimate fitting custom brace.

