

POWER KNEE™

A New Era in Motion

The World's First Actively Powered
Prosthetic Knee.

- Versatile, easy, and intuitive to use.
- Auto-adaptive assistance during stance and swing phase helps improve gait.
- Provides active flexion and extension during walking.
- Supportive stance flexion in level ground walking.
- Provides motor-powered active assistance during stair and ramp ascent and sit-to-stand.
- Controlled resistance when sitting down, descending stairs and walking down ramps.



COMPATIBLE WITH A WIDE VARIETY OF ÖSSUR FEET
AVAILABLE IN ÖSSUR LEGS



POWER KNEE™

PRODUCT OVERVIEW

The Power Knee is unlike anything you or your patients have experienced. A new category of Actively Powered prostheses, you can expect to invest less time in adjustments and more time for helping your patients live their lives fully. Power Knee will improve your patients' short and long term outcomes with the proven benefits of Power – increased symmetry, greater endurance, and new capabilities with powered extension and stance phase control.

Building on proven outcomes, the benefits begin the moment of fitting and extend into whatever activities your patients take on. Power Knee is the most advanced prosthetic knee technology available.

IDEAL FOR

- Transfemoral (bi- or unilateral), knee disarticulation, and hip-disarticulation amputees with a desire to have better symmetry, more natural gait, and participation in the daily activities they enjoy.
- Transfemoral amputees seeking increased assistance from their prosthetic device and who desire the ability to do more.
- Community and workplace ambulators with the ability or potential for variable cadence and ramp/stair walking.
- Low to moderate impact levels.

TECHNICAL SPECIFICATIONS

- Max patient weight: 116 kg
- Product weight (incl. battery pack): 2.65 kg
- Build Height: 27.6 cm
- Knee Flexion: 120°
- Battery Life: 25 hours (during normal use)
- Recharge time: 5h30min for full charge

ORDERING INFORMATION

- PKA10001 Power Knee
Includes 2 replaceable battery packs and charger

WARRANTY INFORMATION

- 2-year limited warranty on knee and charger

EXTENDED WARRANTY OPTIONS

- Additional warranty may be purchased while the device is still under warranty to a total of 6-years of protection. 1, 2, 3 and 4-year warranty extensions are available.

POWER KNEE is currently under investigation for compatibility with osseointegration. It is not presently recommended for use with Osseointegration. Updates will be provided as soon as they become available.



FOLLOW ÖSSUR ON



WWW.OSSUR.COM.AU

Össur Australia

TEL 1300 123 268

FAX +61 2 9475 1114

customercare.au@ossur.com

