

TT SOCKETLESS ÖSSUR LEG ORDER FORM

Send all casts, sockets, files to: 7199 S Conway Rd #100, Orlando, FL 32812
 tel: (888)839-6213, fax: (800)788-9878 | cs@ossur.com | www.ossur.com

BILLING & SHIPPING INFORMATION

Össur Account #:

PO#

Ship To:

Company:

Contact:

Address:

Phone:

Fax:

City/State/Zip:

Email:

REQUIRED PATIENT INFORMATION

Last Name: Impact Level: High Med Low Amputation Level: TT Symes
 First Name: Weight: lbs. K-Level: K1 K2 K3 K4 Left Right Bilat Order

1: LINERS (check one then fill in desired sizes and options)

Comfort Locking 	Dermo Locking 	Iceform Locking 	LINER SIZE Standard Wave 3 mm 6 mm	Seal-In V HP 	Seal-In X5 	LINER SIZE Standard Wave RING SIZE
Comfort w/ Sleeve 	Dermo w/ Sleeve 	Iceform w/ Sleeve 	LINER SIZE Standard Wave 3 mm 6 mm SLEEVE SIZE	Seal-In V 	Seal-In X with Rings 	Seal-In X Locking with Rings (Ring Options)  Classic Volume Grip

2. SUSPENSION & SOCKET ADAPTER (check all that apply)

LOCKING OPTIONS					
Icelock 621 Ratchet & Extra Pin 	Icelock 651 Smooth & Extra Pin 	Icelock 562 Hybrid 	Icelock 600XM 	Icelock 211 Clutch & Extra Pin Icelock 214 4-Hole Clutch & Extra Pin 	Icelock 234 Lanyard 
Works With					
Icelock 672 Male Pyramid Ti Icelock 673 Male Pyramid Al 	Icelock 674 4-Hole Adapter Al Icelock 675 4-Hole Adapter Ti 				
EXPULSION OPTIONS			UNITY ELEVATED VACUUM OPTIONS		
Icelock 544 Expulsion Kit 	Icelock 551 Expulsion Valve 	Unity Valve TT 	Icelock 544 Unity Plate 		
Icelock 544 Expulsion Valve Kit for sleeves 	Works With				
	Icelock 544 Socket Adapter 		3-Prong Socket Adapter 4-Prong Socket Adapter		

Other (write in):

Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

TT SOCKETLESS ÖSSUR LEG ORDER FORM

3: ENDOSKELETAL COMPONENTS (check all that apply)



4-HOLE ADAPTERS

4-Hole Male



4-Hole Female

4-Hole Male w/ Rotation

4-Hole Female w/ Rotation

4-Hole Male w/ 8mm Ax

Offset



4-Hole Female Pyramid 20/20
Offset



4-Hole Female Pyramid 20/20
Offset, 10 degree, Right



4-Hole Female Pyramid 20/20
Offset, 10 degree, Left



INSERTS FOR PRONG

Male Insert for Prong



Female Insert for Prong

LP Female Inserts for Prong



DOUBLE ADAPTERS

Female Double Adapter 32mm

Female Double Adapter 45mm

Female Double Adapter 60mm

Female Double Adapter 75mm

Offset

PYLONS

Height Adjustable Standard

Height Adjustable Long

Female Short Female Kit Short

Female Long Female Kit Long



TUBE CLAMPS

Male Tube Clamp

Female Tube Clamp

4-Hole Tube Clamp



Male Single Adapter Short

Male Single Adapter Long



Female Single Adapter Short

Female Single Adapter Long

Other:






















Össur may select different components to fit your dimensions if necessary. If you do not want Össur to make any changes, check this box:


TT SOCKETLESS ÖSSUR LEG ORDER FORM

PRODUCTS

7: FOOT: WITH OR WITHOUT UNITY (check one and fill out the options below)

No Foot **FOOT OPTIONS:** Foot Size: _____ cm. Category (1-9) _____ Left _____ Right Foot Shell: _____ Beige _____ Brown _____

		Suggested L-Codes†			Suggested L-Codes†															
	Balance Foot S  with Torsion  with Unity	5972, 5986 5972, 5984, 5988 5781*		K2 Sensation with DP Flexion with Unity	5972, 5986 5972, 5968 5781*															
	Balance J with Unity	5976 5781*		Assure with Unity	5981 5781*															
	Pro-Flex ST 	5981, 5986		Vari-Flex with Unity	5981‡ 5781*															
	Talux**	5981‡, 5986‡		Pro-Flex LP Pro-Flex LP Torsion with Unity	5981, 5986 5981, 5984, 5988 5781*															
	Pro-Flex Pivot with Unity	See Össur reimbursement guide		Proprio Foot with Unity	5973‡ 5781*															
	Pro-Flex LP Align	5981, 5986, 5990		ReFlex Rotate with Unity	5987‡, 5984‡ 5781*															
	Pro-Flex XC  Pro-Flex XC Torsion with Unity	5981‡ 5987‡, 5984‡ 5781*		Torsion Shock Adapter with Male Pyramid Adapter with Female Pyramid Adapter	5984, 5988															
	ReFlex Shock with Unity	5987‡ 5781*																		
	Pro-Flex Modular 	5980‡, 5986‡	Desired build height (mm): With foot cover on. No return or exchange after foot is cut.																	
				<table border="1"> <thead> <tr> <th>Foot Size</th> <th colspan="2">Build Height</th> </tr> <tr> <td></td> <th>Min</th> <th>Max</th> </tr> </thead> <tbody> <tr> <td>22-24cm</td> <td>157mm</td> <td>313mm</td> </tr> <tr> <td>25-27cm</td> <td>168mm</td> <td>327mm</td> </tr> <tr> <td>28-30cm</td> <td>176mm</td> <td>340mm</td> </tr> </tbody> </table>	Foot Size	Build Height			Min	Max	22-24cm	157mm	313mm	25-27cm	168mm	327mm	28-30cm	176mm	340mm	
Foot Size	Build Height																			
	Min	Max																		
22-24cm	157mm	313mm																		
25-27cm	168mm	327mm																		
28-30cm	176mm	340mm																		

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

**Discontinued 31 Dec 2022

† The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡ PDAC approved.

KNIT-RITE SOFT-SOCKS

Knit-Rite socks sold separately, not included in Össur Legs pricing.

MEASUREMENT INSTRUCTIONS

Socks are sized in flat measurement, half of the circumference measurement correspond to the size measurements listed in the table. If the sock is used with a liner please measure the circumference with the liner on the patient.



Össur Soft-Sock — Black Polyester/Lycra®							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Long	22"–28" (55–70 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	X-Wide Short	10"–14" (25–35 cm)	10" (25 cm)	7" (17.5 cm)	34, 36, 38, 40
Qty.	Qty.	Qty.	X-Wide Medium	16"–20" (40–50 cm)	10.0–12.0" (25–30 cm)	7–8" (17.5–20 cm)	36, 38, 40, 42, 45

* measured 2 inch (5cm) above distal end

Össur Hole-In-Tpe Soft-Sock - Black Polyester/Lycra®							
Lightweight	3-Ply	5-Ply	Size	Length	Top width	Toe width*	Fits to liner size
Qty.	Qty.	Qty.	Narrow Short	10"–14" (25–35 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Narrow Medium	16"–20" (40–50 cm)	6" (15 cm)	4" (10 cm)	18, 20, 22
Qty.	Qty.	Qty.	Regular X-Short	8"–12" (20–30 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Short	10"–14" (25–35 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Medium	16"–20" (40–50 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Regular Long	22"–28" (55–70 cm)	7–8" (17.5–20 cm)	4–5" (10–12.5 cm)	23.5, 25, 26.5, 28
Qty.	Qty.	Qty.	Wide Short	10"–14" (25–35 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34
Qty.	Qty.	Qty.	Wide Medium	16"–20" (40–50 cm)	9" (23 cm)	6" (15 cm)	28, 30, 32, 34

* Measured 2 inch (5cm) above distal end

CLEAR FORM

6 of 6

FOLLOW ÖSSUR ON



ÖSSUR CUSTOM SOLUTIONS
 7199 S. CONWAY RD #100 | ORLANDO, FL 32812
 TEL (888) 839-6213, FAX (407) 367-3695
 CS@OSSUR.COM | OSSUR.COM

