

AFO DYNAMIC

CARBON AFO



AFO DYNAMIC

AFO Dynamic incorporates Ossur Flex-Foot® technology and the energy storing properties of carbon fibre. The result is a strong, lightweight solution for people with varying degrees of drop foot.

INDICATIONS

Drop foot secondary to CVA, MS, Charcot Marie Tooth disease, other neurological conditions, or mild knee instability

Contra-indicated for severe ankle-foot deformities, severe spasticity, or fluctuation edema

Suggested weight range: up to 120kg

FEATURES & BENEFITS

Carbon fibre design provides excellent strength-to-weight ratio and dynamic motion

The tapered carbon heel stores energy and manages forces at heel strike

Full-length toe lever provides full support while foot is loaded, as well as reduced shock to the contralateral limb

Controlled plantar flexion at initial contact

Creates M-L stability at mid-stance

Reduces energy consumption by reducing hip hike

Medial strut makes it less visible in shoe

Open anterior tibia shell provides superior anatomical fit, and reduces heat and perspiration

ORDERING INFORMATION

PART #	SIDE	SIZE	MENS SHOE	WOMENS SHOE	FOOT PLATE LENGTH
O-102011	Right	X-Small	3 - 5	4 - 6.5	22cm
O-102021	Left	X-Small	3 - 5	4 - 6.5	22cm
O-102012	Right	Small	5.5 - 7.5	7 - 9	24cm
O-102022	Left	Small	5.5 - 7.5	7 - 9	24cm
O-102014	Right	Medium	8 - 10.5	9.5 - 12	26cm
O-102024	Left	Medium	8 - 10.5	9.5 - 12	26cm
O-102016	Right	Large	11 - 14	12+	28cm
O-102026	Left	Large	11 - 14	12+	28cm

ACCESSORIES

PART #	REPLACEMENT SOFTGOODS	SIZE
O-102102	AFO Dynamic Accessory Kit	Small
O-102104	AFO Dynamic Accessory Kit	Medium
O-102106	AFO Dynamic Accessory Kit	Large