

# AFO Light

AFO Light provides support only where it is needed, eliminating unnecessary bulk and providing a cosmetically-appealing look. Using the best of Flex-Foot® technology and the energy storing properties of carbon fiber, AFO Light offers a strong and lightweight solution for people with mild to moderate drop foot.

## Indications

- Mild to moderate drop foot
- Contraindications:
- Severe ankle – foot deformities
- Medial-lateral instability
- Severe spasticity
- Extreme activity
- Fluctuating edema

## Features & Benefits

- Supports the foot for ground clearance during swing phase
- Carbon fiber design provides excellent strength-to-weight ratio and dynamic motion
- Stiffness gradient at heel and the toe correspond to ground reaction forces as they are applied to the foot
- Easy to don and doff
- Weight limit: 100kg



## Ordering Information

Part #	Side	Size	Men's Shoe	Women's Shoe	Foot Plate Length	Calf Height
O-121012	Right	Small	up to 7.5	7 - 9	22cm	32 - 36 cm
O-121022	Left	Small	up to 7.5	7 - 9	22cm	32 - 36 cm
O-121014	Right	Medium	8 - 10.5	9.5 - 12	23.5cm	34 - 38 cm
O-121024	Left	Medium	8 - 10.5	9.5 - 12	23.5cm	34 - 38 cm
O-121016	Right	Large	11 - 14	12+	26cm	36 - 40 cm
O-121026	Left	Large	11 - 14	12+	26cm	36 - 40 cm
O-121018	Right	X-Large	12 - 15	12+	27.5cm	38 - 42 cm
O-121028	Left	X-Large	12 - 15	12+	27.5m	38 - 42 cm

\* AFO Light Accessory Kit includes foam padding, upper hook-and-loop strap, lower hook-and-loop strap, hook-and-loop dots.

## Accessories

Part #	Replacement Softgoods	Size
O-122101	AFO Light Softgoods	Universal