

## BALANCE™ FOOT S



UNITY  
COMPATIBLE

Balance Foot S is a lightweight, waterproof foot specifically designed to meet the needs of low active walkers. For these users, the cushioning heel foam with gradual stiffening and wide sole blade provides enhanced stability perception and the perception of a smooth rollover.

The anatomical foot cover features a steady-grip sole and a sandal clamping toe, for the opportunity to wear different types of footwear.

Össur recommends using Balance Foot S in conjunction with the sleeveless elevated vacuum suspension system Unity®, which adds minimal weight and no build height, and does not interfere with Balance Foot S dynamics. Unity is designed to provide secure suspension in combination with Iceross Seal-In® technology, while stabilising volume fluctuations and maintaining comfortable socket fit throughout the day.



ÖSSUR FEET

### User information

Amputation Level:	Trans-femoral and Trans-tibial
Impact Level:	Low
Maximum Patient Weight:	147kg (324lbs)

### Unity System Information

Maximum Vacuum:	24inHg (609mmHg)
Normal Operating Pressure:	-15inHg to -23inHg
Added System Weight:	130g (0.29lbs) w/ Pump, Socket Valve and Tube

### Foot information

Categories:	1-8
Size:	22-30
Weight of Foot: (Size 27, cat5)	622 g (1,37 lbs) w/ Foot Cover and Sock
Build Height: (Size 27)	126mm (5") w/ Male Pyramid and Foot Cover
Heel Height:	10mm (3/8")
Adapter Options:	Male Pyramid



FULL LENGTH  
TOE LEVER



MULTI-AXIAL  
GROUND  
COMPLIANCE



SPLIT TOE



SANDAL TOE



WATERPROOF

## ORDERING INFORMATION

### FOR COMPLETE FOOT



# BALANCE™ FOOT S


## ORDERING INFORMATION

Please select the appropriate Unity Pump Module kit according to foot module size when placing an order and Unity will come preassembled on the foot.

For extra convenience and system compatibility, Unity Valve Kits are highly recommended as a means of connecting to the socket.

For liner selection charts and item numbers please refer to the appropriate Iceross Seal-In page.


### FOR COMPLETE FOOT WITH UNITY



## CATEGORY SELECTION CHART

Weight kg	45-52	53-59	60-68	69-77	78-88	89-100	101-116	117-130	131-147
Weight lbs	99-115	116-130	131-150	151-170	171-194	195-220	222-256	257-287	288-324
Low Impact Category									
Size 22-24	1	1	2	3	4	5	N/A	N/A	N/A
Size 25-27	1	1	2	3	4	5	6	7	8
Size 28-30	3	3	3	3	4	5	6	7	8


## CLEARANCE



# BALANCE™ FOOT S

## KITS


### BALANCE FOOT S WITH MALE PYRAMID




Includes:

- Foot Module
- Removable Foot Cover\* with a black Spectra Sock

\*When ordering brown covers, add "BR" to the end of the part number.




### BALANCE FOOT S WITH MALE PYRAMID AND UNITY PUMP



Includes:




- Foot Module
- Removable Foot Cover\* with a black Spectra Sock
- Unity Pump

\*Note: Unity pump will come preassembled on the foot. Valve kit needs to be ordered separately.



## UNITY

Picture	Part#	Description	Applicable Foot Kits
	UBS0024	Unity Pump S22-24	Balance Foot S & Balance Foot S Torsion Size 22-24
	UBS0027	Unity Pump S25-27	Balance Foot S & Balance Foot S Torsion Size 25-27
	UBS0030	Unity Pump S28-30	Balance Foot S & Balance Foot S Torsion Size 28-30
	UTT0001	Unity TT Valve Kit	All
	UTF0001	Unity TF Valve Kit	All
	UTF0006	Unity TF Tube Kit	All
	L-544611	Icelock 544 Unity Kit	All
	UTT1103	Unity Valve Air Filter Pack (10 p.)	All



## BALANCE™ FOOT S

### UNITY SPARE COMPONENTS

Part#	Description
UBS00021	Unity BS Components 22-30


Includes:

- Outlet Port
- Inlet Port
- Duckbill Valves
- Tube
- Tube Clamp
- Loctite

### COSMETIC FINISHING

#### FOOT COVER KIT, BEIGE AND BROWN\*

PART#



Includes:

- Foot Cover\*
- Foot Cover Attachment
- Black Spectra Sock

\*For brown foot cover please add a "BR" suffix to the part number.

#### FLEX-FOOT SOCKS

Used to cover the composite foot module to protect the foot and eliminate noise between foot module and foot shell.

Part#	Description	Size	Colour
FCX63022	FF Sock Small 1 pc	22-25	Black
FCX63025	FF Sock Large 1 pc	26-30	Black
FCX63022-15	FF Sock Small 15 pcs	22-25	Black
FCX63026-15	FF Sock Large 15 pcs	26-30	Black

