

CTi®3

CTi®3 brings a new vision to double upright knee braces. Innovative engineering makes CTi3 an excellent option for most of your patient population. Suitable for both everyday activities and sport, CTi3 offers great quality with its enhanced construction and anatomical fit for elevated comfort.



INDICATIONS

- For knee conditions that may benefit from increased AP and ML stability and protection, such as injuries to the ACL, MCL, LCL, PCL, rotary and combined instabilities.

FEATURES AND BENEFITS

- Enhanced construction – Lightweight, high strength carbon composite frames provide the optimal stiffness to weight ratio
- Anatomical fit – Designed to fit a wider range of leg shapes and includes an Accutrak® hinge to closely follow the knee's anatomical movement
- Elevated comfort – Unique liners enhance comfort with their special tread block pattern, which provides a balance of cushion and ventilation
- Easy-to-use accessories and replacement parts, including MX kits, replacement strap sets, liner sets and flexion stops

REQUIRED MEASUREMENTS

- Knee caliper
- Thigh circumference 20 cm (8") above mid patella

ORDERING INFORMATION

CTi3

Size	Caliper	Circumference	Left	Right
S	8-10 cm	40-47 cm	CTI-23002	CTI-22002
M	10-11 cm	45-55 cm	CTI-23003	CTI-22003
L	11-13 cm	52-61 cm	CTI-23004	CTI-22004
XL	13-14 cm	58-67 cm	CTI-23005	CTI-22005

Includes: Condyle Pad Set and Extension Stop Set, 0° Extension Stop and Standard Condyle Pad are pre-installed.

Note: M/L measurements are the preferred and most accurate form of measurement for CTi. Circumferential measurements are approximate.

CTi®3

ACCESSORIES & REPLACEMENT PARTS

CTi3 MX KIT

Size	Left	Right
S	CTI-23132	CTI-22132
M	CTI-23133	CTI-22133
L	CTI-23134	CTI-22134
XL	CTI-23135	CTI-22135

Includes: Impact Guard, Sleeve with Mesh, CTi Decal Sheet and CTi Product Bag.

IMPACT GUARD

Size	Left	Right
S-M	CTI-23112	CTI-22112
L-XL	CTI-23115	CTI-22115

LINER SET

Size	Left	Right
S	CTI-23322	CTI-22322
M	CTI-23323	CTI-22323
L	CTI-23324	CTI-22324
XL	CTI-23325	CTI-22325

HOOK SET

Size	Left	Right
S	CTI-23332	CTI-22332
M	CTI-23333	CTI-22333
L	CTI-23334	CTI-22334
XL	CTI-23335	CTI-22335

STRAP SET

Size	Part #
S	CTI-24222
M	CTI-24223
L	CTI-24224
XL	CTI-24225

MISCELLANEOUS

Description	Part #
Guard Elastic Cord	CTI-24130
Guard Screws	CTI-24140
Tibial Wedge Pad	22077
Condyle Pad Set (Thin, Standard, Thick)	CTI-24320
Extension Stop Set (0°, 5°, 10°, 20°, 30°)	CTI-24430
Flexion Stop Set (0°, 60°, 75°, 90°)	CTI-24440
Hinge Sticker Set	CTI-24450
D-Ring Set	CTI-24510
CTi3 PCL Strap Kit	CTI-24520

SLEEVE WITH MESH

Size	Part #
S	SLV60012
M	SLV60013
L	SLV60014
XL	SLV60015