

CHEETAH® JUNIOR WITH NIKE SOLE



Custom-built, high-performance carbon fibre foot designed primarily for sporting activities.

This is the optimal run and play foot for both Transtibial and Transfemoral users. It attaches posterior to the socket, making it agile, strong and a proven performer for those children who want to practice sports.

The Nike Sole is designed to offer increased stability and smooth natural rollover.



User Information	
Amputation Level:	Transfemoral and Transtibial
Impact Level:	High to Extreme
Maximum Patient Weight:	55Kg



PROPORTIONAL
RESPONSE



ACTIVE TIBIAL
PROGRESSION



WATERPROOF


Foot information	
Categories:	1-4
Sizes:	One Size
Weight:	232g
Weight of Sole:	114g
Build Height:	Min:182 mm Max: 307 mm
Adapter Options:	Lamination Connector

CATEGORY SELECTION CHART

Weight kg	15-24	25-34	35-45	46-55
High to Extreme Impact Level	1	2	3	4

CHEETAH® JUNIOR WITH NIKE SOLE

DIMENSIONS




*Note: Sole adds 10 mm to build height.

KITS

CHEETAH KIT


PART#	CATEGORY (1-4)
CJX0	00



NIKE SOLE JUNIOR VERSION KIT

Part#	Size	Cat.
CM150164	One Size	1-4

- Includes:
- Nike Sole
 - Medallion
 - Loctite
 - Inside Sole
 - Bag for Nike Sole



COMPONENTS

Part#	Description
CHX50001	Connector Kit

