

## EXOFORM® WRIST

### INDICATIONS FOR USE

- Mild to Moderate sprains and strains
- Tendonitis
- Post-cast healing
- Other select soft tissue injuries

### FEATURES AND BENEFITS

- Ergonomically contoured padding
- Rigid plastic supports contour to the hand and wrist
- Adjustable aluminum palmar stay allows for fine-tuning when needed

### REQUIRED MEASUREMENTS

- Wrist circumference



### ORDERING INFORMATION

| Size    | Wrist Circumference         | Left   | Right  |
|---------|-----------------------------|--------|--------|
| X-Small | 5.25 – 6.25" (13 – 16 cm)   | 507082 | 507072 |
| Small   | 6.25 – 7.25" (16 – 18.5 cm) | 507083 | 507073 |
| Medium  | 7.25 – 8" (18.5 – 20.5 cm)  | 507085 | 507075 |
| Large   | 8 – 9" (20.5 – 23.5 cm)     | 507087 | 507077 |
| X-Large | 9 – 10" (23.5 – 25.5 cm)    | 507088 | 507078 |

## EXOFORM® CARPAL TUNNEL WRIST

### INDICATIONS FOR USE

- Carpal tunnel syndrome
- Tendonitis

### FEATURES AND BENEFITS

- Allows for full MCP flexion
- Contoured for anatomical fit
- Provides radial and ulnar styloid relief
- Low profile
- Preshaped into a neutral position for comfortable limited motion/ flexibility

### REQUIRED MEASUREMENTS

- Wrist and palm circumference



### ORDERING INFORMATION

| Size   | Wrist Circumference        | Palm Circumference          | Left   | Right  |
|--------|----------------------------|-----------------------------|--------|--------|
| Small  | 5.25 – 6.5" (13 – 16.5 cm) | 6.75 – 8" (17 – 20.5 cm)    | 517083 | 517073 |
| Medium | 6.5 – 7.75" (16.5 – 19 cm) | 8 – 9.25" (20.5 – 23.5 cm)  | 517085 | 517075 |
| Large  | 7.75 – 10" (19 – 25.5 cm)  | 9.25 – 10.5" (23.5 – 26 cm) | 517087 | 517077 |