

GEL ANKLE STIRRUP

INDICATIONS FOR USE

- Ankle injury rehabilitation / immobilization for grades I & II strains and sprains
- Chronic Instability
- Post Cast healing

FEATURES AND BENEFITS

- Safe hot/cold therapy
- Patented lace-lock feature for added stability
- Available 50% surface or 80%



ORDERING INFORMATION

Size	50% Gel Surface		80% Gel Surface	
	Left	Right	Left	Right
Youth/Pony	E200PL	E200PR		
Adult	E5050L	E5050R	E8020L	E8020R

EXOFORM® DORSAL NIGHT SPLINT

INDICATIONS FOR USE

- Plantar fasciitis
- Post-static pain
- General heel pain
- Drop foot
- Achilles tendonitis
- Lower extremity overuse injuries

FEATURES AND BENEFITS

- Adjustable aluminum stays can be positioned to maintain flexion between 85° – 90°
- Soft, Flex edge over molding improves patient comfort
- Full calf and heel strapping prevent brace from rotating
- Slip resistant padding



REQUIRED MEASUREMENTS

- Shoe size

ORDERING INFORMATION

Size	Men's Shoe Size US	Women's Shoe Size US	Shoe Size UK	Shoe Size EU	Part#
Medium	5 – 9	6 - 10	4.5 – 8	37.5 – 42.5	W-50085
Large	9.5 – 14	10.5 – 14	8.5 – 13.5	43 – 48.5	W-50087