

FLEX-RUN™ JUNIOR WITH NIKE SOLE

The Flex-Run Junior is a custom designed solution for children. Its vertical compliance and efficient energy return are musthaves for high active children.

The Flex-Run Junior offers increased cushioning and lower-frequency response makes this foot suitable for Transfemoral and Transtibial users.

The Nike Sole is designed to offer stability and smooth natural rollover.



| User Information | |
|-------------------------|------------------------------|
| Amputation Level: | Transfemoral and Transtibial |
| Maximum Patient Weight: | 84Kg (185lbs) |



WATERPROOF


| Foot Information | |
|------------------|--------------------------------|
| Categories: | 1-6 |
| Sizes: | One Size |
| Weight of Foot: | 336g (11.9 oz) |
| Build Height: | Min: 200 mm (7 7/ 8") |
| Adapter Options: | Male or Female Bolt on Pyramid |

CATEGORY SELECTION CHART

| | | | | | | |
|------------|-------|-------|-------|---------|---------|---------|
| Weight kg | 15-24 | 25-34 | 35-45 | 46-55 | 57-67 | 68-84 |
| Weight lbs | 33-53 | 54-75 | 76-99 | 100-123 | 124-147 | 148-185 |
| Category | 1 | 2 | 3 | 4 | 5 | 6 |

FLEX-RUN™ JUNIOR WITH NIKE SOLE

CLEARANCE




*Note: Sole adds 10 mm (3/8") to build height.

KITS

FLEX-RUN KIT


| | |
|--------|----------------|
| PART# | CATEGORY (1-6) |
| FSX007 | |



NIKE SOLE JUNIOR VERSION KIT


| Part# | Size | Cat. |
|----------|----------|------|
| CM150164 | One Size | 1-6 |

- Includes:
- Nike Sole
 - Medallion
 - Loctite
 - Inside Sole
 - Bag for Nike Sole



COMPONENTS

| Part# | Description | Cat. |
|----------|------------------|------|
| FWX40001 | Attachment Plate | 1-6 |



FOOT ADAPTER SETS

| Picture | Part# | Description | Weight | Bolt Size | Cat. |
|---------|----------|--|--------|------------|------|
| 1. | FWX00031 | Junior Male Pyramid Foot Adapter Set | 30g | M8 x 30mm | 1-3 |
| 2. | FWX00030 | Adult Male Pyramid Foot Adapter Set* | 80g | M10 x 40mm | 1-6 |
| 3. | CM150161 | Adult Female Pyramid Foot Adapter Set* | 100g | M10 x 35mm | 1-6 |

*Adult componentry is used for Flex-Run Junior

REPLACEMENT BOLT

| Part# | Size | Cat. |
|----------|-------------|------|
| CM130221 | M10 x 40 mm | 1-6 |
| CM130103 | M10 x 35 mm | 1-6 |
| CM130237 | M8 x 30 mm | 1-3 |

