

## FORMFIT® KNEE, HINGED

The Formfit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.

The Formfit Knee brace comes in either a wrap or a sleeve version.



### INDICATIONS FOR USE

For knee conditions that may benefit from increased stability around the knee, such as:

- Mild medial/lateral knee instabilities

### FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Universal buttress (Sleeve only)
- Open popliteal

### REQUIRED MEASUREMENTS

- Measure thigh circumference 6" above mid-patella

### ORDERING INFORMATION

Size	Thigh Circumference	Sleeve	Wrap
XS	13" – 15.5" (33 – 39 cm)	504252	503252
S	15.5 – 18.5" (39 – 47 cm)	504253	503253
M	18.5 – 21" (47 – 53 cm)	504255	503255
L	21 – 23.5" (53 – 60 cm)	504257	503257
XL	23.5 – 26.5" (60 – 67 cm)	504258	503258
XXL	26.5 – 29.5" (67 – 75 cm)	504259	503259
3XL	29.5 – 32" (75 – 81 cm)		503260

### Accessories

Description	Part#
Formfit Knee Brace Buttress (Sleeve only)	503000