

## FORMFIT® KNEE ROM

The Formfit® Knee brace is for patients suffering from mild to moderate strains, sprains or ligament tears (ACL/PCL/MCL/LCL). It provides excellent stability, protection and compression during injury recovery. Made from strong and breathable materials it offers comfort while remaining highly durable. The open popliteal area reduces material bulging and potential irritation, while the malleable rigid hinge arms provide an intimate fit.



The Formfit Knee brace comes in either a wrap or a sleeve version.

### INDICATIONS FOR USE

- Mild to moderate ACL and PCL tears
- MCL and LCL tears
- Combined mild knee sprains and strains
- Instabilities of the knee

### FEATURES AND BENEFITS

- Breathable material
- Easy pull-on
- Malleable hinges
- Open popliteal
- Flexion stops: 0°, 45°, 75°, 90°
- Extension stops: 0°, 20°, 30°, 40°, lockout at full extension

### REQUIRED MEASUREMENTS

- Measure thigh circumference 6" above mid-patella

### ORDERING INFORMATION

Size	Knee Circumference	Sleeve		Wrap	
		Long ROM Sleeve	Short ROM Sleeve	Long ROM Wrap	Short ROM Wrap
XSmall	13 – 15.5" (33 – 39 cm)	508352	508252	507352	507252
Small	15.5 – 18.5" (39 – 47 cm)	508353	508253	507353	507253
Medium	18.5 – 21" (47 – 53 cm)	508355	508255	507355	507255
Large	21 – 23.5" (53 – 60 cm)	508357	508257	507357	507257
XLarge	23.5 – 26.5" (60 – 67 cm)	508358	508258	507358	507258
XXLarge	26.5 – 29.5" (67 – 75 cm)	508359	508259	507359	507259
XXXLarge	29.5 – 32" (75 – 81 cm)			507360	507260

### ACCESSORIES

Description	Part#
Universal Buttress	503000