

FORMFIT® TRACKER

The Formfit® Tracker with micro-adjustable ReflexWing™ buttress realigns the patella, providing dynamic lateral support for common patella conditions. Its FlyWeight™ fabric delivers lightweight breathability and a low-profile design that feels great against the skin and minimises migration. PowerLock™ straps work in tandem with heat moldable CustomFit™ hinges to provide a customisable fit for your patient.



INDICATIONS FOR USE


For knee conditions that may benefit from realignment of the patella and/or patients that would benefit from an added lateral support of the patella, such as:

- lateral patellar subluxation or dislocation
- lateral maltracking

FEATURES AND BENEFITS

- CustomFit™ hinges for customised fit and support
- Anti-Migration Tech™ minimises migration
- PowerLock™ straps for secure lateral support
- ReflexWing™ buttress for comfortable patella realignment
- FlyWeight™ fabric for soft feel and breathability
- Durable and machine washable

REQUIRED MEASUREMENTS



Measure calf circumference 15 cm below mid-patella. Leg in extension.



Measure thigh circumference 15 cm above mid-patella. Leg in extension.

FORMFIT® TRACKER

ORDERING INFORMATION

Size	Calf circumference 15 cm below mid patella	Thigh circumference 15 cm above mid patella	Left	Right
XS	25 - 30	33 - 39	719162	719152
S	30 - 35	39 - 47	719163	719153
M	35 - 40	47 - 53	719164	719154
L	40 - 46	53 - 60	719165	719155
XL	46 - 51	60 - 67	719166	719156
XXL	51 - 56	67 - 75	719167	719157
3XL	56 - 61	75 - 81	719168	719158

ACCESSORIES

Description	Part#
ReflexWing™ Buttress	719100