

# FORMFIT® ANKLE WITH SPEEDLACE

## INDICATIONS FOR USE

- Ankle injury rehabilitation/immobilization (Grade I/II sprains)
- Support chronic instability
- Prophylactic use


## FEATURES AND BENEFITS

- Durable materials
- Quick-Pull closure straps
- Figure 8 straps
- Removable stays
- Circumferential cuff



## REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone. Ankle in 90 degrees.



## ORDERING INFORMATION

Size	Ankle Circ.	Part#
XX-Small	23 – 25.5cm	W-10621
X-Small	25.5 – 28cm	W-10622
Small	28 - 30.5	W-10623
Medium	30.5 – 33cm	W-10625
Large	33 – 35.5cm	W-10627
X-Large	35.5 – 38cm	W-10628
XX-Large	38 – 40.5cm	W-10629