

FORMFIT® ANKLE WITH SPEEDLACE

The Formfit® Ankle Brace with Speedlace provides stability and protection for the ankle – whether used prophylactically, for chronic instability or following an injury. Made from durable materials, its low profile, quick-lace design is comfortable and fits easily into shoes. Removable semi-rigid stays provide additional support and the Figure 8 straps are designed for inversion/eversion control.

INDICATIONS FOR USE

- Grade I-II sprains and strains
- Chronic instability
- Prophylactic use

FEATURES AND BENEFITS

- Durable materials
- Quick-Pull closure straps
- Figure 8 straps
- Removable stays
- Circumferential cuff

REQUIRED MEASUREMENTS

- Measure ankle circumference diagonally over the base of the heel and across the ankle bone.



ORDERING INFORMATION

Size	Ankle Circ.*	Part#
XXSmall	9 – 10" (23 cm – 25.5 cm)	W-10621
XSmall	10 – 11" (25.5 cm – 28 cm)	W-10622
Small	11 – 12" (28 cm – 30.5 cm)	W-10623
Medium	12 – 13" (30.5 cm – 33 cm)	W-10625
Large	13 – 14" (33 cm – 35.5 cm)	W-10627
XLarge	14 – 15" (35.5 cm – 38 cm)	W-10628
XXLarge	15 – 16" (38 cm – 40.5 cm)	W-10629

*Measure ankle circumference diagonally over the base of the heel and across the ankle bone.