

FORMFIT® HINGED LATERAL J

INDICATIONS FOR USE

- Lateral patellar subluxation and dislocations
- Lateral patellofemoral malalignment
- Patellofemoral maltracking and pain
- Chondromalacia due to lateral maltracking

FEATURES AND BENEFITS

- Breathable material
- Removable hinges
- Easy pull-on
- Lateral J buttress
- Open popliteal



REQUIRED MEASUREMENTS

- Thigh Circumference: 6" above mid-patella
- Calf Circumference: 6" below mid-patella

ORDERING INFORMATION

Size	Thigh circumference	Knee circumference	Calf circumference	Right leg	Left leg
X-Small	13 - 15.5" (33 - 39.4cm)	12-13" (30.5-33cm)	10-12" (25.4-30.5cm)	709152	709162
Small	15.5-18.5" (39.4-47cm)	13-14" (33-35.6cm)	12-14" (30.5-35.6cm)	709154	709164
Medium	18.5-21" (47-53.3cm)	14-15" (35.6-38.1cm)	14-16" (35.6-40.6cm)	709155	709165
Large	21-23.5" (53.3-59.7cm)	15-17" (38.1-43.2cm)	16-18" (40.6-45.7cm)	709157	709167
X-Large	23.5-26.5" (59.7-67.3cm)	17-19" (43.2-48.3cm)	18-20" (45.7-50.8cm)	709158	709168
XX-Large	26.5-29.5" (67.3-74.9cm)	19-21" (48.3-53.3cm)	20-22" (50.8-55.9cm)	709159	709169